GC Fuji II LC™ CAPSULE

RADIOPAQUE LIGHT CURED REINFORCED GLASS IONOMER RESTORATIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Restoration of Class III, V and limited Class I cavities.
2. Restoration of primary teeth.
3. Core build up.
4. Cases where a radiopaque restoration is required.
5. Geriatric applications.
6. As a base or liner.

CONTRAINDICATIONS

1. Pulp capping.
2. Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0.33 / 0.10
Mixing Time (sec.)	10"
Working Time (min., sec.)	3'15"
Light Curing Time (sec.)	20"
Depth of Cure (A2) (mm)	1.8

Test conditions : Temperature (23 ±1°C)
Relative humidity (50 ±10%)

1. CAVITY PREPARATION

- a) Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- b) Wash and dry but DO NOT DESICCATE. CAVITY CONDITIONER is recommended to remove the smear layer (Fig. 1).

2. CAPSULE ACTIVATION AND MIXING

- a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- b) To activate the capsule, push the plunger until it is flush with the main body and hold it down for 2 seconds (Fig. 3).
Note:
1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
2) The capsule should be activated just before mixing and used immediately.

- c) Immediately remove the capsule and set it into a mixer (an amalgamator) and mix for 10 seconds at high speed (~4,000 RPM) (Fig. 4).

3. RESTORATIVE TECHNIQUE

- a) Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- b) Click twice to prime the capsule then syringe (Fig. 5).
The working time is 3 minutes 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
Note:
1) To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
2) To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.
- c) Remove surface moisture but DO NOT DESICCATE.
- d) Extrude cement directly into preparation (Fig. 6). Avoid air bubbles.
- e) Form the contour and place a matrix if required.
- f) Light-cure for 20 seconds using a suitable visible light curing device (470nm wavelength) (Fig. 7). Place light source as close as possible to the cement surface.
Note:
For cavities deeper than 1.8 mm, use a layering technique.

4. FINISHING

- 1) Finish under water spray using superfine diamond bur, silicone point and polishing strips (Fig. 8).

Note:
Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 9).

STORAGE

Recommended for optimal performance store the capsules in the original aluminium foil in a cool and dark place (4-25°C) (39.2-77°F).

SHADE

A1, A2, A3, A3.5, A4, B2, B3, B4, C2 and C4 (10 shades based on Vita® shades).

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

1. Single shade package : 48 capsules.
2. Dark Shade Assortment package : 24 capsules.
5 pcs each of A3.5, A4, B3 and B4; 4 pcs each of C4.
3. Light Shade Assortment package : 24 capsules.
5 pcs each of A1, A2, A3 and B2; 4 pcs each of C2.
*Average contents per capsule : 0.33g powder and 0.10g (0.085mL) liquid.
*Minimum net volume of mixed cement per capsule : 0.10mL.
4. Single Shade Export package: 200 capsules
Option : GC Capsule Applier (1 piece).

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at: <http://www.gcamerica.com>
They can also be obtained from your supplier.

Last revised: 11/2019

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8568, Japan
RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
DISTRIBUTED BY
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
www.gcamerica.com
TEL: +1-708-597-0900

GC SOUTH AMERICA
Rua Helicóptero, 399, Santana - São Paulo, SP - Brasil
CEP: 02022-051 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP. TEC. Mayara De Santis Ribeiro - CRO/SP: 105.982

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