GC FIBER POST RADIOPAQUE ESTHETIC POST

Primo use, carefully read the instructions for use. (**EN**)

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

GC FIBER POST DRILL is used to prepare post space. GC FIBER POST, composed of glass fiber reinforced composite, is intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

CONTRAINDICATIONS

- A curved or twisted root canal where the post cannot be seated straight.
- In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Fiber Post					
	Φ 0.8 mm	Φ 1.0 mm	Φ 1.2 mm	Φ 1.4 mm	Φ 1.6 mm
Color code	Gray	Purple	Red	Yellow	Blue
Length	22 mm				
Length of the tapered part	2.2 mm	3.5 mm	4 mm	5 mm	6 mm
Tapered angle	4°				
Greatest dimension	0.8 mm	1.0 mm	1.2 mm	1.4 mm	1.6 mm
Smallest dimension	0.5 mm	0.5 mm	0.5 mm	0.6 mm	0.7 mm

Drill			
	1.2	1.4	1.6
Color code	Red	Yellow	Blue
Drill blade length	10 mm		
Drill blade diameter	1.3 mm	1.5 mm	1.7 mm

DIRECTIONS FOR USE

- Select the appropriate post and drill size by radiographic examination. Determine the required preparation depth.
- Remove any remaining root filling material (gutta percha) using e.g. a peeso reamer leaving at least 4mm at the apical end.
- Prepare the post space to the pre-planned depth using the selected drills (2000-5000rpm) with moderate pressure using water spray, incrementally starting with the narrowest drill. For post size 0.8 and 1.0, use a peeso reamer or the narrowest drill to prepare the post space.
- The preparation should include at least 1.5 mm in height of sound tooth structure circumferentially. **The thickness of the sound tooth structure should be at least 1.0 mm**.
 - Clean and autoclave drills before the first use. After use, follow the instruction below; Keep the water pouring into a container and soak drills in water to remove the blood or debris from the surface of the drills. Make sure to keep the drills away from direct water stream. Wash all blood or debris from drills under water, using a nylon brush. Clean drills using an ultrasonic cleaner for 5 minutes in chlorhexidine gluconate solution or glutaric aldehyde solution. Remove the cleaning agent thoroughly using a brush under running water. Soak drills in purified (distilled) water in a clean container and clean with ultrasonic cleaner for 5 minutes. After cleaning, immediately dry drills thoroughly. Sterilize by autoclaving at 121°C for 15 minutes or 126°C for 10 minutes under moist condition. **Dry heat sterilization or chlorine disinfectant such as sodium hypochlorite cannot be used.**
 - Do not use a drill which is damaged, deformed, rusted or contaminated.
 - Ensure that the drill is securely attached to the hand piece. Operate the hand piece outside the patient’s mouth to check the drill is attached correctly.
 - Use protective eye wear while drilling.
 - Take care not to put inappropriate force on the drill.
 - Clean and disinfect the post space according to standard techniques. Rinse and dry thoroughly with oil free air and use paperpoints to remove any remaining liquid

NOTE: Make sure no residues of temporary materials or chemical disinfectant remains in the post space. .
 - Remove the post from the blister and place in the post space to check the fit. Mark the final length.

NOTE: About 2/3 of the post should be located in the root space and 1/3 in the core structure.
 - Remove the post and cut the post to the desired length using a diamond disc or bur. Scissors or similar tools cannot be used. This may damage the post.
 - Clean the post with cotton soaked with alcohol. Dry with oil-free air.
 - Treat the surface of the post with a silane coupling agent such as GC Ceramic Primer.** Refer to the respective manufacturer’s instructions for use.

NOTE: After silanization, use forceps to handle the post and avoid direct contact with fingers. If the surface is contaminated, clean with cotton soaked with alcohol and treat with silane coupling agent again.
 - Prepare the post space according to the respective manufacturer’s instructions of luting or core build-up material. Make sure the post space is clean and dry.
 - Fill the post space with the luting or core build-up material. Seat and cement the post according to the respective manufacturer’s instructions.
 - Complete the core with core build-up material.

CAUTION

- GC FIBER POST is for single use only. Do not reuse.
- Take care not to drop the post in patient’s mouth.
- Use of a rubber dam is recommended to isolate tooth and to prevent accidental ingestion of the post.
- GC FIBER POST and GC FIBER POST DRILL are not compatible with any other post or drill.
- Do not leave drills wet. Store after thorough drying.
- Use a dust collector, wear a dust mask and protective eyewear when using discs or burs.

STORAGE

GC FIBER POST: Store away from direct sunlight.

GC FIBER POST DRILL: Store in a clean and dry area.

PACKAGES

- 400001** - GC FIBER POST ASSORTMENT KIT
Post: 5 pieces each of 1.2 mm, 1.4 mm, 1.6 mm (in single-dose blister packaging)
Drill: 1 piece each of 1.2 mm, 1.4 mm, 1.6 mm
- 400101** - GC FIBER POST ASSORTMENT KIT
Post: 5 pieces each of 1.0 mm, 1.2 mm, 1.4 mm (in single-dose blister packaging)
Drill: 1 piece each of 1.2 mm, 1.4 mm
- Refill of GC FIBER POST 10 pieces (available in all sizes, 0.8 mm, 1.0 mm, 1.2 mm, 1.4 mm, 1.6 mm) (in single-dose blister packaging)
- Refill of GC FIBER POST DRILL 1 piece (available in 3 sizes - 1.2mm, 1.4mm, 1.6mm)

- Fill the post space with the luting or core build-up material. Seat and cement the post according to the respective manufacturer’s instructions.
- Complete the core with core build-up material.

VORSICHT

- GC FIBER POST is for single use only. Do not reuse.
- Take care not to drop the post in patient’s mouth.
- Use of a rubber dam is recommended to isolate tooth and to prevent accidental ingestion of the post.
- GC FIBER POST and GC FIBER POST DRILL are not compatible with any other post or drill.
- Do not leave drills wet. Store after thorough drying.
- Use a dust collector, wear a dust mask and protective eyewear when using discs or burs.

STORAGE

GC FIBER POST: Store away from direct sunlight.

GC FIBER POST DRILL: Store in a clean and dry area.

PACKAGES

- 400001** - GC FIBER POST ASSORTMENT KIT
Post: 5 pieces each of 1.2 mm, 1.4 mm, 1.6 mm (in single-dose blister packaging)
Drill: 1 piece each of 1.2 mm, 1.4 mm, 1.6 mm
- 400101** - GC FIBER POST ASSORTMENT KIT
Post: 5 pieces each of 1.0 mm, 1.2 mm, 1.4 mm (in single-dose blister packaging)
Drill: 1 piece each of 1.2 mm, 1.4 mm
- Refill of GC FIBER POST 10 pieces (available in all sizes, 0.8 mm, 1.0 mm, 1.2 mm, 1.4 mm, 1.6 mm) (in single-dose blister packaging)
- Refill of GC FIBER POST DRILL 1 piece (available in 3 sizes - 1.2mm, 1.4mm, 1.6mm)

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