GC Fuji ORTHO™ LC

LIGHT-CURED ORTHODONTIC BONDING ADHESIVE

For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

- Bonding metal brackets and attachments.
- Bonding ceramic brackets.
- Bonding acrylic appliances.
- Band cementation only when extra band retention is desired.

CONTRAINDICATIONS

- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.
- Bonding polycarbonate brackets.
- Cementation of stainless steel crown-retained appliances (e.g., Herbst Appliances).

DIRECTIONS FOR USE

- Enamel Preparation
 - Using a prophyl cup or a brush, clean the bonding surfaces of the teeth with plain (non-fluorinated) pumice and water.
 - Rinse thoroughly with water.
 - Using a sponge or a cotton pellet, apply ORTHO CONDITIONER to the bonding surfaces of the teeth for 20 seconds.
 - Rinse the surfaces. Although successful results have been obtained with no enamel pretreatment, conditioning of the enamel bonding surfaces will increase the adhesive's bond strength. Rinse thoroughly with water before use.
- Enamel Preparation
 - Using a prophyl cup or a brush, clean the bonding surfaces of the teeth with plain (non-fluorinated) pumice and water.
 - Rinse thoroughly with water.
 - Using a sponge or a cotton pellet, apply ORTHO CONDITIONER to the bonding surfaces of the teeth for 20 seconds.
 - Rinse the surfaces. Although successful results have been obtained with no enamel pretreatment, conditioning of the enamel bonding surfaces will increase the adhesive's bond strength. Rinse thoroughly with water before use.

NOTE: The enamel bonding surfaces must be moist. An overly dry or desiccated enamel surface will adversely affect the bond strength. The optimum level of surface moisture can be obtained by wiping the bonding surfaces of the teeth with a moistened cotton roll immediately prior to bracket bonding. If the bonding surfaces become desiccated during the bonding procedure, rehydrate them with the moistened cotton roll before bonding with GC Fuji ORTHO LC.

- Powder and Liquid Dispensing
 - The standard powder to liquid ratio is 3.0 g / 1.0 g. 1 level large scoop of powder to 2 drops of liquid.
 - For accurate dispensing of powder, tap the bottle gently. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close the bottles immediately after use.

- Mixing
 - For rebonds or a small mix, use the small scoop and one drop of liquid.
 - Divide the powder into 2 equal parts. Mix the first portion with ALL the liquid and mix for about 10 seconds. Incorporate the remaining powder and mix thoroughly for an additional 10-15 seconds (total 20-25 seconds).

NOTE: The working time is approximately 3 minutes from start of mixing at 23° (73°F). Higher temperatures will shorten the working time, and lower temperatures will extend it.

- Bonding Procedure
 - Brackets
 - Coat the bonding surface of the bracket completely with the mixed adhesive.
 - Position the coated bracket on the tooth.
 - Press the bracket firmly against the enamel surface. Using an explorer or a scaler, remove the excess adhesive at this time. Press the bracket again. Continue placing additional brackets. Place all brackets in a quadrant or in the full arch. More than one mix of adhesive may be necessary. If bracket drifts in a corner, "lock" the bracket with light curing for approximately 5 seconds.

AFTER REMOVING EXCESS ADHESIVE, BE CAREFUL NOT TO DISTURB THE BRACKETS BEFORE THE ADHESIVE IS CURED.

MOVEMENT OF THE BRACKETS AT THIS POINT MIGHT DECREASE THE BOND STRENGTH.

- Using a curing light (470nm wavelength), cure each bracket for 10 seconds each from the occlusal, mesial, distal and gingival aspects. It is very important that the curing light be tested periodically with a light meter to ensure adequate light output.

- Insert leveling wire and complete initial procedures. A light force wire (e.g., 016 NITI or equivalent) is recommended at the initial bonding.

- Bonding ceramic brackets: Chemically retentive bases must be etched and silanated. If not etched and silanated by the manufacturer, etch and silanate before bonding with GC Fuji ORTHO LC. If bracket bases become contaminated, follow the manufacturer's directions for cleaning, and/or etching and silanating. No treatment is required for mechanically-retained, lasso type brackets or porcelain restorations.

- Prepare porcelain surfaces by etching and silanating according to manufacturer's directions. Bond with GC Fuji ORTHO LC as per above.

- Bonding brackets to amalgam/alloy restorations: Lightly roughen the metal surface with a disc or a fine diamond point before placing GC Fuji ORTHO LC.

- Acrylic Appliances
 - Using a microetcher or a rotary instrument, roughen and place retention holes in the internal surfaces of the appliance.
 - Prepare the appliance surfaces by etching and silanating according to manufacturer's directions. Bond with GC Fuji ORTHO LC as per above.

- Removal of old appliance: Use a scaler or rotary instrument to remove old appliance, flex the adhesive in multiple sites with band removing pliers. If this does not remove the appliance, section the appliance, desiccate the adhesive, and remove the appliance.

- Storage
 - For optimal performance, store in a cool and dark place (4-25°C / 39-77.0°F).

- Packages
 - Introductory Package: 15 g powder (1), 8 g (6.8 mL) liquid (1) with a double-ended powder scoop, a plastic spatula, a set of mixing pads (No. 22)
 - Standard Package: 40 g powder (1), 8 g (6.8 mL) liquid (2) with a double-ended powder scoop, a plastic spatula, a set of mixing pads (No. 23)
 - Replacement Packages: a) 40 g powder with a double-ended powder scoop b) 8 g (6.8 mL) liquid

- CAUTION
 - In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
 - In case of contact with eyes, flush immediately with water and seek medical attention.
 - Do not mix powder or liquid with any other glass ionomer materials.
 - In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
 - Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gcgroup.com> or for The Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection, immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERS. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EDUMED database (<https://ec.europa.eu/tools/edumed>) or contact us at Regulatory.gcc@dc.dental

Undesired effects-Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gcc.dental. In this way you will contribute to improve the safety of this product.

UK Responsible Person GC UNITED KINGDOM LTD Cooper Court Newport Pagnell, Buckinghamshire, MK16 8US, United Kingdom Last revised: 10/2022

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- Acrylische Apparate
 - Using a microetcher or a rotary instrument, roughen and place retention holes in the internal surfaces of the appliance.
 - Prepare the appliance surfaces by etching and silanating according to manufacturer's directions. Bond with GC Fuji ORTHO LC as per above.

- Entfernung des alten Appliance: Use a scaler or rotary instrument to remove old appliance, flex the adhesive in multiple sites with band removing pliers. If this does not remove the appliance, section the appliance, desiccate the adhesive, and remove the appliance.

- Speicherung
 - For optimal performance, store in a cool and dark place (4-25°C / 39-77.0°F).

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