Prior to use, carefully read the instructions for use.

**GB**

## COMPOSITE PRIMER

PRIMER FOR COMPOSITE BONDING

COMPOSITE PRIMER

COMPOSITE PRIMER is a light-cured priming agent designed to be used in bonding micro-filled hybrid resin prostheses to tooth using resin adhesive cement or in repairing them.

COMPOSITE PRIMER

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**  
Priming the bonding surfaces of micro-filled (MFR) hybrid resin prostheses (made from GC GRADIA, etc.).

**CONTRAINDICATIONS**  
In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

- DIRECTIONS FOR USE**
- Surface preparation of resin prosthesis
    - For bonding resin prosthesis

Treat the bonding surface of the MFR hybrid resin prosthesis using a mild sandblaster (aluminum oxide). Clean with an air syringe.
    - For repairing resin prosthesis

Roughen the bonding surface of the MFR hybrid resin to be repaired using a sandblaster (aluminum oxide) with a pressure of less than 3kg/cm<sup>2</sup> (4,3 psi) or other instrument. Clean with an air gun or other device.

Note :

After cleaning, do not contaminate the bonding surface with extraneous matter or moisture.

- Application and Light curing of COMPOSITE PRIMER
  - Dispense one or two drops of COMPOSITE PRIMER into a dispensing dish. Apply a thin layer to the bonding surface of the MFR hybrid resin prosthesis using a Micro-tip applicator, etc.

Note :

    - If the primer is stored in a refrigerator, remove from refrigerator and leave to stand for at least 5 minutes at room temperature before use.
    - After dispensing, use as soon as possible. Replace the bottle cap immediately after use.
  - Light cure for 20 seconds using a visible light curing unit (470 nm wavelength). After light curing, do not remove the uncured layer. It will serve as a bonding layer between the resin prosthesis and resin adhesive cement.

- Cementation

Apply resin adhesive cement such as GC LINKMAX to the primed surface of the prosthesis or apply new MFR resin material such as GC GRADIA to the primed resin prosthesis to be repaired. For details, refer to manufacture's instructions for use.

**STORAGE**  
Store in a cool, dark place (4-25°C) (39,2-77,0°F).  
(Shelf life : 2 years from date of manufacture)

**PACKAGE**  
Bottle of 3mL

- CAUTION**
- COMPOSITE PRIMER is flammable. Do not use near naked flame. Keep away from sources of ignition. Keep away from direct sunlight. Use in a well ventilated place.
  - In case of contact with oral tissues, immediately remove with a cotton pellet or a sponge. Flush immediately with water after cementing.
  - In case of contact with skin, remove immediately with a cotton pellet or a sponge. Rinse with soap and water.
  - In case of contact with eyes, flush immediately with water and seek medical attention.
  - Take care to avoid ingestion of the material.
  - Do not mix with any other primers.

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.