
# PERFORMANCE TEST KIT

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## KIT DE TESTE DE DESEMPENHO

## KIT DE TEST DE PERFORMANCES

MANUFACTURED BY  
GC EUROPE N.V.  
Interleuvenlaan 13  
B-3001 Leuven, Belgium  
DISTRIBUTED BY  
GC AMERICA INC.  
3737 West 127th Street  
Alsip, Illinois 60803 U.S.A.  
GC ASIA DENTAL PTE. LTD.  
19 Lyons Way, #06-27  
Singapore 508724  
PRINTED IN USA  
5950802

# PERFORMANCE TEST KIT

Prior to use, carefully read the instructions for use.

## IN GENERAL

The Test Performance Kit is designed to measure the effectiveness of a curing light over time by comparing the measured length of cured samples made out of the same VLC material. Ideally, a test should be performed upon receipt of a curing light and subsequently at regular time intervals. When using the same VLC material, cured samples should have similar lengths.

## DIRECTIONS FOR USE

1. Assemble the two parts
2. Select the cavity you want to fill
3. Fill the cavity with the VLC material
4. Light cure the material with your curing light in the preferred mode of cure / curing time
5. Remove the cured material from the cavity by separating the two parts
6. Remove any uncured material from the bottom of the sample using alcohol
7. Measure the length of the cured sample

## IMPORTANT NOTE

When using VLC materials in the mouth, maximum layer thickness recommended by the manufacturer of the VLC material used should be respected. In general, VLC materials should be cured in layers no thicker than 2 to 3 mm.

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