LIGHT-CURED FLOWABLE COMPOSITE

UniFil® Flow is a multi-purpose light-cured, fluoride releasing, micro-filled hybrid type composite resin with excellent flowability. UniFil® Flow allows for effective direct injection using a special dispensing tip, resulting in aesthetic restorations.

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Restoration of Class I, II, III, IV, V cavities (particularly for small Class I cavities / shallow Class V cavities / other small cavities).
- Restoration of root surface caries.
- Restorations in deciduous teeth.
- Filling tunnel shaped cavities.
- Sealing hypersensitive areas.
- Liner / base / filling in cavity undercuts.
- Sealant.
- Fixing loose teeth.
- Additions to composite restorations. (Fig. 1)

**CONTRAINDICATIONS**

- Pulp capping.
- In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of the product and refer to a physician.

**DIRECTIONS FOR USE**

**1. Preparations**

- Hold the syringe upright and remove the wing cap by turning counterclockwise. Take care not to expose material to direct light from the dental lamp or natural light (Fig. 2).
- Promptly and securely attach the dispensing tip (plastic or needle type) to the syringe by turning clockwise (Fig. 2)(Fig. 3).  
Note :  
The plastic type is recommended for use with children.
- After attaching the dispensing tip, protect it with the cover until ready for use in order to avoid exposure to light (Fig. 2).

**2. Cavity Preparation**

Prepare cavity using standard techniques.  
Note :  
For pulp capping, use calcium hydroxide.

**3. Shade Selection**

Select shade from 5 shades of A1,A2,A3,A3.5,A4 based on Vita® Shade and CV (Cervical Color).

\* Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

**4. Use of Light-cured Bonding System**

For bonding UniFil® Flow to tooth structure, use a light-cured bonding system (Fig. 4). GC UniFil® Bond which contains adhesive monomer (4-MET) for superior bondability to the structure is recommended.  
Note :  
When using light-cured bonding system, follow manufacturer's instructions for use.

**5. Placement of UniFil® Flow**

- Remove cover from the dispensing tip on the syringe (Fig. 5).
- To remove any air from the dispensing tip, with the tip pointing upwards gently push forward the syringe plunger until material reaches the mouth of the tip (Fig. 6).  
Note :  
If there is air inside the dispensing tip, air bubbles may be formed at the time of injection.
- Place the dispensing tip as close as possible to the cavity, and slowly push the plunger to inject material into it (Fig. 7). Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.  
Note :  
a. When attaching the dispensing tip, make sure that no material is sticking to the tip and the syringe in order to ensure a tight connection.  
b. If the syringe does not extrude smoothly, remove the dispensing tip and extrude material directly from the syringe to make sure that material is coming out.  
c. The material will start to harden if exposed to the dental light or ambient light.  
Be sure to protect it from light when working from a mixing pad.  
d. After use, immediately remove the dispensing tip and tightly close the syringe with the wing cap.

**Clinical Hint 1**

In order to inject effectively, use the surface tension of the material to ensure uniformly across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

**Clinical Hint 2**

When filling a large cavity, it is recommended to place material incrementally into the cavity. Another effective method is to use UniFil® Flow for filling in undercuts or as a liner/base, and then to place composite resin (GC UniFil S, etc.) on top.

**6. Light Curing**

Light cure the UniFil® Flow using a light curing unit (Fig. 8). Refer to the following chart for Irradiation Time and Effective Depth of Cure.  
Note :  
When light curing material, wear protective glasses.

**Irradiation Time and Effective Depth of Cure**

Irradiation time	Shade	
	A1, A2, A3	A3.5, A4, CV
3 sec. * (Plasma arc)	2.0 mm	1.5 mm
20 sec. ** (Halogen)(LED)		
5 sec. * (Plasma arc)	3.0 mm	2.5 mm
40 sec. ** (Halogen)(LED)		

\* Plasma arc: GC Filpo  
\*\* Halogen : GC New Light VL-II, Co-bee or Coe Lunar TA LED : GC e-Light (Mode : Turbo)

**7. Shaping and Polishing**

Shape and polish using standard techniques.

**SHADES**

A1, A2, A3, A3.5, A4, CV (cervical color)

**STORAGE**

Store in a cool and dark place (4 - 25°C / 39.2 - 77.0°F) (Shelf life : 2 years from date of manufacture)

**PACKAGES**

- Starter Package  
Syringe 1.5g (0.8mL) x 4 (shades of A2, A3, A3.5, CV)  
10 dispensing tips (5 each of plastic / needle types)  
4 light protective covers
- Replacement packages  
a. Refill package (one shade each of A1, A2, A3, A3.5, A4 or CV)  
Syringe 1.5g (0.8mL) x 2  
with 4 dispensing tips (2 each of plastic / needle types), 2 light protective covers.  
b. Dispensing tips package (plastic or needle type)  
Pack of 20 dispensing tips, 2 light protective covers.

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take special care that the patient should not swallow material.
- Do not mix with other similar products.
- The dispensing tip cannot be sterilized in an autoclave or chemically.


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7


Fig. 8

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