G-CERA ORBIT VEST

Revised Sept. 2002

S3300

Refractory die material

Feuerfestes Stumpfmateriale

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Ildfast modelmateriale

GC

Prior to use, carefully read
the instructions.

INSTRUCTIONS FOR USE

G-CERA ORBIT VEST is a phosphate-bonded refractory die material developed for production of all-ceramic crowns, inlays, onlays and veneers. G-Cera Orbit Vest allows various porcelains to be built directly onto the refractory model; has outstanding dimensional stability and provides sufficiently high strength to withstand repeated firings.

FEATURES

1. TWO TYPES OF POWDER OFFERED

G-Cera Orbit Vest has two types of powder with different expansions for crowns/veneers, and for inlays. This provides for a wide range of use in producing both external ceramic restorations such as crowns or veneers, and internal ceramic restorations such as inlays.

2. SUPERIOR STRENGTH

G-Cera Orbit Vest has sufficiently high strength to endure repeated firings. It has been formulated to curb the expansion rate to a minimum in a given range of temperature rise, causing no crack or damage to its model.

3. SMOOTH SURFACE

Thanks to its unique granular formula, G-Cera Orbit Vest provides a smooth model surface.

4. EXCELLENT COMPATIBILITY WITH VARIOUS PORCELAINS

G-Cera Orbit Vest is designed for use with most metal bonding porcelains which have average coefficients of thermal expansion of $11.5-13.5 \times 10^{-6}$ at $100-500^{\circ}\text{C}$ ($212-932^{\circ}\text{F}$). It will not cause the porcelain to crack or separate from the model during firing.

PHYSICAL DATA

| Item | Value |
|----------------------------------------------------|--------------------|
| Powder / Liquid ratio | 33 g / 5 ml (6.5g) |
| Pouring time (min. sec.) (after completion of mix) | 2'00" |
| Setting time (min. sec.) | 7'00" |
| Compressive strength (MPa) (after degassing) | 30.6 |

Test conditions: Temperature (23±2°C) (73.4±3.6°F) - Relative humidity (50±5%).

DIRECTIONS FOR USE:

1. THE POWDER / LIQUID RATIO

The powder/liquid ratio is 33g (one pack of powder) / 5ml (6.5g).

Note:

Powder and liquid should be correctly dispensed. Otherwise proper physical properties such as thermal and setting expansions can not be obtained.

2. MIXING

After hand-mixing thoroughly, vacuum mix uniformly for 60 seconds.

Note:

The smallest possible mixing container (5-6cm internal diameter) should be used. The use of a larger container will not allow for a uniform consistency, resulting in improper physical properties such as thermal and setting expansion.

3. POURING THE IMPRESSION

Prior to pouring the refractory, a wetting agent can be applied to the impression in order to obtain a bubble-free refractory model. After applying the wetting agent, lightly rinse under running water and dry well. After mixing the refractory, immediately pour the mixture into the impression while using a gypsum vibrator (maximum pouring time after mix: about 2 minutes).

Note:

Some wetting agents may roughen the surface of the refractory model.

4. REMOVING FROM THE IMPRESSION

Remove the model from the impression 1-1.5 hour after pouring.

Note:

Avoid early removal as it will adversely effect the strength, dimensional precision and surface conditions of the model. Also do not allow the model to stay in the impression for over 2 hours, or you will have difficulty with separation.

5. DEGASSING THE MODEL

Gradually heat the model in the ring furnace from room temperature up to 700°C (1,292°F) over 1 hour. Maintain (700°C / 1,292°F) temperature for 10 minutes.

Note:

Too rapid a rise in temperature may fracture the model due to too rapid evaporation of moisture inside the refractory model.

Next, place the model in the preheated (700°C / 1,292°F) porcelain furnace and raise to 1,050°C (1,922°F) at a rate of 55°C (131°F) per minute under atmospheric pressure. Maintain this temperature for 5 minutes. Remove the model and bench-cool to room temperature. It is not necessary to bench-cool in the porcelain furnace.

Notes:

- 1) Do not degass the model at a temperature higher than 1,050°C (1,922°F). Otherwise it may affect its dimensional precision.*
- 2) Degassing the model can also be performed in a porcelain furnace, but phosphate-bonded refractory materials give off small amounts of ammonia gas, which may affect the porcelain furnace. Therefore a ring furnace is recommended to heat the model up to 700°C (1,292°F).*

NOTES:

1. When using G-Cera Orbit Vest, use a dust collector and dust mask to avoid inhaling powder dust.
2. Avoid contact of the material with eyes. In case of contact, immediately flush with water and seek medical attention.
3. When cutting and trimming the hardened refractory model, wear goggles for eye protection.
4. When heating the material, work in a well-ventilated room.
5. Do not mix the powder or liquid of G-Cera Orbit Vest with other products.
6. This product is to be used only by a dental professional.
7. This product should be used only for the applications described in the instructions.

Other notes:

1. If the liquid should be contaminated with extraneous matter such as plaster, acid, table salt, etc., or should be placed below 5°C (41°F), it will turn into a gel and can not be used again.
2. Keep both powder and liquid away from direct sunlight. Store at room temperatures (20-25°C / 68-77°F).
3. After opening the powder pack, use material up as soon as possible.
4. Some components contained in the powder will delay the setting time of gypsums or gypsum bonded investments. Designate a mixing bowl to be used exclusively with G-Cera Orbit Vest.

PACKAGES:

1. Box of 30 x 33g powder for crowns/veneers.
2. Box of 30 x 33g powder for inlays.
3. Bottle of 200ml liquid.

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
TAISEI DENTAL MFG. CO., LTD.