GC EUROPE N.V.**Head Office**

Interleuvenlaan 13
B - 3001 Leuven
Tel. +32.16.39.80.50
Fax. +32.16.40.02.14
info@gceurope.com
www.gceurope.com

GC GERMANY GmbH

Paul-Gerhardt-Allee 50
D - 81245 München
Tel. +49.89.89.66.74.0
Fax. +49.89.89.66.74.29
info@germany.gceurope.com
www.germany.gceurope.com

GC ITALIA S.r.l.

Via Calabria 1
I - 20098 San Giuliano Milanese
Tel. +39.02.98.28.20.68
Fax. +39.02.98.28.21.00
info@italy.gceurope.com
www.italy.gceurope.com

GC UNITED KINGDOM Ltd.

22-23, Coopers Court
Newport Pagnell
UK - Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info@uk.gceurope.com
www.uk.gceurope.com

GC FRANCE s.a.s

9 bis, Avenue du Bouton d'Or – BP 166
F - 94384 Bonneuil sur Marne Cedex
Tel. +33.1.49.80.37.91
Fax. +33.1.49.80.37.90
info@france.gceurope.com
www.france.gceurope.com

GC EUROPE N.V.**Austrian Office**

Tallak 124
A - 8103 Rein bei Graz
Tel. +43.3124.54020
Fax. +43.3124.54020.40
info@austria.gceurope.com
www.austria.gceurope.com

GC EUROPE N.V.**Benelux Office**

Edisonbaan 12
NL - 3439 MN Nieuwegein
Tel. +31.30.630.85.00
Fax. +31.30.605.59.86
info@benelux.gceurope.com
www.benelux.gceurope.com

GC EUROPE N.V.**East European Office**

Čazmanska 8
HR - 10000 Zagreb
Tel. +385.1.46.77.251
Tel. +385.1.46.78.474
Fax. +385.1.46.78.473
info@eeo.gceurope.com
www.eeo.gceurope.com

GC EUROPE N.V.**Finnish and Baltic States Office**

Vanha Hommaksentie 11B
FIN - 02430 Masala
Tel. & Fax. +358.9.221.82.59
info@finland.gceurope.com
www.finland.gceurope.com

GC EUROPE N.V.**Scandinavian Office**

Kungssporten 4 A
S - 427 50 Billdal
Tel. +46 31 939553
Fax. +46 31 914246
info@scandinavia.gceurope.com
www.scandinavia.gceurope.com

GC EUROPE N.V.**Scandinavian Office – Denmark**

Tværdiget 22
DK - 2730 Herlev
Tel. +45 44 53 53 82
Fax. +45 44 53 53 87
info@denmark.gceurope.com
www.denmark.gceurope.com

GC EUROPE N.V.**Swiss Office**

Wilerstrasse 3
CH - 9545 Wängi
Tel. +41.52.366.46.46
Fax. +41.52.366.46.26
info@switzerland.gceurope.com
www.switzerland.gceurope.com

GC America Inc.


3737 West 127th Street
USA - Alsip, Illinois 60803
Tel. +1.800.323.7063
Fax. +1.708.371.5103
sales@gcamerica.com
www.gcamerica.com

GC Asia Dental Pte. Ltd.

19 Loyang Way #06-27
SINGAPORE 508724
Tel. +65.6546.7588
Fax. +65.6546.7577
gcasia@mbox3.signet.com.sg
www.gcasia.info

GC FUJIVEST® PLATINUM

PHOSPHATE BONDED INVESTMENT (TYPE I)
PHOSPHATGEBUNDENE EINBETTMASSE (TYP I)
REVETEMENT A LIANT PHOSPHATE (TYPE I)
RIVESTIMENTO A LEGANTE FOSFATICO (TIPO I)
FOSFAATGEBONDEN INBEDMASSA (TYPE I)
REVESTIMIENTO DE BASE DE FOSFATO (TIPO I)
FOSFATBUNDEN INBÄDDNINGSMASSA (TYP I)
FOSFATBUNDEN INDSTØBNINGSMASSE (TYPE I)
FOSFAATTISIDONNAINEN VALUMASSA (TYYPPI I)

REVISED OCTOBER 2006

GC FUJIVEST® PLATINUM

GB	PHOSPHATE BONDED INVESTMENT	4
D	PHOSPHATGEBUNDENE EINBETTMASSE	7
FR	REVETEMENT A LIANT PHOSPHATE	10
IT	RIVESTIMENTO A LEGANTE FOSFATICO	13
NL	FOSFAATGEBONDEN INBEDMASSA	16
ES	REVESTIMIENTO DE BASE DE FOSFATO	19
SW	FOSFATBUNDEN INBÄDDNINGSMASSA	22
DK	FOSFATBUNDEN INDSTØBNINGSMASSE	25
FIN	FOSFAATTISIDONNAINEN VALUMASSA	28

Dear customer,

We thank you for purchasing a high quality GC product.

This Instruction For Use gives you a clear summary of each working stage and is based on the results of numerous lab tests carried out by our Research and Development Department.

If you have any further questions on the use of this product, please contact your local GC representative.

For more information on GC products please visit our website www.gc-europe.com.

GENERAL REMARK

All information included in this Instruction For Use is based on extensive testing and an extended range of casting trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.), different end-results may be obtained.

GC FUJIVEST® PLATINUM is a carbon-free phosphate bonded precision casting investment for quick heating as well as conventional step heating procedure.


Specially developed for precision castings of precious, semi-precious and Pd-base alloys for crown and bridgework.

GC FUJIVEST® PLATINUM offers you a number of advantages:

- Can be used in a quick-heating method or a progressive heating-up schedule
- Extremely smooth casting surfaces
- Perfect fit of the castings even for long span restorations like implant works
- Very easy devesting
- Easy mixing because of optimal wettability

PHYSICAL DATA (typical data)

ISO/DIN 9694.2: Measured at 23°C with 100% Expansion liquid	20' setting	120' setting
Heating up / Burn out procedure	Quick heating	Conventional step heating
Setting expansion	0.50%	0.90%
Thermal expansion	0.80%	0.60%
Total expansion (Linear)	1.30%	1.50%
Initial setting time	6'30"	6'30"
Compressive strength	5 MPa	5 MPa
Working time	5'	5'
Flow	13,5 cm	13,5 cm

INSTRUCTIONS FOR USE

1. Storage

Store powder and liquid at normal room temperature (23°C).

Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the castings.

- ① Protect the liquid against low temperatures!
Do not store below 5°C.
Once frozen, the liquid cannot be used anymore.
Be careful with winter deliveries.

2. Working temperature

Use at ± 23°C room temperature (19°C minimum).

- ① Higher working temperatures slightly reduce working time but improve the surface smoothness.
Lower temperatures prolong setting time and could result in surface roughness.

3. Preparations before investing

Use GC Multi Sep as a wax separator, which leaves no oily film on the die surface.

Use a 1 mm thick dry casting liner, GC New Casting Liner.

Seal the borders of the liner with a thin layer of vaseline.

- ① Open or wet liners can change the liquid content of the mixture and thus affect the expansion.
Be sure that the total inner surface of the metal ring is covered with casting liner.
Ensure that the liner is flush with the upper rim of the metal ring, so no direct contact of investment with the metal ring.
X3 ringsize = 1 layer of GC New Casting Liner
X6 ringsize = 2 layers of GC New Casting Liner
X9 ringsize = 2 layers of GC New Casting Liner

GC FUJIVEST® PLATINUM is a very fluid investment, which can be used without any wetting agents. Nevertheless, if a wetting agent is used, be sure to totally dry the surface prior to investing.

4. Ring sizes

Can be used with a ring and ringless technique.

- ① Possible to use with: X1, X3, X6, X9 for ring technique and X1, X3, X6 for ringless technique
- ① With ringsize X1 care should be taken. Due to smaller amount of material, meaning less exothermic reaction during setting, smaller castings could result. Temperatures below 23°C may result in rough surfaces.
Ringsize X3 gives the best results.

5. Powder/liquid ratio

Ringsize	Powder	Liquid
X1	60 g	13,2 ml
X3	150 g	33,0 ml
X6	300 g	66,0 ml
X9	420 g	92,4 ml

- ① Standard P/L ratio = 100 g / 22 ml
Exact powder/liquid measurement is necessary to obtain stable results.
Use adapted measure equipment such as electric balance, liquid measuring cup, pipette.

6. Expansion

Basic info on liquid concentration.

The used concentration of the expansion liquid depends on various parameters:

- storage temperature of the investment (powder and liquid)
- mixing time
- setting under pressure
- the amount of mixture

- ① Pure liquid gives the maximum setting expansion, meaning the largest casting.
Pure water gives the minimum setting expansion, meaning the smallest casting.
Use only distilled water to dilute.

Liquid dilution chart

Type Alloy	Ring-size	High Precious >70% Au		Semi Precious < 55% Au		Pd -Base		High Precious ceramic alloy		Semi Precious ceramic alloy		Pd-Base ceramic alloy		
		ml liquid	ml dist. water	ml liquid	ml dist. water	ml liquid	ml dist. water	ml liquid	ml dist. water	ml liquid	ml dist. water	ml liquid	ml dist. water	
Inlays/onlays Partial Crowns	Wax patern	x1	6	72	6,6	6,6	9	4,2						
		x3	15	18	15	18	23	10						
		x6	30	36	30	36	46	20						
		x9	41,4	51	41,4	51	64	28,4						
Crowns Bridges Primary Crowns	Wax Pattern	x1	6,6	6,6	7,2	6	10	3,2	9,9	3,3	9,2	4	9,2	4
		x3	16,5	16,5	18	15	25	8	25	8	23	10	23	10
		x6	33	33	36	30	49	17	50	16	46	20	46	20
	GC Pattern Resin	x9	46,2	46,2	50,4	42	69	23,4	69,4	23	64,4	28	64,4	28
		x1	8	5,2	8,6	4,6	10,6	2,6						
		x3	20	13	21	12	26	7						
Secondary Units	GC Pattern Resin	x6	40	26	43	23	53	13						
		x9	55,5	37	60	32,4	74	18,4						
		x1	6,4	6,8	7,2	6	9,2	4	9,2	4	8	5,2	8	5,2
		x3	16	17	18	15	23	10	23	10	20	13	20	13
		x6	32	34	34	32	46	20	46	20	40	26	40	26
		x9	44,8	47,6	47	45,5	64,4	28	64,4	28	56	36,4	56	36,4

NOTE

All information included in this Instruction For Use is based on extensive testing and a whole series of casting trials.

However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) different end-results may be obtained.

All fit tests have been carried out with the quick heating method (20 min. setting)

7. Mixing

1. Pre-mix powder and liquid thoroughly by hand with a spatula.

- ① Be sure all powder is wetted out by the liquid to give a uniform mixture.

2. Mix for 60 seconds under vacuum (420 rpm).

- ① Always use clean mixing bowl and check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles in the castings.

8. Working time

At least 4 minutes pouring time at room temperature (23°C).

- ① Working time / pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.

9. Investing

Invest under gentle vibration (low frequency).

- ① GC FUJIVEST® PLATINUM is very fluid, so strong vibration is not necessary and not advisable.

From the moment the ring is totally filled, **stop vibration immediately** and do not touch the investment until set.

10. Setting time

Leave to set for 20 min. from **start of mixing**.

- ① **Best results are obtained by putting immediately into a preheated furnace after 20 min.**

Scrape the top surface of the investment ring with a sharp knife.

- ① Diamond disks for trimmers are damaged by the quartz and cristobalite particles of the investment.

11. Heating Up – Burn Out Procedure

	Quick heating	Conventional Step heating
Insertion temperature	Pre-heated furnace at 750°C-800°C/1380°F-1470°F for Au-alloys 800-850°C/1470-1560°F* for Ceramic alloys	Room temperature
Step 1		Room temperature (23°C/73°F) to 260°C/500°F Heat rate 3°C/37°F min.
Step 2		Holding time at 260°C/500°F 60 - 90 Min
Step 3		Rise temperature from 260°C/500°F to 580°C/1076°F 5°C/41°F per Min
Step 4		Holding time at 580°C/1076°F 20 - 50 Min.
Step 5		Rise temperature from 580°C/1076°F to 750°C-800°C/1380°F-1470°F for Au-alloys Rise temperature from 580°C/1076°F to 800-850°C/1470-1560°F* for Ceramic alloys Heat rate 5°C/41°F per Min
Holding time	X 1 40 min. at end temperature	X 1 30 min. at end temperature
	X 3 50 min. at end temperature	X 3 40 min. at end temperature
	X 6 60 min. at end temperature	X 6 50 min. at end temperature
	X 9 90 min. at end temperature	X 9 60 min. at end temperature

Due to aggressive burning out, do not open the furnace during heating-up.

In case of vacuum pressure casting, raise the end temperature by 50°C/122°F.

* 850°C/1560°F: for casting-on technique of prefabricated metal parts.

When several investment rings are put into the furnace at the same time, prolong the heating period.

750°C/1380°F is the minimum pre-heating temperature

12. Casting

Cast in the usual manner: centrifugal casting, vacuum-pressure casting, etc.

Take care with the positioning of the casting ring in the casting device.

Cast as soon as possible after removing the ring from the furnace.

13. Cooling

Cool down the casting as slowly as possible.

e.g. insert into cold furnace and close the door.

Place the ring upside down for cooling down.

Notes

- Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC FUJIVEST® PLATINUM will delay the setting time of gypsum products.
- It is recommended that mixing bowls are stored in water between uses.
- Store powder & liquid at room temperature (+/- 23°C / 73°F).

Safety recommendations and hazard warnings

- Investment materials contain quartz.

Do not inhale dust!

Risk of lung damage (Silicosis, lung cancer).

Advise : wear protection mask Type FFP 2 – EN 149:2001.

Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.

Rinse the empty investment material bag with water before disposal.

- Remove dust from your working place only when it is wet.
- To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
- When sandblasting the cast object, always use a fine-dust filter extraction system.

Packages

Powder: 6 kg box (60 g pack x 100)
6 kg box (90 g pack x 67)
6 kg box (150 g pack x 40)

Liquid: 900 ml bottle

Expiry date: 2 years from the manufacturing date

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