# GC FUJIVEST® PREMIUM

PHOSPHATE BONDED INVESTMENT (TYPE I)  
PHOSPHATGEBUNDENE EINBETTMASSE (TYP I)  
REVETEMENT A LIANT PHOSPHATE (TYPE I)  
RIVESTIMENTO A LEGANTE FOSFATICO (TIPO I)  
FOSFAATGEBONDEN INBEDMASSA (TYPE I)  
REVESTIMIENTO DE BASE DE FOSFATO (TIPO I)  
FOSFATBUNDEN INBÄDDNINGSMASSA (TYP I)  
FOSFATBUNDEN INDSTØBNINGSMASSE (TYPE I)  
FOSFAATTISIDONNAINEN VALUMASSA (TYYPPI I)

# GC FUJIVEST® PREMIUM

GB	PHOSPHATE BONDED INVESTMENT	4
D	PHOSPHATGEBUNDENE EINBETTMASSE	7
FR	REVETEMENT A LIANT PHOSPHATE	10
IT	RIVESTIMENTO A LEGANTE FOSFATICO	13
NL	FOSFAATGEBONDEN INBEDMASSA	16
ES	REVESTIMIENTO DE BASE DE FOSFATO	19
SW	FOSFATBUNDEN INBÄDDNINGSMASSA	22
DK	FOSFATBUNDEN INDSTØBNINGSMASSE	25
FIN	FOSFAATTISIDONNAINEN VALUMASSA	28

Dear customer,

We thank you for purchasing a high quality GC product.

This Instruction For Use gives you a clear summary of each working stage and is based on the results of numerous lab tests carried out by our Research and Development Department.

If you have any further questions on the use of this product, please contact your local GC representative.

For more information on GC products please visit our website [www.gceurope.com](http://www.gceurope.com).

### General remark

All information included in this Instruction For Use is based on extensive testing and an extended range of casting trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.), different end-results may be obtained.

GC FUJIVEST® PREMIUM is a carbon-free phosphate bonded investment for Crowns & Bridges, for use in both quick and step-heating procedure.


Suitable for all dental alloys with special attention to non-precious alloys.

GC FUJIVEST® PREMIUM offers you a number of advantages:

- Suitable for all dental alloys with special attention to non-precious castings
- Ringless techniques possible in both heating-up procedures (quick & slow)
- Total expansion up to 3.5%
- Perfect expansion control by liquid dilution
- High fluidity
- Robust & strong investment but easy devesting

### PHYSICAL DATA (typical data)

ISO/DIN 9694.2; Measured at 23°C with 100% Expansion liquid	20' setting	120' setting
Setting expansion	2.00%	2.50%
Thermal expansion	1.00%	1.00%
Total expansion (Linear)	3.00%	3.50%
Initial setting time	9 min.	9 min.
Compressive strength	7,5 MPa	7,5 MPa
Working time	6 min.	6 min.
Flow	13,4 cm	13,4 cm

## INSTRUCTIONS FOR USE

### 1. Storage

Store powder and liquid at normal room temperature (23°C). Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the castings.

- ① Protect the liquid against low temperatures!  
Do not store below 5°C.  
Once frozen, the liquid can not be used anymore.  
Be careful with winter deliveries.

### 2. Working temperature

Use at ± 23°C room temperature (19°C minimum).

- ① Higher working temperatures slightly reduce working time but improve the surface smoothness.  
Lower temperatures prolong setting time and could result in surface roughness.

### 3. Preparations before investing

Use GC Multi Sep as a wax separator, which leaves no oily film on the die surface.

Use a 1 mm thick dry casting liner, GC New Casting Liner.

Seal the borders of the liner with a thin layer of vaseline.

- ① Open or wet liners can change the liquid content of the mixture and thus affect the expansion.  
Be sure that the total inner surface of the metal ring is covered with casting liner.  
Ensure that the liner is flush with the upper rim of the metal ring, so no direct contact of investment with the metal ring.  
X3 ringsize = 1 layer of GC New Casting Liner  
X6 ringsize = 2 layers of GC New Casting Liner  
X9 ringsize = 2 layers of GC New Casting Liner
- ① GC FUJIVEST® PREMIUM is a very fluid investment, which can be used without any wetting agents.  
Nevertheless, if a wetting agent is used, be sure to totally dry the wax surface prior to investing.

### 4. Ring sizes

Quick-heating & step-heating method.

Quick heating and slow heating method can be used with ring and ringless techniques.

- ① Possible to use with: X1, X3, X6, X9 for ring technique and X1, X3, X6 for ringless technique
- ① With ringsize X1 care should be taken. Due to smaller amount of material, meaning less exothermic reaction during setting, smaller castings could result. Temperatures below 23°C may result in rough surfaces.  
Ringsize X3 gives the best results.

## 5. Powder/liquid ratio

Ringsize	Powder	Liquid
X1	60 g	13,2 ml
X3	150 g	33,0 ml
X6	300 g	66,0 ml
X9	420 g	92,4 ml

- ① Standard P/L ratio = 100 g / 22 ml  
Exact powder/liquid measurement is necessary to obtain stable results.  
Use adapted measure equipment such as electric balance, liquid measuring cup, pipette.

## 6. Expansion

Basic info on liquid concentration:

The used concentration of the expansion liquid depends on various parameters:

- storage temperature of the investment (powder and liquid)
  - mixing time
  - setting under pressure
  - the amount of mixture
- ① Pure liquid gives the maximum setting expansion, meaning the largest casting.  
Pure water gives the minimum setting expansion, meaning the smallest casting.  
Use only distilled water to dilute

## Liquid dilution chart

Type of Alloy	Liquid / water ratio %	Ring size : Liquid / water				
		X1 / 60g	90g	X3 / 150g	X6 / 300g	X9 / 420g
High Precious >70% Au	50/50	6.6ml/6.6ml	9.9ml/9.9ml	16.5ml/16.5ml	33ml/33ml	42.2ml/42.2ml
Semi Precious < 55% Au	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Pd-base	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Non-precious alloy	NiCr 75/25	10ml/3.2ml	15ml/4.8ml	25ml/8ml	50ml/16ml	69.4ml/23ml
	CoCr 100%	13.2ml	19.8ml	33ml	66ml	92.4ml
High Precious ceramic alloy	55/45	7.3ml/5.9ml	11.8ml/8ml	18ml/15ml	36ml/30ml	51ml/41.4ml
Semi Precious ceramic alloy	55/45	7.3ml/5.9ml	11.8ml/8ml	18ml/15ml	36ml/30ml	51ml/41.4ml
Pd-base ceramic alloy	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Non-precious ceramic alloy	NiCr 75/25	10ml/3.2ml	15ml/4.8ml	25ml/8ml	50ml/16ml	69.4ml/23ml
	CoCr 100%	13.2ml	19.8ml	25ml/8ml	66ml	92.4ml

### NOTE

All information included in this Instruction For Use is based on extensive testing and a whole series of casting trials.

However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) different end-results may be obtained.

All fit tests have been carried out with the quick heating method (20 min. setting).

120 min. setting will result in a somewhat wider casting.

## 7. Mixing

1. Pre-mix powder and liquid thoroughly by hand with a spatula.

- ① Be sure all powder is wetted out by the liquid to give a uniform mixture.

2. Place the mixture under vacuum for 15 seconds without mixing.

- ① To reach optimal degassing.

3. Mix for 60 seconds under vacuum (320-420 rpm).

- ① Always use clean mixing bowl and check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles in the castings.

## 8. Working time

At least 6 minutes pouring time at room temperature (23°C).

- ① Working time / pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.

## 9. Investing

Invest under gentle vibration (low frequency).

- ① GC FUJIVEST® PREMIUM is very fluid, so strong vibration is not necessary and not advisable. From the moment the ring is totally filled, **stop vibration immediately** and do not touch the investment until set.

## 10. Setting time

1. Leave to set for **20 min from start of mixing**.

- ① Longer setting is possible.  
Best results are obtained by putting immediately into a preheated furnace after 20 min.

2. Scrape the top surface of the investment ring with a sharp knife.

- ① Diamond disks for trimmers are damaged by the quartz and cristobalite particles of the investment.

## 11. Heating up - Burn Out Procedure

	Quick heating	Conventional Step heating
<b>Insertion temperature</b>	Pre-heated furnace at 700-750°C/1290-1380°F for Au - alloys 800-850°C/1470-1560°F* for Ceramic alloys 900°C/1650°F for Non-precious alloys	Room temperature
<b>Step 1</b>		Room temperature (23°C) to 260°C/500°F Heat rate 3°C/37°F per min.
<b>Step 2</b>		Holding time at 260°C/500°F 40 - 90 min.
<b>Step 3</b>		Rise temperature from 260°C/500°F to 580°C/1076°F Heat rate 6°C/43°F per min.
<b>Step 4</b>		Holding time at 580°C/1076°F 20 - 50 min.
<b>Step 5</b>		Rise temperature from 580°C/1076°F to 750°C/1380°F for Au-Alloys Rise temperature from 580°C/1076°F to 800-850°C/1470-1560°F* for Ceramic Alloys Rise temperature from 580°C/1076°F to 900°C/1650°F for Non-precious Alloys Heat rate 9°C/48°F per min.
<b>Holding time</b>	X 1 40 min. at end temperature	X 1 30 min. at end temperature
	X 3 50 min. at end temperature	X 3 40 min. at end temperature
	X 6 60 min. at end temperature	X 6 50 min. at end temperature
	X 9 90 min. at end temperature	X 9 60 min. at end temperature

Due to aggressive burning out, do not open the furnace during heating-up.

In case of vacuum pressure casting, raise the end temperature by 50°C/122°F.

\*850°C/1560°F: for casting-on technique of pre-fabricated metal parts.

When several investment rings are put into the furnace at the same time, prolong the heating period.

## 12. Casting

Cast in the usual manner: centrifugal casting, vacuum-pressure casting, etc.

Take care with the positioning of the casting ring in the casting device.

Cast as soon as possible after removing the ring from the furnace.

## 13. Cooling

Cool down the casting as slowly as possible.

e.g. insert into cold furnace and close the door.

Place the ring upside down for cooling down.

## Notes

1. Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC FUJIVEST® PREMIUM will delay the setting time of gypsum products.
2. It is recommended that mixing bowls are stored in water between uses.
3. Store powder & liquid at room temperature (+/- 23°C / 73°F).

## Safety recommendations and hazard warnings

1. Investment materials contain quartz.  
Do not inhale dust!  
Risk of lung damage (Silicosis, lung cancer).  
Advise: wear protection mask Type FFP 2 – EN 149:2001.  
Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.  
Rinse the empty investment material bag with water before disposal.
2. Remove dust from your working place only when it is wet.
3. To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
4. When sandblasting the cast object, always use a fine-dust filter extraction system.

## Packages

Powder: 6 kg box (60 g pack x 100)  
6 kg box (90 g pack x 67)  
6 kg box (150 g pack x 40)

Liquid: 900 ml bottle

Expiry date : 2 years from the manufacturing date

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### Tableau ratio liquide / eau distillée

Alliage	Ratio standard liquide / eau distillée %	Taille du cylindre: liquide / eau				
		X1 / 60g	90g	X3 / 150g	X6 / 300g	X9 / 420g
Précieux >70% Au	50/50	6.6ml/6.6ml	9.9ml/9.9ml	16.5ml/16.5ml	33ml/33ml	42.2ml/42.2ml
Semi-précieux < 55% Au	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Alliage Pd-base	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Alliage non-précieux	NiCr 75/25	10ml/3.2ml	15ml/4.8ml	25ml/8ml	50ml/16ml	69.4ml/23ml
	CoCr 100%	13.2ml	19.8ml	33ml	66ml	92.4ml
Alliage céramique précieux	55/45	7.3ml/5.9ml	11.8ml/8ml	18ml/15ml	36ml/30ml	51ml/41.4ml
Alliage céramique semi-précieux	55/45	7.3ml/5.9ml	11.8ml/8ml	18ml/15ml	36ml/30ml	51ml/41.4ml
Alliage céramique à base de Pd	60/40	8ml/5.2ml	11.8ml/8ml	20ml/13ml	40ml/26ml	55.4ml/37ml
Alliage céramique non-précieux	NiCr 75/25	10ml/3.2ml	15ml/4.8ml	25ml/8ml	50ml/16ml	69.4ml/23ml
	CoCr 100%	13.2ml	19.8ml	25ml/8ml	66ml	92.4ml

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**GC America Inc.**
**GC Asia Dental Pte. Ltd.**