Initial^{AL}

BASIC SET

ADVANCED SET

IN VIVO / IN SITU SET

Technical Manual/**Technische Arbeitsanleitung**
Manuel Technique/**Manuale tecnico**

GC Initial AL

GC

REVISED MAY 2008

*True-to-nature ...
Naturidentisch ...*

*Proche du naturel ...
Estetica naturale ...*


By means of a simple standardized build up technique, lifelike aesthetics can be reproduced ...

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Standardized build up procedure Standardisierte Schichttechnik

The aluminumoxide (Al_2O_3) framework.
GC Initial AL can be used as veneering ceramics for all aluminumoxide (Al_2O_3) frame materials.

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GC Initial AL Liner (Lin-1 – Lin-6) is applied thinly on the framework.

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
GC Initial AL Liner equally applied on the coping.

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
Standardized build up procedure Standardisierte Schichttechnik


INside (IN) material of the chosen colour is applied on the cervical as well as on the approximal part (see colour combination charts).

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Fluo Dentin (FD-91 – FD-93) is applied in the incisal area for optimal light transmission, as well as thinly over the entire labial surface. (see colour combination chart)

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
Dentin material in its correct shade is applied.

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
Standardized build up procedure Standardisierte Schichttechnik

The complete Dentin layering in the desired inner structure.

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
Application of a thin layer of Clear Fluorescence (CL-F) over the entire Dentin layering.

This layer is 0,2 mm in depth and is a copy of human dentition.

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Standardized build up procedure Standardisierte Schichttechnik

Enamel (E) or Enamel Opal (EOP) is applied on top of the Clear Fluorescence (CL-F) layer up to the final shape.

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Final labial build up. To create the milky white incisal contour, Enamel Opal (EO-15) should be used.

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
Standardized build up procedure Standardisierte Schichttechnik


The palatal surface of the restoration is covered with the high chroma INside (IN) material. (see colour combination chart)

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The marginal ridges are lined with the same Dentin mixture that has been used on the labial side.

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Final oral build up with the same Enamel mixture that has been used on the labial side.

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Standardized build up procedure Standardisierte Schichttechnik

The restoration is fired following the firing parameters of the 1st dentin firing (see firing instructions). The fired ceramic surface shows a slight gloss.

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The proximal areas are supplemented with the same material used for the Dentin / Enamel layering. The restoration is finally adjusted using Enamel Transparent materials. Cervical Translucent (CT) can be used as an alternative for the Dentin / Enamel layering.

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Completed "GC Initial" restoration after glazing. Application of individualized characterizations with INvivo / INSitu Stains is possible.

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Multi-chromatic build up procedure

Polychrome Schichttechnik

The aluminumoxide (Al_2O_3) framework.

GC Initial AL can be used as veneering ceramics for all aluminumoxide (Al_2O_3) frame materials.

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The GC Initial AL Liner is applied equally on the coping.

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Multi-chromatic build up procedure Polychrome Schichttechnik


INside (IN) of the chosen colour is applied on the cervical part as well as on the approximal part (see colour combination chart).

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
Fluo Dentin (FD-91 – FD-93) is applied in the incisal area for optimal light transmission and thinly over the entire labial surface (see colour combination chart).

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Multi-chromatic build up procedure

Polychrome Schichttechnik

A thin layer of INside (IN) Primary Dentin covers the entire labial frame surface (see colour combination chart).

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Standard Dentin powder is mixed with 10% of INside powder (see colour combination chart) and very thinly applied (about 50% of the usual Dentin layering thickness).

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This way of working allows a true-to-nature build up procedure of the Transpa / Enamel layering. Dentin is applied in its final form.

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
Multi-chromatic build up procedure Polychrome Schichttechnik


The build up of an incisal shield with Enamel (E/EOP) and / or Transpa (T) mixtures.

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By means of a „Cut - Back“, the incisal part is reduced up to the internal Fluo-Dentin layering.

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
The reduced area is reconstructed with Fluo Dentin (FD-91) in the desired Mamelon structure. This is a very important point for colour synchronization of different base structures placed next to each another.

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Multi-chromatic build up procedure

Polychrome Schichttechnik

The level of colour intensity of the Mamelon structure is defined by INside colours or INvivo stains (see colour combination chart).

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This base is now covered with Clear Fluorescence (CL-F), the layer is 0,2 mm in depth and is a copy of human dentition.

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
Multi-chromatic build up procedure Polychrome Schichttechnik

The Enamel / Transpa layering, where special incisal and transparent mixtures are alternatively placed next to another (e.g. E58, TO, EOP2, CL-F).

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The final contouring is made out of a mixture of Enamel Opalescence (EOP) and standard Enamel (E). Depending on the age of the patient, the translucent level of the Enamel layer can be modified with TN, TO & CL-F. To create the milky white incisal contour, EO-15 should be used.

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The palatal surface of the restoration is covered with a high chroma INside colour (see colour combination chart).

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Multi-chromatic build up procedure

Polychrome Schichttechnik

The oral marginal ridges are lined with standard Dentin material.

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Final oral build up and functional contouring is done with the same Enamel mixture that has been used on the labial side and finished with Enamel Occlusal (EO-15).

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The restoration is fired following the parameters of the 1st Dentin firing (see firing instructions). The fired ceramic surface shows a slight gloss.

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
Multi-chromatic build up procedure Polychrome Schichttechnik

*Labial view after the first Dentin firing.
The very low shrinkage makes the following correction firing more easy.*

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The completed "GC Initial" restoration after glazing. Application of individualized characterizations with INvivo / INsitu stains is possible.

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Firing Instructions Brennanleitung

Nature is the ultimate, GC Initial an equal partner.

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GC Initial – Colour Combination Chart

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GC Initial INside – Field of application and reference chart

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GC Initial INSide – Field of application and reference chart

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- Approximal / Cervical / Oral colours
Approximal / Cervical / Oral Farben

Shade Group Farbgruppe	Shades Farben	“INSide” reference “INSide” Referenz
A-shades A-Farben	A1, A2, A3	IN-42 Terracotta
		IN-44 Sand
	A3.5, A4	IN-42 Terracotta
		IN-45 Havanna
		IN-46 Brasil
IN-49 Maracuja		
B-shades B-Farben	B1, B2	IN-43 Sun
		IN-47 Sienna
	B3, B4	IN-47 Sienna
		IN-48 Kurkuma
IN-50 Curry		
C-shades C-Farben	C1, C2	IN-51 Olive
	C3, C4	IN-51 Olive
		IN-45 Havanna
D-shades D-Farben	D2, D3, D4	IN-44 Sand
		IN-51 Olive

High Fluorescence Dentin – Field of application and reference chart

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GC Initial Shoulder Porcelain – Reference chart

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The Shoulder Opaques can be used as a first basic layering at the border frame / tooth shoulder. The percentages used in this list can vary depending on the thickness of the total ceramic shoulder.

BASIC RULE: The thinner the total ceramic shoulder, the more high chroma shoulder powders should be used. Furthermore, all shoulder powders can be mixed with the INside powders to intensify or to individualize the shoulder shades until a maximum of 20% INside powders of this total mixture.

In this case the final temperature of shoulder firing should be reduced by 10°C to 20°C, depending on the amount of INside powders mixed with the shoulder porcelain.

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Firing charts for various furnaces

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
AUSTROMAT 3001

1st Opaque Firing Washbrand	C580 T0 T120.L9 V9 T60 T080.C950 V0 T60 C0 L0 T2 C580
1st and 2nd Shoulder Firing 1. und 2. Schulterbrand	C600 T180 T120.L9 V9 T60 T080.C950 V0 T60 C0 L0 T2 C600
1st Dentin Firing 1. Dentinbrand	C580 T180 T180.L9 V9 T60 T055.C910 V0 T60 C0 L0 T2 C580
2nd Dentin Firing 2. Dentinbrand	C580 T180 T180.L9 V9 T60 T055.C900 V0 T60 C0 L0 T2 C580
Glaze Firing Glanzbrand	C600 T120.L9 T60 T055.C910 C0 L0 T2 C580
Glaze firing glaze powder Glanzbrand Glasur	C480 T120.L9 T60 T045.C880 T60 C0 L0 T2 C580
Correction Powder Firing Korrekturmasse-Brand	C450 T120 T120.L9 V9 T60 T045.C790 V0 T60 C0 L0 T2 C580

We wish to strongly emphasise that all tests with these furnaces were in-house studies only and no long-term tests were performed. Therefore we would appreciate your understanding that this info might only be taken as a guidance.

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Vita[®] is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.
GC NORDIC AB