GC Initial™ LiSi Press LITHIUM DISILICATE GLASS CERAMICS

For use only by a dental professional in the recommended indications.

- INDICATIONS FOR USE
• Occlusal veneers
• Thin veneers
• Veneers
• Onlays
• Crowns in the anterior and posterior region
• 3-unit bridges in the anterior region
• 3-unit bridges in the premolar region up to the second premolar as the terminal unit
• Crown or splinted crown top on top of an implant abutment
• Severe reduction of occlusal vertical dimension
• Parafuncions (e.g. suffering from bruxism)


Table with 3 columns: Trans, Processing technique, Indications. Rows include HT, MT, LT, LIHQ, MO.

CONTRAINDICATIONS
Not suitable for patients with:
• Severely reduced residual dentition
• Severely reduced occlusal vertical dimension
• Parafuncions (e.g. suffering from bruxism)


PREPARATION DESIGN

When designing restorations, the following preparation guidelines and minimum dimensions should be maintained:


- Preparation guidelines for all-ceramic restorations
• Prepare margins with deep chamfer or rounded shoulder (fig. 1)
• Do not prepare undercut
• Ensure that the cavity walls form an angle of 6° with the long axis of the tooth
• All internal edges and angles should be rounded.


Occlusal veneer
• Reduce the incisal crown third - in the occlusal area by approx. 1.0 mm.


Thin veneer
• If possible, the preparation should be located in the enamel
• Reduction in the cervical and/or labial area by 0.3 mm, and the incisal edge by 0.4 mm.


Inlays
• A preparation depth of at least 1.0 mm and an isthmus width of at least 1.0 mm must be observed in the fissure area.

Crowns in the anterior region
• Width of the shoulder/chamfer should be at least 1.0 mm
• Reduce the incisal crown third - in the occlusal area by approx. 1.5 mm

Crowns in the posterior region
• Reduce the anatomical shape and observe the stipulated minimum
• Width of the shoulder/chamfer should be at least 1.0 mm

3-unit bridges
• The maximum acceptable pontic width is different in the anterior and posterior region. The pontic width is determined on the unprepared tooth.

WAXING UP
The investment is carried out with GC Initial LiSi Press Framework (fig. 2)
Fabricate the working model and apply the a press (fig. 2)

INVESTING
GC Initial LiSi Press
1. The investing is carried out with GC Initial LiSi Press, a specially developed phosphate bonded investment, offering easy removal of the reaction layer.

PREHEATING
1. After setting, remove the mould base former and ring base former.
2. Simultaneously, place the investment ring on the investment ring.

PRESSING
1. Remove the investment ring from the preheating furnace immediately after completion of the preheat step.

PRESSING PROGRAM
1. DEKEMA AUSTROMAT 64 (PANAMAT Press)

Table with 4 columns: Ring size, Start temperature, Hold time, Pressure. Rows include HT, MT, LT, LIHQ, MO.

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