## METAL PRIMER II

PRIMER FOR BONDING RESIN TO METAL

METALPRIMER II contains a special monomer (MEPS : thiophosphoric methacrylate) which promotes bonding of acrylic resins and composites to all types of metal.

#### RECOMMENDED USE

Bonding of acrylic and composite resins to all types of metal as set out below in the DIRECTIONS FOR USE.

#### NOTES AND CONTRAINDICATIONS

- For use only by dental professionals for the applications recommended.
- Contains methacrylate. Not to be used by operators, or on patients, with a history of allergy to methacrylates. In cases of allergic response, seek medical attention.
- In case of contact with skin or eyes, immediately flush with water and seek medical attention.
- Flammable. Avoid high temperatures, direct sunlight and sources of ignition. Use in a well ventilated area.
- Volatile. Seal container immediately after use.
- Assign the brush provided (No. 7) for use ONLY with METALPRIMER II.

#### DIRECTIONS FOR USE

#### GENERAL

In all cases

- Sandblast metal surfaces with alumina (aluminium oxide), and blow clean with air gun.

Note : Avoid contamination with oil, separator or moisture from fingers, etc.
- Dispense 2-3 drops of METALPRIMER II into a dappen dish or similar container. Using the brush provided, apply a thin layer to the bonding surface. Allow to dry.

Note : A second coat may be applied if necessary.
- Apply resin to metal as indicated below.

Note : Apply resin to metal immediately after the surface is prepared.

#### RECOMMENDED APPLICATIONS

- Resin cement to metal for metal inlays, crowns and bridges
  - Prepare metal surface as above and bond restoration with resin cement.
- Opaque resin to metal for composite resin veneers and crowns
  - Prepare surface as above and apply opaque (GC GRADIA) in usual manner.

Cure as recommended by the manufacturer.
- Denture base to magnetic attachment
  - Prepare surface as above and bond attachment in usual manner using self curing acrylic resin (GC UNIFAST Trad).
- Denture base resin to metal base, clasp or attachment
  - Prepare bonding surface with alumina sandblast and blow clean with air gun.
  - Form a wax denture.
  - After flushing out the wax, clean the bonding surface using a steam cleaner. Remove wax completely using solvent.
  - Apply METALPRIMER II.
  - Pack and cure denture base in usual manner.
- In the mouth repair of damaged crown and bridge veneers
  - Place rubber dam.
  - Clean metal surface with carbide bur and prepare as above.

Special note : Be sure to use rubber dam and avoid direct contact with oral tissue. In case of contact, immediately flush with water and seek medical attention.
- Repair of denture with metal base, clasp or attachment
  - Roughen metal bonding surface using carborundum point.
  - Prepare bonding surface as above.
  - Build up using self-cure resin in the usual manner.
- Bonding resin to metal orthodontic appliances
  - Prepare bonding surfaces as above.
  - Apply orthodontic resin to the appliance in the usual manner.

#### PACKAGE

1 bottle (5 mL) METALPRIMER II
1 brush (No. 7)

#### STORAGE

Store in a cool place (4-25°C) out of direct sunlight. Shelf-life - 2 years from date of manufacture.

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.