**Prior to use, carefully read the instructions for use**

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For use only by a dental professional in the recommended indications.

### DESCRIPTION:

NATURE-CRYL HI-20 ET is a premium high impact denture base material that is recommended for use in the fabrication of full dentures and partial dentures. The NATURE-CRYL HI-20 ET formulation includes a unique methyl methacrylate copolymer and a specially formulated liquid that yields a denture base system that exhibits a number of exciting features and benefits.

Features	Benefits
High Impact Strength – Exceeds ISO1567 requirement for High Impact denture base material	Reduced risk of breakage
20 Minute Boiling Water Cure	Speed and convenience
Extended Working Time at Room Temperature	Improves workflow, productivity, and compatible with popular injection system
Extended Working Time when refrigerated	Improves workflow, productivity, and eliminates waste
Cadmium Free Pigments	Eliminate cadmium toxicity concerns

Type 1- Class 1 (Heat-polymerizable polymers - Powder and Liquid) per ISO1567

### RECOMMENDED INDICATIONS:

For use in the fabrication of full and partial dentures using conventional heat cured dough pack and or injection processing technique.

### WARNING:

Contains methyl methacrylate. In susceptible individuals, sensitization to the product cannot be excluded. Use of the product should be discontinued if allergic reactions are observed.

### CAUTION:

- 1.Flammable Liquid. Do not use near open flames or sources of ignition. Close the container immediately after use.
- 2.Avoid contact with skin and eyes. In case of contact with eyes, flush immediately with water and seek medical attention.
- 3.Avoid prolonged contact with the liquid and unpolymerized dough and inhalation of monomer (liquid) fumes/vapors.
- 4.Refer to the MSDS for additional health and safety information.
- 5.Do not mix with other resins.

### STORAGE:

NATURE-CRYL HI-20 ET Powder and Liquid should be stored in a cool (60-80°F, 15-27°C), dark, and dry place away from heat and sources of ignition or from sunlight.

### SHELF LIFE:

NATURE-CRYL HI-20ET Powder: 3 years from the date of manufacture  
NATURE-CRYL HI-20ET Liquid: 3 years from the date of manufacture.

### INSTRUCTIONS FOR USE:

#### Preparatory Work:

Wax-up the appliance using conventional materials and techniques. Invest the case using gypsum, i.e. COE-CAL™ Type III Laboratory Stone. Boil out in the usual manner by immersing flask in boiling water for about 4 minutes to soften the wax so that it can easily be removed from the mold. Flush and then rinse the mold thoroughly with clean boiling water. The mold surfaces must be wax-free in order for the tinfoil substitute to form a protective film. *Do not use solvents such as acetone or chloroform.* When the mold stops steaming, apply COE-SEP™ tinfoil substitute.

#### Mixing:

Equipment required: mixing jar with cover, spatula, and a balance or graduated cylinder. Be sure all measuring and mixing utensils are clean, and dry.

**Proportions by weight: 10g liquid (monomer) to 20g powder (polymer)**

**Proportions by volume: 10mL liquid (monomer) to 27mL powder (polymer)**

- Carefully measure the desired amount of liquid and pour it into the mixing jar.
- *Warning: The liquid is highly flammable. Keep it away from open flames or sources of ignition.*
- Shake powder bottle to redistribute fibers and then measure the desired amount of powder.
- Stir the polymer powder into the liquid monomer until powder is wet (about 30 seconds).

#### Packing:

- If using a conventional 2 hour or 9 hour curing technique, the mix will be ready to start packing after it reaches the dough stage or when it is no longer tacky and doesn't stick to the sides of the mixing jar. Under ordinary conditions of 23°C/73°F this will be between 7-10 minutes.

- If using the 20 minute cure technique, the mix will be ready to start packing after it reaches the snap stage or when it repeatedly breaks cleanly with a snapping sound. *It is very important that the dough reaches this snap stage.* Under ordinary conditions of 23°C/73°F this will be between 12-15 minutes.

#### Note:

*In warmer temperatures less time will be required for the mix to reach desired consistency, in cooler temperatures more time may be required.*

One mix can be used for approximately 1.0-1.5 hours if the closed/sealed container is kept at room temperature of 23°C/73°F. The mix can be used for approximately 24 hours if the closed/sealed container is placed in the refrigerator immediately, and if the resin is packed within a reasonable time after removal from the refrigerator. Allow the material to return to room temperature before packing.

The flask should be at room temperature and the resin at the dough or snap stage depending on cure technique. Remove enough dough from the mixing jar to slightly over-fill the mold and re-cover the mixing jar to protect remaining contents from the effects of evaporation. The trial packing technique is recommended. Pack the dough into the mold and cover with a dry sheet of polyethylene. Slowly close the flask halves together until they are just shy of metal-to-metal contact. Open the flask and examine the dough in the mold immediately. If it has a velvety gloss appearance and the land area is covered with flash, the contents of the mold have been adequately compressed (otherwise add more material and repeat). Remove the flash, separate with a polyethylene sheet or other separator, and close the flask halves attempting to obtain metal-to-metal closure of the flask with moderate pressure. Open the flask again and remove any remaining flash and separator sheets. Obtain final metal-to-metal closure of the flask, using no more pressure than will be maintained by the curing compress.

#### Curing Schedule:

- **20 Minute Cure\*** - Immerse flask in boiling water for 20 minutes.
- **2 Hour Cure** - Immerse the flask in a water bath at 162-165°F(72-74°C) to cure for 1-1/2 hours: then transfer to boiling water bath for 1/2 hour.
- **9 Hour Cure** - 9 hours at 165°F(74°C)

*\*Note: The 20 min. boiling water curing method is not recommended for cases exceeding 8.5mm in thickness.*

#### Cooling:

Bench cool to room temperature. An alternate and acceptable method is to bench cool the flask for 30 minutes at room temperature and then accomplish final cooling under running tap water at 56-91°F(13-33°C) for 15 minutes.

#### Finishing and Polishing:

Finish and polish using standard technique.

#### Note:

1. Do not generate excessive frictional heat when finishing and polishing. High temperatures can cause distortion of the denture base.
2. Important- Keep the finished denture in water until it is inserted in the patients mouth. Instruct the patient to keep denture in water during any period when it is not being worn to prevent it from drying out. This will ensure better conformation to the mouth structures.

#### Repairs:

It is recommended that repairs be made using auto polymerizing (cold cure) resins such as NATURE-CRYL™ POUR in order to keep dimensional changes in the repaired denture to a minimum. Be sure that plaster surfaces underlying the repair area are protected by a tinfoil substitute. For maximum strength and density, repairs made with cold cure resins should be cured under pressure in an ACRI-DENSE™ Curing Unit according to the instructions furnished with the unit. If it is necessary to heat cure repairs, the curing temperature should never exceed 165°F(74°C). DO NOT BOIL.

#### SHADES:

##### **Powder available in 7 Shades**

(COE Light Fibered Plus, COE-LOR™ Mild, COE-LOR™ Moderate, COE-LOR™ Heavy, Original, Light Reddish Pink, Light)

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MADE IN THE U.S.A.
PRINTED IN U.S.A.