Prior to use, carefully read the instructions for use

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Self-Cured Denture Base & Repair Resin

For use only by a dental professional in the recommended indications.

DESCRIPTION:

NATURE-CRYL POUR is a versatile denture base material that may be used to either construct full and partial dentures using the pour resin technique or for hard relines and repair.

NATURE-CRYL POUR yields a denture base system that exhibits a number of exciting features and benefits.

Features	Benefits
Room Temperature Acrylic	Ideal for Repairs
Fluid Resin Technique	Labor saving and simple to process
Dimensionally Accurate	Dentures retain original fit, fewer Relines and Rebases
Color Stable	Denture will not change or fade in color
Cadmium Free Pigments	Eliminate cadmium toxicity concerns

Type 2 - Class 2 (Autopolymerizable polymers - Powder and Liquid) per ISO1567

RECOMMENDED INDICATIONS:

For use in the fabrication of full and partial dentures using conventional pour resin processing technique and for hard reline and repair of existing dentures. *Not recommended for chairside denture relining.*

WARNING:

Contains methyl methacrylate. In susceptible individuals, sensitization to the product cannot be excluded. (Use of the product should be discontinued if allergic reactions are observed.)

CAUTION:

- Flammable Liquid. Do not use near open flames or sources of ignition. Close container immediately after use.
- Avoid contact with the liquid on skin and in your eyes. In case of contact with eyes, flush immediately with water and seek medical attention.
- Avoid prolonged contact with the liquid and unpolymerized dough and inhalation of monomer (liquid) fumes/vapors.
- Do not mix with other resins.
- Refer to the MSDS for additional health and safety information.

STORAGE:

NATURE-CRYL POUR Powder and Liquid should be stored in a cool (60-80°F, 15-27°C), dark, and dry place away from heat, sources of ignition, and direct sunlight.

SHELF LIFE:

NATURE-CRYL POUR Powder: 3 years from the date of manufacture.
NATURE-CRYL POUR Liquid: 3 years from the date of manufacture.

I. INSTRUCTIONS FOR USE: FULL or PARTIAL DENTURE FABRICATION

Recommended Materials & Equipment:

- Investing materials: Alginate (COE ALGINATE™) or Hydrocolloid duplicating material (DUPLI-COE-LOID™)
- Flask: (POUR-n-CURE FLASK)
- Curing Unit: ACRI-DENSE™ CURING UNIT or equivalent.
- Disposable plastic mixing cups, spatula, balance or graduated cylinder.

Preparatory Work:

- Wax-up the full or partial denture using conventional materials and techniques.
- Trim the sides of the cast to provide a slight convergence towards the land areas to permit easy withdrawal. Land areas should be completely flat and 1/8" wide.
- Soak the waxed up cast in warm water for 10 minutes to eliminate air from the cast.
- Place the waxed-up cast in the POUR-n-CURE FLASK at an angle so that one sprue will be slightly longer. Be sure the cast is flat on the base of the flask.
- Investing:
 - When using **DUPLI-COE-LOID** hydrocolloid:
 - Follow the instructions for melt down of the duplicating material and cool it down to the dispensing temperature of 140°F.
 - Slowly pour the DUPLI-COE-LOID into the assembled flask between the cast and the flask wall, avoid pouring the duplicating material directly on the cast. Continue pouring to cover the teeth and completely fill the flask and reservoir
 - Carefully place the flask in cool circulating water for 30 minutes. *The depth of the water should cover approximately the bottom 3/4 height of the flask.*
 - When using **COE ALGINATE**:
 - Mix 16 oz. of water to one 8 oz. cup of powder. Use cold water to obtain longer working time.
 - Pour the mixed alginate over the wax-up cast.
 - Immediately place the flask in the ACRI-DENSE Pneumatic Curing Unit. Apply 20 psig air pressure.
 - Remove the flask from the ACRI-DENSE unit one minute after the alginate residue in the mixing bowl has set.
- After the investing material has gelled or set, remove the reservoir and trim the material flush with the bottom surface of the cover plate.
- Turn the POUR-n-CURE FLASK upside down, remove the insert to expose the base of the cast and gently tease the waxed-up cast out of the investment material.
- Remove the wax from the cast and clean the mold with hot (130°F – 55°C) water and dry it. Painting the surface of the mold with COE-SEP tinfoil substitute is recommended to seal the model surface and thereby make a smoother denture.
- For complete dentures use a large sprue cutter to make holes of proper size for rapid filling and venting of the mold. Cases are sprued and vented to the highest portion of the wax up. This prevents void caused by trapped air. For the saddle areas of partial dentures use a small sprue cutter.
- Remove the wax and teeth from the cast and thoroughly clean the teeth and the cast with clean boiling water. Replace dry teeth in the mold. *Note: To improve retention of plastic denture teeth mechanical retention grooves may be cut into the ridge lap of the denture teeth.*
- Paint the surface of the cast with a thin coat of COE-SEP that has been diluted with equal parts of water. And then allow to dry.
- Carefully reinsert the cast into the mold, seating it firmly on the land areas.
- Reassemble the POUR-n-CURE FLASK so the sprue holes are facing up.

Mixing, Pouring, and Curing:

Warning: The liquid is highly flammable. Keep it away from open flames or sources of ignition. Close container immediately after use. Be sure all measuring and mixing utensils are clean, and dry.

*	Powder (P)	Liquid (L)	P/L Ratio
By Weight	20 g	14 g	1.43
By Volume	28 cc	15 mL	1.87

*Quantities correspond to material required for a typical denture. Use more or less material for larger or smaller dentures.

- Shake powder container to redistribute fibers prior to measuring the desired amount.
- Powder and liquid must be measured accurately. If too much liquid is used, excessive shrinkage will occur. If too much powder is used, the viscosity of the mix will be too heavy to pour, resulting in voids or incomplete denture.
- The liquid can be refrigerated if additional pouring time is desired.

- Carefully measure the recommended proportions of powder and liquid.
- Pour the liquid into disposable mixing plastic cup.
- Slowly add the powder to the liquid. Stir thoroughly with a metal spatula for ten seconds. Pour the mixture between two cups at least twice to fully incorporate powder and liquid.
- One minute from the start of mixing pour the NATURE-CRYL POUR fluid resin into the sprue hole. Pour in a continuous stream to avoid entrapping air and continue pouring until the vent hole is filled.
- Place the flask in the ACRI-DENSE Pneumatic Curing Unit with the sprues upright.
- Add enough warm water (100°F – 38°C) to the ACRI-DENSE Unit to submerge the flask to 3/4 of its total height.
- Engage and completely lock the cover and pressurize the unit to 20 psig.
- Cure the case under 20-psig pressure for 30 minutes.
- Remove the flask from the ACRI-DENSE Unit and de-flask the denture.
- Go to Finishing and Polishing section of the instructions

II. INSTRUCTIONS FOR USE: “LABORATORY” RELINE USING NATURE-CRYL POUR DENTURE BASE

Caution: NATURE-CRYL POUR is not intended for use in chairside relines and repair. Processing must be done on and accurate COE-CAL model not in the patient's mouth.

- Pour a COE-CAL stone cast in the impression.
- Take a COE-CAL index of the teeth, allow to set.

- Mount the cast on an articulator (anterior pin stop preferred) or a reline jig.
- Separate the denture from the cast. Completely remove all impression material.
- Remove 1/16” to 1/8” (1.5 mm to 3 mm) around the entire peripheral border of the denture. With a pear shaped bur grind the surfaces around the ridge. Labial, buccal, and palatal vault areas to permit the even flow of reline material.
- Apply COE-SEP to the cast. Also paint the labial and buccal surfaces of the new denture with COE-SEP to prevent the reline material from adhering.

Warning: The liquid is highly flammable. Keep it away from open flames or sources of ignition. Close container immediately after use.

Powder / Liquid Mixing Proportions

	Powder	Liquid	Powder/Liquid Ratio
By Weight	10 g	7 g	1.43
By Volume	14 cc	7.5 mL	1.87

Be sure all measuring and mixing utensils are clean, and dry. Shake powder container to redistribute fibers prior to measuring the desired amount. Powder and liquid must be measured accurately. If too much liquid is used, excessive shrinkage will occur. If too much powder is used, the viscosity of the mix will be too heavy to pour, resulting in voids or incomplete denture. The liquid can be refrigerated if additional working time is desired.

- Using a brush or eyedropper moisten the tissue surface of the denture with NATURE-CRYL POUR Liquid.
- Mix 10g (14 cc) NATURE-CRYL POUR Powder and 7 g (7.5 mL) NATURE-CRYL POUR Liquid. Stir for 15 seconds. In approximately 4 minutes the mixture will thicken and flow slowly.
- Spread the mixture on the tissue side of the denture in the same manner you apply impression paste.
- Carefully place and seat the denture on the cast. Close the articulator into the index, making sure of the pin stop closure. Hold the articulator in place with a large rubber band. If a reline jig is used seat the denture in place and carefully close the jig using the wing nuts.
- Curing: The use of the ACRI-DENSE Curing Unit is recommended to facilitate a stronger and denser repair with less porosity. The case & jig should be immersed in warm water in the ACRI-DENSE Curing Unit and cured for 30 minutes under 20 psig pressure.
- Remove the flask from the ACRI-DENSE Unit and de-flask the denture.
- Go to Finishing and Polishing section of the instructions.

Finishing and Polishing:

Finish and polish using standard techniques.

Note:

- Do not generate excessive frictional heat when finishing and polishing. High temperatures can cause distortion of the denture base.
- Important- Keep the finished denture in water until it is inserted in the patient's mouth. Instruct the patient to keep denture in water during any period when it is not being worn to prevent it from drying out. This will ensure better conformation to the mouth structures.

III. INSTRUCTIONS FOR USE: “DENTURE REPAIR USING NATURE-CRYL POUR DENTURE BASE

Caution: NATURE-CRYL POUR is not intended for use in chairside relines and repair. Processing must be done on an accurate COE-CAL model not in the patient's mouth.

Powder / Liquid Mixing Proportions

For POUR MIXTURE: Powder / Liquid Mixing Proportions

	Powder	Liquid	Powder/Liquid Ratio
By Weight	10 g	7 g	1.43
By Volume	14 cc	7.5 mL	1.87

FOR PACKING MIXTURE: Powder / Liquid Mixing Proportions (Wait until mixtures reaches dough state)

	Powder	Liquid	Powder/Liquid Ratio
By Weight	10 g	5 g	2.0
By Volume	14 cc	5.4 mL	12.59

FOR DUSTING: Carefully moisten the fracture with liquid and apply a layer of powder. Repeat until the desired thickness is achieved. Be sure all measuring and mixing utensils are clean, and dry. Shake powder container to redistribute fibers prior to measuring the desired amount.

Powder and liquid must be measured accurately. If too much liquid is used, excessive shrinkage will occur. If too much powder is used, the viscosity of the mix will be too heavy to pour, resulting in voids or incomplete denture. The liquid can be refrigerated if additional pouring time is desired.

REPLACING A TOOTH:

- Make a box preparation on the palatal or lingual side of the denture. The preparation should be straight and even.
- Select a replacement tooth of the proper mold and shade. Hold it in place with sticky wax.
- Apply a lubricant to the labial surface of the denture and make a COE-CAL matrix of the labial surface. Remove the matrix.
- Mix NATURE-CRYL POUR powder and liquid and pack it into the box preparation, or dust in alternative layers of powder and liquid. Build up slightly over the margins of the box preparation.
- Cover the resin with a piece of cellophane and apply pressure to both the matrix and resin until the resin is polymerized (set).

REPAIRING A FRACTURE:

- Brush-Bead Build-up or Brush on technique – This technique use a camel's hair brush and two dappen dishes, one for the NATURE-CRYL POUR powder and one for the NATURE-CRYL POUR Liquid. The brush is first dipped in the liquid and then into the powder (only the tip of the brush should be dipped in the powder and then extracted very quickly) the result being a saturated bead of polymer on the brush tip. Build up the repair area slightly higher than the level of the denture. Case is now ready for curing placement in the ACRI-DENSE Curing Unit or it can be placed in warm water until cured
- Dispenser (Salt & Pepper) Technique: This technique uses a camel's hair brush and a small polyethylene dispenser. Paint the area with NATURE-CRYL POUR Liquid and then apply NATURE-CRYL POUR Powder by depressing or squeezing the flexible bottle slowly to control the amount of powder. Add liquid with the brush until the powder is completely saturated. Build up the repair area slightly higher than the level of the denture. Case is now ready for curing placement in the ACRI-DENSE Curing Unit or it can be placed in warm water until cured.

SHADES:

Powder Available in 3 Shades
(COE Light Fibered +, COE-LOR™ Mild, COE-LOR™ Moderate, COE-LOR™ Heavy, Original, Light Reddish Pink, Light)

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