G-CEM CAPSULE

SELF ADHESIVE LUTING CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges (max. 2 abutments)
- Cementation of metal, ceramic, fiber posts, and cast post and cores.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	0,35 / 0,15
Standard Mixing Time (sec.) (4000 rpm at 23°C)	10"
Working Time (min., sec.) (at 23°C)	2'30"

Test conditions: Relative humidity (50+/-10%).

A. CEMENTATION TECHNIQUE for ceramics, indirect composites and metal based restorations

1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
- Clean the prepared tooth with pumice and water.
- Rinse thoroughly with water. Dry by gently blowing with oil free air. Prepared surfaces should appear dry.

- Do not use substances like H₂O₂ (hydrogen peroxide) and/or disinfectants as these may lower the bond strength.
 - Special care should also be taken to remove so called imaging powders very thoroughly with pumice and water.

2. RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

3. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To remove excess cement, plunge the applicator until it flush with the main body (Fig. 2).

c) Immediately place the capsule into a GC CAPSULE APPLIER and click the lever once (Fig. 3). The capsule is now activated.

Note: The capsule should be activated just before mixing and used immediately.

- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000RPM) (Fig. 4). More details on mixing times can be found in following chart.

Mixing time chart

Mixer Speed	Mixing Time
Standard mixer (approx 4,000rpm)	10 seconds
High-speed mixer (over 4,500rpm)	
Rotational mixer *No centrifuge option	

Note: Refer to the respective manufacturer's instructions for use.

4. DISPENSING

- Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.

- Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- Extrude directly into the restoration (Fig. 6).

Note: Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.

- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will show working time. Refrigeration of G-CEM CAPSULE is recommended to extend the working time.

- Maintain moderate pressure.

6. EXCESS CEMENT REMOVAL

- To remove excess cement, light cure using a visible light curing unit (Halogen, GC G-Light, LED) for 2-4 seconds (Fig. 7) or allow the cement to self-cure until the excess cement feels rubbery.

- Remove excess cement with a scaler (Fig. 8).

- Light cure each surface for 10 seconds (single surface from occlusal:10 seconds, any other surfaces:10 seconds each) or allow the cement to self-cure for 4 minutes after seating of the restoration.

NOTE:

G-CEM is a dual cure cement. To insure proper curing involving cementation of the following restorations, G-CEM must be used in self cure mode for 4 minutes.

- Metal, metal coping restorations
- Metal-free dark and opaque restorations
- Translucent metal-free restorations greater than 2mm in thickness
- Finish and polish margins with appropriate instruments.

B. CEMENTATION TECHNIQUE for metal, fiberglass and ceramic posts

1. TOOTH PREPARATION

Prepare the post space according to manufacturer's instructions. The post space should be cleaned, rinsed and dried thoroughly using paper points. The post space should be dry.

NOTE:

Sodium hypochlorite solution is recommended to use for the chemical cleaning of the post space.

- Do not use H₂O₂ (hydrogen peroxide) to chemically clean the post space as it may lower the bond strength.

2. POST PREPARATION

Make sure that the post is pretreated and handled according to the manufacturer's instructions.

3. CAPSULE ACTIVATION AND MIXING

Refer to A.3.

4. DISPENSING & CEMENTATION

NOTE:

Refrigerate G-CEM Capsules to extend working time.

- Immediately remove the mixed capsule from the mixer and attach a GC Elongation Tip on the capsule nozzle firmly until it clicks. GC Elongation Tip can only be attached to G-CEM CAPSULE with a retention ring on the nozzle (Fig. 9). Do not attach the tip on other capsules. Immediately place the capsule with the elongation tip into a GC CAPSULE APPLIER.

- Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body. If necessary, the capsule nozzle can be bended.

- Place the tip deep in the post space and start extruding the cement within 45 seconds from start of mix. Keep the tip immersed in the cement whilst dispensing to avoid trapping air. (Fig. 10).

NOTE:

Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.

- Inset the post immediately into the post space within 1 minute after cement application.

Note:

1) The oral temperature accelerates the setting reaction of the cement. 2)DO NOT use lenuito spirales to place the cement in the post space as this will shorten working time.

- Maintain moderate pressure.

- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

5. EXCESS CEMENT REMOVAL

- Continue to maintain moderate pressure while removing excess cement with appropriate instruments.

- Light cure all surfaces for 20 seconds or allow the cement to self cure for 4 minutes after seating of the post (Fig. 11).

STORAGE

Store the capsules in the original aluminium foil in a cool and dark place (4-25°C) (39.2-77.0°F).

(Shelf life : 2 years from date of manufacture)

SHADES

Vita® (Vita® shade), A03 (opaque), B01 (opaque) and Translucent Vita® are registered trademarks of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE

- Single shade package : 50 capsules.
- Assorted package : 50 capsules (20 each of A2 and Translucent, 5 each of A03 and B01).

Average contents per capsule : 0.35g powder and 0.15g (0.14mL) liquid. Minimum net volume of mixed cement per capsule : 0.17mL.

- GC CAPSULE APPLIER (1 piece) is available separately.
- GC Capsule Elongation Tip: 50 pieces

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.

- In case of contact with eyes, flush immediately with water and seek medical attention.


- If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks.

- Do not use G-CEM in combination with eugenol containing materials as eugenol may hinder G-CEM from setting or bonding properly.

- Use a protective light shield or similar protective eye wear during light curing. This product is not indicated for filling or core build-up.

- When wiping elongation tip with a cotton soaked with alcohol, prevent the alcohol from entering the inside of the tip. Elongation tip is for single use only and cannot be autoclaved or chemically.

Last revised : 02/2011

MANUFACTURED BY GC CORPORATION 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8685, Japan

EUROPEAN REPRESENTATIVE GC CORPORATION 33, Rue de la Woluwe, 1200 Brussels, Belgium

GC AMERICA INC. 3737 West 127th Street, Minneapolis, MN 55438, U.S.A.

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