Prior to use, carefully read the instructions for use. EN

G-aenial™ Flo X

LIGHT-CURED RADIOPAQUE FLOWABLE COMPOSITE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Blocking out undercuts
2. Repair of (in) direct aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel
4. Sealing hypersensitive areas
5. Fissure sealant on sensitive Bereiche
6. Direct restorative for Class II, III, IV, I and V cavities including small cavities which require high flowability.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

Operation

- 1) Hold the syringe upright and remove the cap by turning counterclockwise.
- 2) Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
- 3) Place the light protective cap until use.

Note:
Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnant material.

2. Shade Selection

Clean the tooth with pumice and water. Shade selection should be made prior to isolation. Select the appropriate shade from 8 shades of A1, A2, A3, A3.5, A4, A4.2, AO2, AO3 and CV. A shade guide such as the Vita® classical shade guide. *Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

3. Cavity Preparation Prepare cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate prepared tooth from contamination with saliva, blood or sulcus fluid.

Note:
For pulp capping, use calcium hydroxide.

4. Bonding Treatment For bonding G-aenial Flo X to enamel and/or dentine, use a light-cured bonding system such as G-Premio BOND, G-aenial Bond, G-BOND, GC Fuji Bond LC, UniFil® Bond. Follow manufacturer's instructions.

5. Placement of G-aenial Flo X

- 1) Remove cap from the dispensing tip.

Note:
Prior to syringe material into the cavity, make sure to check the following. 1. Dispensing tip is securely attached to the syringe. 2. Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the tip pointing upwards until material reaches the opening of the tip (Fig. 2).

- 2) Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to syringe material.

Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.

- a. If plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
- b. Minize exposure to ambient light. Ambient light can shorten the manipulation time.
- c. After use, immediately remove and dispose the dispensing tip and tightly close the syringe with the cap.

Clinical Hint

In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

6. Light Curing

Light cure the G-aenial Flo X using a light curing unit. Keep light guide as close as possible using the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure.

Shade	Irradiation time	10 sec. (High power LED) (1200mW/cm ²)	20 sec. (Halogen/LED) (700 mW/cm ²)
A1, A2, A3		2,0mm	
A3.5, A4, CV, AO2, AO3		1,5mm	

- Note:
1. Material should be placed and light cured in layers. For maximum layer thickness, please consult above table.
 2. Lower light intensity may cause insufficient curing or discoloration of the material.

7. Finishing and Polishing Finish and polish using standard techniques.

SHADES

A1, A2, A3, A3.5, A4, AO2, AO3 and CV

STORAGE

Recommended for optimal performance, store in a cool and dark place (-25°C / 39,2-77,0°F) away from high temperatures or direct sunlight.

PACKAGES

1. Syringe 3.6g (2,0mL), 20 dispensing tip III (plastic type), 1 light protective cap
- Dispensing tip package:
 2. 30 dispensing tip II (plastic type), 2 light protective caps
 3. 30 dispensing tip III (needle type), 2 light protective caps

CAUTION


1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Take care to avoid ingestion of the material.
4. Do not mix with other similar products.
5. Dispensing tips are for single use only. To prevent cross contamination between patients, do not reuse the dispensing tip. The dispensing tip cannot be sterilized in an autoclave or chemically.
6. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.
7. Wear protective eye glasses during light curing.
8. When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
9. Avoid getting material on clothing.
10. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
11. Do not use G-aenial Flo X in combination with eugenol containing materials as eugenol may hinder G-aenial Flo X from setting.
12. Do not wipe the dispensing tip with excessive alcohol as this may cause contamination of the paste.
13. Dispose of all wastes according to local regulations.
14. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
15. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

http://www.gceurope.com
or for The Americas
http://www.gcamerica.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use remove and label for deterioration. Discard device if damaged.
DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised: 06/2018

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

EU : **GC EUROPE N.V.**
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

GC SOUTH AMERICA
Rua Heitor de Sá, Santana - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-51
RESP: TEG: Mayara de Sá Ribeiro - CRO/SP 105.982

GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05 Singapore 528729
TEL: +65 6548 7588

GC AUSTRALASIA DENTAL PTE. LTD.
1753 Botolph Claydon, Mackay, NSW 2019, Australia
TEL: +61 2 9301 8200

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