Prior to use, carefully read the instructions for use.

EN

G-ænial Universal Flo

UNIVERSAL LIGHT-CURED RADIOPAQUE FLOWABLE COMPOSITE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Direct restorative for Class I, IV, III, II and V cavities.
- Fissure sealant.
- Sealing hypersensitive areas.
- Repair of (in) direct aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel
- Blocking out undercuts
- Liner or base

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

Operation

- Preparations
 - Hold the syringe upright and remove the cap by turning counterclockwise.
 - Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
 - Place the light protective cover until use.

Note : Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnant material.
- Shade selection

Clean the tooth with pumice and water. Shade selection should be made prior to isolation.

Select the appropriate shade from 16 shades of A1, A2, A3, A3.5, A4, B1, B2, B3, C3, AO1, AO2, AO3, BW, CV, JE, AE, A, B, C shades are based on the Vita® classical shade guide. *Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
- Cavity Preparation

Prepare cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate prepared tooth from contamination with saliva, blood or sulcus fluid.


- Preparation

For pulp capping, use calcium hydroxide.

Bonding treatment
For bonding G-ænial Universal Flo to enamel and/or dentine, use a light-cure bonding system such as G-Premio BOND, G-ænial Bond, G-BOND, GC Fuji Bond LC, UniFil® Bond. Follow manufacturer's instructions.
- Placement of G-ænial Universal Flo
 - Remove cover from the dispensing tip.

Note: Prior to syringe material into the cavity, make sure to check the following:
 - Dispensing tip is securely attached to the syringe.
 - Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the tip pointing upwards until material reaches the opening of the tip (Fig. 2).

- Placing the syringe
 - Dispensing tip is securely attached to the syringe.
 - Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the tip pointing upwards until material reaches the opening of the tip (Fig. 2).
- Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to syringe material. Alternatively, dispense material onto a mixing pad and transfer to the cavity using a suitable instrument.
 - if plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
 - Minimize exposure to ambient light. Ambient light can shorten the manipulation time.
 - After use, immediately remove and dispose the dispensing tip and tightly close the syringe with the cap.

- if plugged, remove the dispensing tip and extrude a small amount of material directly from the syringe.
- Minimize exposure to ambient light. Ambient light can shorten the manipulation time.
- After use, immediately remove and dispose the dispensing tip and tightly close the syringe with the cap.

Clinical Hint
In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provide a smooth surface over the restoration.

- Light Curing
Light cure the G-ænial Universal Flo using a light curing unit. Keep light curing as close as possible to the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure.

Irradiation time	9 sec. (Plasma arc) (2000mW/cm²)
Shade	10 sec. (GC G-Light) (1200mW/cm²)
	20 sec. (Halogen/LED) (700 mW/cm²)
A1, A2, A3, B1, B2, JE, AE	2.0mm
A3.5, A4, B3, C3, BW, CV, AO1, AO2, AO3	1.5mm

- Note :
 - Material should be placed and light cured in layers. For maximum layer thickness, please consult above table.
 - Lower light intensity may cause insufficient curing or discoloration of the material.

- Finishing and Polishing

Finish and polish using standard techniques.

SHADES

A1, A2, A3, A3.5, A4, B1, B2, B3, C3, AO1, AO2, AO3, BW, CV, JE, AE

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F) away from high temperatures or direct sunlight. (Shelf life : 3 years from date of manufacture)

PACKAGES

- Syringe 3.4g (2.0mL), 20 dispensing tips (plastic type), 1 light protective covers
- Dispensing tip package
- 30 dispensing tips (plastic type), 2 light protective covers
- 30 dispensing tips (needle type), 2 light protective covers

CAUTION

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingestion of the material.
- Do not mix with other similar products.
- Dispensing tips are for single use only. To prevent cross contamination between patients, do not reuse the dispensing tip. The dispensing tip cannot be sterilized in an autoclave or chemically.
- Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.
- Wear protective eye glasses during light curing.
- When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
- Avoid getting material on clothing.
- In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
- Do not use G-ænial Universal Flo in combination with eugenol containing materials as eugenol may hinder G-ænial Universal Flo from setting.
- Do not wipe the dispensing tip with excessive alcohol as this may cause contamination of the paste.
- Dispose of all wastes according to local regulations.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com or for The Americas
http://www.gcamerica.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised: 05/2017

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

RESPECTUABLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY
GC CORPORATION
1-1 Harasumi-cho, Ichikawa-shi, Tokyo 174-8565, Japan

EU : GC EUROPE N.V.
Heersdijk 1, Heersdijk-Leuven 1340, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00
TEL: +1 708-597-0900
www.gcamerica.com

ASIA : GC ASIA DENTAL PTE. LTD.
11, Riverside Court, #03-06, Singapore 528729
TEL: +65 6546 7588

AUSTRALASIA : GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Bankstown NSW, NSW 2010, Australia
TEL: +61 2 9301 8200

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