GC Fuji® Automix LC

RADIOPAQUE LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE

For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

- Class III and V restorations
- Restoration of primary teeth
- Core Build-up
- Cases where radiopacity is required
- Base material for Class I and Class II cavities using a sandwich laminate technique.

CONTRAINDICATIONS

- Direct pulp capping
- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

PHYSICAL PROPERTIES

Paste/Paste Ratio (g/g)	1.6 / 1.0
Working Time (min., sec.)	1'50"
Light Curing Time* (sec.)	20"
Depth of Cure (mm)	2

*Light curing required.
Test conditions: Temperature (23 +/-1°C) (73 +/-2°F), Relative humidity (50 +/-10%)

IDENTIFICATION OF PARTS (Fig.1)

- A. Paste Pak Cartridge SL
1. Cover 2.Inner plunger 3.Cartridge tip 4. Sleeve
B. Paste Pak Dispenser II
5. Piston 6.Sliding block 7. Piston release 8. Lever
C. GC Fuji Mixing Tip SL with intraoral tip
9. Intraoral tip

DIRECTIONS FOR USE

1. PASTE PAK CARTRIDGE SL LOADING

Make sure that the piston of Paste Pak Dispenser II is completely retracted into the dispenser (Fig.2). To retract the piston, push the piston release and at the same time pull the sliding block (on the upper side of the dispenser) with fingers. Load the cartridge into the dispenser. Ensure the ▲ mark on the cartridge lines up with the ▲ mark on the dispenser as illustrated (Fig.3). Then turn the cartridge fully to the right position. After the cartridge loading, push the piston release forward until it stops (Fig.4). Do not push the piston on the upper side.

NOTE:

To remove the cartridge, retract the piston fully and rotate the cartridge back to the marked position to release from the dispenser.

2. CAVITY PREPARATION

- Prepare the cavity using standard techniques. Extensive mechanical retention is not necessary. For direct pulp capping, use calcium hydroxide. For a non cavitated Class V, clean the lesion with pumice and water, rinse thoroughly, and dry gently with an oil free air syringe
- Apply CAVITY CONDITIONER (10 seconds) or DENTIN CONDITIONER (20 seconds) to the prepared surfaces to remove the smear layer.

NOTE:

- Do not use phosphoric acid.
- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an oil-free air syringe. DO NOT DESICCATE.

3A. DISPENSING IN AUTOMIX MODE

- Remove the cartridge cover.
- Before every use, check the two openings of the cartridge tip and confirm the pastes are at the same level to ensure even flow from the cartridge.

Note: Before dispensing for the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.

- Attach the GC Fuji Mixing Tip SL on the cartridge. Push the mixing tip completely until the end. (Fig. 5-1) Hold the sides of sleeve and slide it forward to lock the mixing tip (Fig.5-2). The direction of the intraoral tip can be adjusted by rotating it in the desired direction.

NOTE:

- Make sure the GC Fuji Mixing Tip SL is securely attached to the cartridge. If not, remove the mixing tip and reattach.
- If the material has been refrigerated, leave the cartridge at room temperature for 10 minutes.

- Before dispensing the material in the cavity, bleed a small amount of the mixed paste on a pad.
- SLOWLY depress the lever to dispense the material into the prepared cavity.

- After dispensing, release the lever. Slide the sleeve backward and remove the mixing tip (Fig. 6). Clean the cartridge tip and always put the cover back on.

NOTE:

- Do not leave the used mixing tip on the cartridge for a long period of time.
- The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the upper side of the dispenser. When only a small amount of paste remains in the cartridge, the paste may not be dispensed from the mixing tip.

(3) To avoid damage, do not drop the dispenser with the cartridge and mixing tip.

3B. DISPENSING IN HANDMIX MODE

- Remove the cartridge cover.
- Before every use, check the two openings of the cartridge tip and confirm the pastes are at the same level to ensure even flow from the cartridge.

NOTE:

- Before dispensing for the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.
- Depress the lever to dispense the required amounts of paste onto the mixing tip.

- While depressing the lever, level and cut off the extruded material from the cartridge tips by pushing the cartridge tips on the mixing pad at the right angle.

- When the lever is released after dispensing, residual pastes at the cartridge tips are retracted into the cartridge. Clean the cartridge tip and replace the cartridge cover.

NOTE:

- The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the upper side of the dispenser.

- After dispensing, spread the two pastes thinly on the mixing pad using a plastic spatula or a mixing tip, with tapping strokes, for 10 seconds. Be careful not to incorporate air bubbles.

NOTE:

- Larger amounts of paste require a longer mixing time (15 seconds).

4. APPLICATION

- Place a matrix system if required.
- Remove surface moisture but DO NOT DESICCATE.
- Extrude the material directly into the preparation from the mixing tip or using an appropriate instrument in handmix mode. For deep cavities, use layering technique with increments of 2mm. The working time is 1 minute 50 seconds from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time.
- Form the contour.

- Light-cure* for 20 seconds using a visible light curing device (LED/Halogen >500mW/cm²). Place light source as closely as possible to the cement surface.

*Light curing required.

- Finish under water spray using superfine diamond bur, silicone point and polishing strips. Surface can be polished following the routine procedure using a composite polishing system.
- Apply a final coat of G-COAT PLUS, GC Fuji VARNISH or GC Fuji COAT LC. Refer to respective instructions for use.

STORAGE

Recommended for optimal performance, store at room temperature away from heat, moisture and direct sunlight (4-25°C) (39.2-77.0°F).

SHADE

- A1, A2, A3 based on Vita® shades
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE

- 1 x Paste Pak Cartridge SL (14.9g / 7.2mL)
- 20 x GC Fuji Mixing Tip SL with intraoral tip
- 1 x Paste Pak Dispenser II

CAUTION

- In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.

- Do not use GC Fuji Automix LC in combination with eugenol containing materials as eugenol may hinder GC Fuji Automix LC from setting or bonding properly.
- Take care to avoid ingestion of the material.
- Do not mix with other materials.

- GC Fuji Mixing Tip SL is for single use only. To prevent cross contamination between patients, do not reuse the mixing tip. The mixing tip cannot be sterilized in an autoclave or chemically.
- Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.

- Wear protective eye glasses during light curing.
- When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
- Avoid getting material on clothing.
- Do not wipe the GC Fuji Mixing Tip SL with excessive alcohol as this may cause contamination of the paste.
- Dispose of all wastes according to local regulations.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <http://www.gcamerica.com>. They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

CLEANING / STERILIZATION for Paste Pak Dispenser II

For cleaning and sterilization of Paste Pak Dispenser II, refer to the respective instructions for use.

Last revised: 07/2018

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Tabatah-ku, Tokyo 174-8585, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Tabatah-ku, Tokyo 174-8585, Japan
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-589-0500
www.gcamerica.com
GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP, BRASIL
RFP: 02029-091 - TEL: +55-11-2925-0961
CNPJ: 08.279.999/0001-61
RESP. TEC: Mayara de Santos Ribeiro - CRO/SP 105.982

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