# GC FujiCEM GC FujiCEM Automix

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

**RECOMMENDED INDICATIONS**

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

**CONTRAINDICATIONS**

- Pulp capping
- In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

**IDENTIFICATION OF PARTS (Fig. 1)**

A. Paste Pak Cartridge

- Cover 2 Inner plunger 3. Cartridge tip

B. Paste Pak Dispenser

- Piston 5. Adjusting guide 6. Lever 7. Sliding block 8. Piston release

C. GC FujiCEM Mixing Tip

- Release tab 10. Grip 11. Guide pin

**DIRECTIONS FOR USE**

Paste / Paste ratio (g/g)	2,0 / 1,0
Mixing time (sec.)	10"-15"
Working time (min., sec.) (at 23°C)	2'15"
(from start of mixing)	2'15"
Final finishing commencing time (min., sec.) (after seating the restoration)	4'30"

**PASTE PAK CARTRIDGE LOADING**

- Make sure that the piston is completely retracted into the dispenser (Fig. 2).  
Note:  
To retract the piston, push the piston release and at the same time pull the sliding block (on the underside of the dispenser) with fingers.
- Load the cartridge into the Paste Pak Dispenser (Fig. 3). Ensure the ▲ mark on the cartridge lines up with the ▲ mark on the Paste Pak Dispenser as illustrated. Then turn the cartridge fully to the right position.
- After cartridge loading, push the piston release forward until it stops (Fig. 4).  
a) Push the piston release forward.  
b) Do not push the upper side.  
Note:  
To remove the cartridge, retract the piston fully and turn the cartridge back to the marked positions.

**CEMENTATION TECHNIQUE – AUTOMIXING WITH GC FujiCEM Mixing Tip**

- TOOTH PREPARATION**
  - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
  - Clean the prepared tooth with pumice and water.  
Note:  
Use GC Fuji Plus Conditioner on the prepared tooth surface for 20 seconds for optimal adhesion.
  - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE.  
Prepared surfaces should appear moist (glistening).
- RESTORATION PREPARATION**  
Make sure that the restoration is pretreated and handled according to the manufacturers instructions.
- DISPENSING AND AUTOMIXING**
  - Remove the cartridge cover.
  - Before every use, check the two openings of the cartridge tip and confirm the pastes are at the same level to ensure even flow from the cartridge.  
Note:  
Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.
  - Attach the GC FujiCEM Mixing Tip correctly on the cartridge by aligning the guide pin on the GC FujiCEM Mixing Tip with the hole on the cartridge, simultaneously pushing the tip until it clicks on the cartridge (Fig. 5-1).  
Note:  
1) Make sure GC FujiCEM Mixing Tip is securely attached to the cartridge. If not, remove the GC FujiCEM Mixing Tip and reattach (Fig. 5-2).  
2) If the material has been refrigerated, leave the cartridge in room temperature for 10 minutes.
  - Set the adjusting guide to the front of the lever (Fig. 6).  
Note:  
Before dispensing the material into the appliances, bleed a small amount of the mixed paste on a pad.
  - Slowly depress the lever to dispense the material into the appliance.
  - After dispensing, release the lever. To remove the GC FujiCEM Mixing Tip, push the GC FujiCEM Mixing Tip forward while depressing the release tab with one hand, and move the front part of the GC FujiCEM Mixing Tip from side to side with the other hand (Fig. 7). Clean the cartridge tip and replace the cartridge cover.  
Note:  
1) Do not leave the used GC FujiCEM Mixing Tip on the cartridge for a long period of time.  
2) The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser. As the cartridge cylinders are translucent, the residual amount can also be confirmed by the position of the inner plunger seen through the cylinders. When only a small amount of paste remains, the amount of paste may not be dispensed from the GC FujiCEM Mixing Tip.

**4. CEMENTATION**

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 15 seconds from start of mixing at 23°C (74°F). Higher temperature will shorten working time. Maintain moderate pressure.
- Start removing excess cement with a scaler when the excess cement feels rubbery.
- Finishing can be started 4 minutes 30 seconds after seating the restoration.

**CEMENTATION TECHNIQUE – HANDMIXING**

- TOOTH PREPARATION**
  - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
  - Clean the prepared tooth with pumice and water.  
Note:  
Use GC Fuji Plus Conditioner on the prepared tooth surface for 20 seconds for optimal adhesion.
  - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE.  
Prepared surfaces should appear moist (glistening).
- RESTORATION PREPARATION**  
Make sure that the restoration is pretreated and handled according to the manufacturers instructions.
- DISPENSING**
  - Remove the cartridge cover.
  - Depress the lever to dispense the required amounts of paste onto the mixing pad.  
Note:  
1) Before dispensing the first time from a new cartridge or after a long interval in between use, bleed the pastes from the cartridge.
  - The adjusting guide can be moved forwards or backwards to change the quantity of material to be dispensed (Fig. 8).  
Note:  
1. Position of adjusting guide for one inlay.  
2. Position of adjusting guide for one crown.
  - While depressing the lever (Fig. 9-1), level and cut off the extruded material from the cartridge tips by moving the cartridge and dispenser to an upright position on the mixing pad (Fig. 9-2).
  - When the lever is released after dispensing, residual pastes at the cartridge tips are retracted into the cartridge. Clean the cartridge tip and replace the cartridge cover.  
Note:  
The amount of residual material in the cartridge can be confirmed by the position of the sliding block on the underside of the dispenser. As the cartridge cylinders are translucent, the residual amount can also be confirmed by the position of the inner plunger seen through the cylinders.
- MIXING**  
After dispensing, spread the two pastes thinly on the mixing pad using a plastic spatula. Mix thoroughly with lagging strokes, for 10 seconds. Be careful not to incorporate air bubbles.  
Note:  
1) The two pastes are similar in color, so mix with care. Larger amount of paste requires a longer mixing time (15 seconds).  
2) Do not mix GC FujiCEM pastes with other materials.

**5. CEMENTATION**

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 15 seconds from start of mixing at 23°C (74°F). Higher temperature will shorten working time. Refrigeration of the cartridge is recommended to extend the working time.
- Maintain moderate pressure.
- Start removing excess cement when the excess cement feels rubbery.
- Finishing can be started 4 minutes 30 seconds after seating the restoration.

**CLEANING / STERILIZATION OF PASTE PAK DISPENSER**  
For cleaning and sterilization of PASTE PAK DISPENSER, refer to the respective instructions for use.

**SHADE**  
Yellow

**STORAGE**  
Store in a cool and dark place (4-25°C) (39.2-77.0°F). If not in use for a prolonged period of time, store in refrigerator. (Shelf life : 2 years from date of manufacture)

**PACKAGES**

- GC FujiCEM Refill (Automix compatible)  
Paste Pak Cartridge (13.3g / 7.2mL) (2) with mixing pad (No.22)
- GC FujiCEM Automix  
Paste Pak Cartridge (13.3g / 7.2mL) (2), GC FujiCEM Mixing Tip (44)
- Paste Pak Dispenser (1)

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- In rare cases metal may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of Paste Pak Dispenser (nickel-plated aluminum) and refer to a physician.
- Store Paste Pak Dispenser away from chemicals such as iodine tincture containing halogen element. Otherwise the dispenser may discolor.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC FujiCEM in combination with eugenol containing materials as eugenol may hinder GC FujiCEM from setting or bonding properly.
- This product is not indicated for filling or core build-up.

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