Prior to use, carefully read the instructions for use.

EN

## GC Fuji II CAPSULE

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

#### RECOMMENDED INDICATIONS

- Restoration of primary teeth.
- Core build-up.
- Restorations of Class III, V and limited Class I cavities.

#### CONTRAINdicATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

#### DIRECTIONS FOR USE

Verhältnis Powder / Liquid Ratio (g/g)	0,30 / 0,11
Mixing Time (sec.)	10"
Working Time (min., sec.)	2'00"
Net Setting Time (min., sec.)	2'20"

Test conditions : Temperature (23 +/-1°C), Relative humidity (50 +/-5%)

ISO 9917 (Glass polyalkenoate cement)

#### 1. CAPSULE ACTIVATION AND MIXING

- Before activating, shake the capsule or tap on a hard surface to loosen the powder (Fig. 1).
  - To activate the capsule push plunger until it is flush with the main body (Fig. 2).
  - Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 3). The capsule is now activated.







Note :

The capsule should be activated just before mixing and used immediately.
  - Immediately set the capsule into the mixer and mix for 10 seconds (+/- 4,000 RPM) (Fig. 4).
- RESTORATIVE TECHNIQUE
    - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary.
    - Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer. For pulp capping use calcium hydroxide.
    - Mix the cement. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
    - Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
    - Make two clicks to prime the capsule then syringe (Fig. 5).
    - Remove surface moisture, but DO NOT DESICCATE.
    - Extrude cement directly into preparation. Avoid air bubbles.

Note :

      - To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
      - To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.
    - Form the preliminary contour and cover with a matrix if required.
    - Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC FUJI COAT LC (light cure).
  - FINISHING
    - One visit treatment
      - After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.
      - Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.
      - Final finishing and polishing under water spray can now be carried out.
      - Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.)
    - Standard method

24 hours after placement finish and polish the restoration under water spray using standard techniques.

#### STORAGE

Store at temperature of 4-25°C (39.2-77.0°F). (Shelf life : 2 years from date of manufacture).

#### SHADE

No. 21	pale yellow	No. 23	dark grey
No. 22	yellow brown	No. 23M	deep grey
No. 22.5	grey brown	No. 26	deep brown

#### PACKAGES

- Assortment package
    - 50 capsules (5 pcs. No. 21, 26, 10 pcs. No. 22, 22.5, 23, 23M).
  - Single shade package : 50 capsules
    - \* Average contents per capsule : 0.30g powder and 0.11g (0.09mL) liquid.
    - \* Minimum net volume of mixed cement per capsule : 0.12mL.
- Option : GC Capsule Applier (1 piece).

#### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

Last revised : 12/2008

**CE** 0086

**'GC'**

MANUFACTURED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-Ku, Tokyo 174-8585, Japan

DISTRIBUTED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-Ku, Tokyo 174-8585, Japan

**EU** : GC EUROPE N.V.  
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,  
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900.

GC ASIA DENTAL PTE. LTD.  
19 Loyang Way, #06-27 Singapore 508724  
TEL: +65 6546 7389

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