GC Fuji IX GP™ CAPSULE

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0.40 / 0.11
Mixing Time (sec.)*	10"
Time of initial extrusion after mixing	10 seconds maximum
Working Time (min., sec.)*	2'00"
Net Setting Time (min., sec.)*	2'20"
Final Finishing Commencing Time (37°C, after starting mix) (min.,sec.)	6'00"

Test conditions : Temperature (23 +/-1°C),
Relative humidity (50 +/-5%)

1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).
- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- To activate the capsule, push the plunger until it is flush with the main body and hold it down for 2 seconds (Fig. 3).

Note :

- Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
- The capsule should be activated just before mixing and used immediately.

- Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000RPM) (Fig. 4).

3. RESTORATIVE TECHNIQUE

- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
 - Make two clicks to prime the capsule then syringe (Fig. 5). The working time is 2 minutes from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
 - Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 6).
 - Form the preliminary contour, and cover with a matrix if required.
 - When set, immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure).
- Note :
- To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

4. FINISHING

- Final finishing under water spray using standard techniques can begin approximately 6 minutes after start of mixing (Fig. 7).
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 8,9).
- Instruct the patient not to apply pressure for 1 hour.

STORAGE

Recommended for optimal performance, store the capsules in the original aluminum foil in a cool and dark place (4-25°C)(39.2-77.0°F).

SHADES

A2, A3, A3.5, B2, B3 and C4.

Shade numbers according to Vita® shade guide.

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package : 48 capsules.
- Assortment package : 48 capsules (10 each of A2, A3 and A3.5; 9 each of B2 and B3).
Average contents per capsule : 0.40g powder and 0.11g (0.09mL) liquid.
Minimum net volume of mixed cement per capsule : 0.14mL.
Option : GC CAPSULE APPLIER (1 piece).
- Single Shade Export package: 200 capsules.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at:
<http://www.gcamerica.com>

They can also be obtained from your supplier.

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