GC Fuji IX GP FAST

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT IN CAPSULES

Prior to use, carefully read the instructions for use.

GC Fuji IX GP FAST is a fast-setting version of GC Fuji IX GP CAPSULE.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE


	Fast	Original
Powder / Liquid Ratio (g/g)	0,40 / 0,11	0,40 / 0,11
Mixing Time (sec.)*	10"	10"
Time of initial extrusion after mixing	10 seconds maximum	
Working Time (min., sec.)*	1'15"	2'00"
Net Setting Time (min., sec.)*	2'00"	2'20"
Final Finishing Commencing Time (37°C, after starting mix) (min.,sec.)	3'00"	6'00"

Test conditions : Temperature (23 +/-1°C), Relative humidity (50 +/-5%)

*ISO 9917 : 1991 (E) (Glass polyalkenoate cement)




1. CAVITY PREPARATION

- Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).
- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).




2. CAPSULE ACTIVATION AND MIXING

- Before activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- To activate the capsule, push the plunger until it is flush with the main body (Fig. 3).
- Immediately place the capsule into a metal GC Capsule Applier and click the lever once (Fig. 4). The capsule is now activated.
Note :
The capsule should be activated just before mixing and used immediately.
- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000RPM) (Fig. 5).

3. RESTORATIVE TECHNIQUE



- Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.
- Make two clicks to prime the capsule then syringe (Fig. 6). The working time is 1 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 7).
- Form the preliminary contour, and cover with a matrix if required.
- When set, immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) (Fig. 8).

- Note :
- To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
 - To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

4. FINISHING

- Final finishing under water spray using standard techniques can begin approximately 3 minutes after start of mixing (Fig. 9).
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) (Fig. 10).
- Instruct the patient not to apply pressure for 1 hour.

STORAGE

Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life : 2 years from date of manufacture).

SHADES

A1, A2, A3, A3.5, B2, B3, C4
Shade numbers according to Vita® shade guide.

Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- Single shade package : 50 capsules.
- Assortment package : 50 capsules (10 each of A2, A3, A3.5, B2, B3).
Average contents per capsule : 0.40g powder and 0.11g (0.09mL) liquid.
Minimum net volume of mixed cement per capsule : 0.14mL.

Option : GC CAPSULE APPLIER (1 piece).

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

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