Prior to use, carefully read the instructions for use. **EN**

GRADIA® CORE

DUAL-CURED RADIOPAQUE COMPOSITE FOR CORE BUILD-UP & POST CEMENTATION

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Core build-up.
- Post cementation.


CONTRAINDICATIONS

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Prior to use, take out the GRADIA CORE SELF-ETCHING BOND and GRADIA CORE cartridge from the refrigerator and allow to stand for at least 5 minutes at room temperature.

1. LOADING CARTRIDGE INTO DISPENSER

Insert a cartridge into GRADIA CORE DISPENSER GUN and push the plunger of the dispenser until it engages with the syringe. Load the cartridge upright, and remove the cap by rotating 1/4-tum counterclockwise.

Before every use, blend a small amount of pastes from the cartridge to ensure the even flow of both pastes. Attach a GC AUTOMIX TIP For Endo by aligning the V shaped notch on the rim of the mixing tip with the V shaped notch between the cartridge barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 tum clockwise. If desired, the extension Tip can be attached to the mixing tip. Push the extension tip firmly to the mixing tip until it clicks.

Note:
a) GRADIA CORE Cartridge is to be used always with GRADIA CORE DISPENSER GUN.
b) Use GC AUTOMIX TIP For Endo. Do not use any other mixing tips to avoid the inadequate mixing.
c) Hold the base part of the extension tip when attaching to the mixing tip.

2. CAVITY PREPARATION
Isolate the tooth with a rubber dam. Prepare root canal (s) using a standard technique. The cleaned root canal (s) should be rinsed and dried thoroughly. Place matrix around the prepared tooth if necessary.

In case of a core build up in vital teeth, make sure that the pulp is appropriately protected.

3. APPLICATION GRADIA CORE SELF-ETCHING BOND
a) Shake the bottles of GRADIA CORE SELF-ETCHING BOND LIQUID A and B thoroughly. Dispense one drop of each into the dispensing dish and mix thoroughly for 5 seconds using the micro-tip applicator. Remove any residual material around the nozzles and replace bottle caps immediately after use.

Note:
1) Always use GRADIA CORE SELF-ETCHING BOND together with GRADIA CORE.
2) Use the dispensing dish designated for GRADIA CORE SELF-ETCHING BOND. After use, clean the dish immediately using a cotton roll moistened with alcohol.

3) GRADIA CORE SELF-ETCHING BOND is volatile. Dispense immediately before use. After mixing, immediately proceed to the next step. Be sure to use self-etching Bond within one minute after mixing.

4) GRADIA CORE SELF-ETCHING BOND is dual-cured material. It starts to cure under ambient light and by chemical reaction.

5) Make sure that the micro-tip applicator is firmly pushed into the application holder.

b) Apply the mixture on the prepared enamel and/or dentin surfaces.

c) Leave for 30 seconds. Blot any excess material away with a paper point. Dry thoroughly for 10 seconds with oil-free air until the adhesive film no longer moves.

4) unit (Table 1).

Light curing device (Light intensity >450mW/cm²)	Time
Halogen	10 seconds
GC G-Light LED (Wavelength range 460-480nm)	10 seconds
Plasma Arc (470 Tip)	3 seconds

Note:
1) Light curing of GRADIA CORE SELF-ETCHING BOND prior to the application of GRADIA CORE is essential to achieve the optimal bonding strength.

2) If there is more than one root canal involved, light cure in the direction of each root canal for the recommended time.

4. POST CEMENTATION & CORE BUILD-UP

a) After having light cured the bonding agent, dispense GRADIA CORE into the prepared root canal. Keep the tip immersed in the cement while dispensing to avoid trapping air.

Do not use a lentulo spiral to place the material in the root canal.

b) Insert post which is prepared and pre-treated according to the manufacturer's instructions. Light cure for a few seconds (tack cure) to temporarily fix the position of the post.

Note:
Post length should be adjusted prior to cementation to prevent dislodging.

c) Continue to dispense the paste around the post to form the core. Working time is 3 minutes from start of mixing at 23°C (73.4°F). Higher temperature will shorten the working time.

Note:
1) For easier manipulation of a large amount of material, it is recommended to build up incrementally and light each time for a few seconds. Refer to the table 2 for the light curing time for tack cure.

2) The used mixing tip can be left on the syringe to serve as a storage cap until next use. Prior to the next use, place a new mixing tip.

d) Light cure for final setting from labial (buccal) and lingual side using a visible light curing unit or allow the material to self-cure for 5 minutes. Refer to Table 2 for the light curing times.

Table 2	Light curing device (Light intensity >450mW/cm²)	Light curing time (Incremental build-up) (Depth of cure: 2.5 mm)	Final Light-curing Time (Depth of cure: 2.5 mm)
Halogen	5 seconds	20 seconds	20 seconds
GC G-Light LED (Wavelength range 460-480nm)	5 seconds	10 seconds	10 seconds
Plasma Arc (470 Tip)	3 seconds	3 seconds	3 seconds

GRADIA CORE can be layered in thicknesses of 2.5 mm.

Note:
1) Longer curing time is required if the root canal is deep or there is more than one root canal. (Recommended curing time x the number of root canals is recommended.)

2) GRADIA CORE is dual cured material. Allow the material to self cure for 5 minutes if the light cannot reach the material.

5. CONTOURING

Once fully cured, prepare the core in a standard technique.

6. ADDITIONAL BUILD-UP

Roughen the surface of the cured composite. Apply freshly mixed GRADIA CORE SELF-ETCHING BOND on the build-up with GRADIA CORE following the procedure described above.

STORAGE

Store in refrigerator (2-8°C / 35.6-46.4°F). (Shelf life: 2 years from date of manufacture)

SHADE

Universal shade

PACKAGES

1. GRADIA CORE KIT
Cartridge 10 mL (20g) (1), GC AUTOMIX TIP For Endo with extension tip (20), GRADIA CORE SELF-ETCHING BOND LIQUID A 3.0 mL (1), GRADIA CORE SELF-ETCHING BOND LIQUID B 1.5 mL (1), Micro-Tip Applicator (25), Applicator Holder (1), Dispensing Dish (1)

2. GRADIA CORE REFILL
Cartridge 10 mL (20g) (1), GC AUTOMIX TIP For Endo with extension tip (20)

3. GRADIA CORE SELF-ETCHING BOND LIQUID A 3.0 mL (1)

4. GRADIA CORE SELF-ETCHING BOND LIQUID B 1.5 mL (1)

5. GC AUTOMIX TIP For Endo including extension tip (10)

6. GRADIA CORE DISPENSER GUN (1)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cotton butter can be used to isolate the operation field from the skin or oral tissue. Use gloves to avoid contact.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

3. If the tissue contacted by GRADIA CORE SELF-ETCHING BOND turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks.

4. Avoid inhalation of vapor or digestion of GRADIA CORE SELF-ETCHING BOND or GRADIA CORE.

5. Do not use GRADIA CORE and GRADIA CORE SELF-ETCHING BOND in combination with eugenol containing materials as eugenol may hinder GRADIA CORE and GRADIA CORE SELF-ETCHING BOND from setting or bonding properly.

6. Do not mix GRADIA CORE or GRADIA CORE SELF-ETCHING BOND with other products. Do not use GRADIA CORE SELF-ETCHING BOND with other composites or cementitious materials.

7. GRADIA CORE SELF-ETCHING BOND liquids are flammable. Do not use near naked flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.

8. GRADIA CORE SELF-ETCHING BOND liquids are volatile. Use in a well ventilated place. Close bottles immediately after use.

9. When wiping GC AUTOMIX TIP or extension tip with a cotton soaked with alcohol, prevent the alcohol from entering the inside of the tip.

10. Micro-Tip Applicator, GC AUTOMIX TIP and extension tip are for single use only and cannot be autoclaved or chemically.

11. GRADIA CORE DISPENSER GUN can be autoclaved. The operating temperature should not exceed 135°C (275°F). For chemical cleaning, do not use sodium hypochlorite, sodium percarbonate, phenol, or cresol.

12. Wear protective eyeglasses during light curing.

13. Do not mix with other similar products.

14. Avoid getting material on clothing or equipment.

15. When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of curing dust.

16. Dispose of all waste according to local regulations.

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