Prior to use, carefully read the instructions for use.

EN

GC RELINE BONDING AGENT

GC RELINE BONDING AGENT is a specially designed methacrylate base bonding agent that enhances the bond between GC Reline Material and denture base materials.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS:

For use as a bonding agent when processing GC RELINE Hard Denture Chairside Reline Material.

CONTRAINDICATIONS:

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE:

Note: Consult the instructions for use for the GC RELINE Hard Denture Chairside Reline Material for complete processing information.

1. Roughen the area of the denture to be relined. Clean and dry the denture thoroughly.
2. Coat the artificial teeth, labial and buccal surfaces of the denture with COE LUBRICANT or Vaseline using a cotton swab being careful to stay at least 3 mm (1/8 inch) from the peripheral border. *Do not contaminate the tissue side of the denture with lubricant.*
3. Using the brush provided, coat the area to be relined with a uniform coating of the GC RELINE BONDING AGENT. Air-dry the bonding agent for approximately 10 seconds to allow the solvent to evaporate.
Note: GC RELINE BONDING AGENT is highly volatile; therefore be sure to replace and tightly close the cap after each use.
4. Continue processing the GC RELINE Hard Denture Chairside Reline Material following the recommended procedures.

CAUTIONS:

1. The solvent in GC RELINE BONDING AGENT is highly volatile and flammable. Avoid sources of ignition and keep tight capped when not in use.
2. In case of contact with skin, wash the affected area thoroughly.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. Wash the brush with GC RELINE liquid or ethanol after using the GC RELINE BONDING AGENT.

STORAGE:

Store at normal temperature 15-25 °C (60 - 77 °F)
(Shelf life: 3 years from date of manufacture)

Last revised: 09/2011

MANUFACTURED by
GC AMERICA INC.
3737 West 127th Street
Alsip, Illinois 60803 U.S.A.

DISTRIBUTED by
EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240,
Interleuvenlaan 33,
B-3001 Leuven, Belgium
TEL: +32 16 74 10 00

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