GC RELINE™ II Soft GC RELINE™ II ExtraSoft GC RELINE™ II ExtraExtraSoft

VINYL POLYSILOXANE SOFT LINING MATERIAL
GC RELINE™ II PRIMER FOR RESIN
RESIN PRIMER/ADHESIVE FOR VINYL POLYSILOXANE SOFT LINING MATERIAL
GC RELINE™ II FINISHING MATERIAL
COATING MATERIAL FOR VINYL POLYSILOXANE SOFT LINING MATERIAL
GC RELINE™ II POINT FOR TRIMMING
POINT FOR TRIMMING SOFT LINING MATERIAL
GC RELINE™ II WHEEL FOR FINISHING
WHEEL FOR FINISHING SOFT LINING MATERIAL
GC RELINE™ II REMOVER FOR RESIN
REMOVER FOR SOFT LINING MATERIAL
For use only by a dental professional in the recommended indications.

- INDICATIONS FOR USE**
- Relining of dentures in case of thin sensitive mucosa and resorbed ridges.
 - Relining of dentures in case of sharp and bony ridges.
 - Relining of dentures in case of bony undercuts (with difficulties involving insertion and/or removal of the denture).
 - Relining of implant-supported dentures


CONTRAINDICATIONS
Avoid using this product with patients having known silicone allergies.

DIRECTIONS FOR USE
A) DIRECT RELINING METHOD

- Check the resin denture for misfit and sensitive areas of the mucosa (e.g. with Fit Checker Advanced or Fit Checker Advanced Blue).
- Disinfect with a non-alcohol based product and dry the denture thoroughly. Relieve the area to be relined and roughen it, ideally by microabrasion.
- Clean in an ultrasonic cleaning machine containing an ultrasonic or enzymatic solution and dry thoroughly. Apply GC RELINE II PRIMER FOR RESIN to the bonding surface with a clean and dry brush and gently dry with clean, dry and oil-free air. Do not contaminate the surface with latex gloves, powder from gloves, etc. Note: Close the bottle with the cap immediately after use. Do not contaminate the primer bottle with moisture, in order to avoid that the primer liquid becomes blurry.
- Dispense GC RELINE II Soft, ExtraSoft or ExtraExtraSoft directly onto the denture. The working time is 2 minutes at room temperature 23°C (73.4°F). Note: 1) Do not apply a small amount of onset material (1-2 mm of mixed material) as this is slow to set and you will need to retain the denture in the mouth for over 5 minutes. 2) Do not use a solvent like alcohol on the area between GC RELINE II and the denture as this will cause separation.
- Seat the denture in the mouth. Ask the patient to gently close into centric occlusion and then perform muscle trimming.
- Retain the denture in place for at least 5 minutes and remove it from the mouth.
- Using a scalpel or scissors, remove excess material on the edge of GC RELINE II Soft, ExtraSoft or ExtraExtraSoft. Using light pressure, progressively trim with GC RELINE II POINT FOR TRIMMING using a speed of less than 20,000 RPM, and finish with GC RELINE II WHEEL FOR FINISHING using a speed of less than 6,000 RPM. Before use, the wheel must be set on a mandrel with a fixation cap or screw to hold it securely in place.
- If it is considered necessary to achieve an even smoother finish, follow the procedure below: (NB: GC RELINE II FINISHING MATERIAL A and B are sold separately) a) Apply GC RELINE II PRIMER FOR RESIN to the peripheral borders and margins of the denture in order to bond GC RELINE II FINISHING MATERIAL to it. Do not apply the primer to the relining material. Gently dry the primer (as in step 3 above). b) Dispense equal amounts of GC RELINE II FINISHING MATERIAL A and B onto the mixing pad. Mix with a plastic spatula and apply the mixture to the edge of the previously placed GC RELINE II Soft, ExtraSoft or ExtraExtraSoft material and to the denture using a sponge. The working time is 2 minutes and the setting time is 4 minutes at room temperature 23°C (73.4°F). If it is necessary to check the fit of the relined denture, use an uncured fit checking material.

B) INDIRECT RELINING METHOD

- Proceed according to Steps No. 1-4 for Direct Relining Method.
- Mount the denture on the gypsum model and press it on an articulator.
- Remove the denture from the model after about 30 minutes.
- If required, proceed according to Steps No. 7- 8 for Direct Relining Method.

MAINTENANCE OF THE RELINED DENTURE

- Advise the patient to clean the relined denture under running water, with a soft brush and soap, at least once daily.
- Within 3 months of applying GC RELINE II Soft, ExtraSoft or ExtraExtraSoft check the denture for any accumulation of calculus and check the condition of the oral tissues. If necessary, relieve the pressure spots and apply additional GC RELINE II Soft, ExtraSoft or ExtraExtraSoft material to the relieved area.
- If an additional layer is required to be applied to an existing relining, relieve the original relined area with GC RELINE II POINT FOR TRIMMING and apply material within 3 months of the original application.
- To replace an existing relined material, follow the below described procedure: a) Expose the acrylic surface by cutting the silicon material with a scalpel. b) Put a few drops of GC RELINE II REMOVER FOR RESIN on the surface of the acrylic denture. c) Insert the instrument or spatula between the silicone and acrylic resin interface and remove the silicone material from the acrylic resin with the instrument or spatula. d) After removing the silicone material, relieve the fitting surface of the denture to provide a completely clean surface. Afterwards, follow steps No. 2-8 in the instructions above.

STORAGE
Recommended for optimal performance, store at room temperature (15-25°C) (59-77°F) and out of direct sunlight. (Shelf life: 2 years from the date of manufacture)

PACKAGES
1. INTRODUCTORY KIT
GC RELINE II CARTRIDGE (One of the following cartridges: Soft 61g (48mL), ExtraSoft 57g (48mL), or ExtraExtraSoft 55g (48mL) (1)
GC RELINE II PRIMER FOR RESIN 8g (10mL) (1)
GC RELINE II POINT FOR TRIMMING (1)
GC RELINE II WHEEL FOR FINISHING (1)
GC RELINE II REMOVER FOR RESIN 5g (5.5mL) (1)
MIXING TIP II L (6 pieces)
Disposable brush (3)
Plastic spatula (1)
2. REFILL PACKAGES
a) GC RELINE II CARTRIDGE (One of the following cartridges: Soft 61g (48mL), ExtraSoft 57g (48mL), or ExtraExtraSoft 55g (48mL) (1)
b) GC RELINE II PRIMER FOR RESIN 8g (10mL) (1)
c) GC RELINE II FINISHING MATERIAL A and B 10g (9mL) (1 each)
d) GC RELINE II POINT FOR TRIMMING (3)
e) GC RELINE II WHEEL FOR FINISHING (12)
f) GC RELINE II REMOVER (FOR RESIN) 5g (5.5mL) (1)
g) Disposable brush (6)

CAUTION

- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- The primer and the remover are volatile and flammable. Do not use them near fire. Close the bottle immediately after use and keep away from sources of ignition.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Avoid getting the relining material on clothing. It is difficult to remove once it has set and penetrated the fabric.
- When applying the relining material, care should be taken to avoid mixing or contact with the following materials. They may delay setting time or result in poor adhesion. * Latex gloves, powder from gloves * Condensation type silicone impression material * Polysulfide impression material * Eugenol type material * Soft relining material from other companies * Uncured acrylic resin
- Before using, check that the shaft of GC RELINE II POINT FOR TRIMMING or GC RELINE II WHEEL FOR FINISHING rotates correctly without distortion.
- Rotate the point or wheel from the relining material toward the denture base.
- Do not use excessive pressure when using the point or wheel.
- When trimming or finishing the relining material, wear a dust mask to avoid inhaling silicone dust, and wear protective eye glasses.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gcamerica.com> They can also be obtained from your supplier.

CLEANING AND DISINFECTING MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. **DO NOT IMMERS**e. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised: 02/2017
Rx Only
GC
MANUFACTURED BY GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC.
3737 West 127th Street, Atsip, IL 60803 U.S.A.
TEL: +1-708-907-0900
www.gcamerica.com
DISTRIBUTED BY GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
GC AMERICA INC.
3737 West 127th Street, Atsip, IL 60803 U.S.A.
TEL: +1-708-907-0900
www.gcamerica.com
PRINTED IN JAPAN

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INDICATIONS FOR USE

- Relining of dentures in case of thin sensitive mucosa and resorbed ridges.
- Relining of dentures in case of sharp and bony ridges.
- Relining of dentures in case of bony undercuts (with difficulties involving insertion and/or removal of the denture).
- Relining of implant-supported dentures

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
DIRECTION FOR USE
A) DIRECT RELINING METHOD

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B) INDIRECT RELINING METHOD

- Proceed according to Steps No.1-4 for Direct Relining Method.
- Mount the denture on the gypsum model and press it on an articulator.
- Remove the denture from the model after about 30 minutes.
- If required, proceed according to Steps No. 7- 8 for Direct Relining Method.

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PACKAGES
1. INTRODUCTORY KIT
GC RELINE II CARTRIDGE (One of the following cartridges: Soft 61g (48mL), ExtraSoft 57g (48mL), or ExtraExtraSoft 55g (48mL) (1)
GC RELINE II PRIMER FOR RESIN 8g (10mL) (1)
GC RELINE II POINT FOR TRIMMING (1)
GC RELINE II WHEEL FOR FINISHING (1)
GC RELINE II REMOVER FOR RESIN 5g (5.5mL) (1)
MIXING TIP II L (6 pieces)
Disposable brush (3)
Plastic spatula (1)
2. REFILL PACKAGES
a) GC RELINE II CARTRIDGE (One of the following cartridges: Soft 61g (48mL), ExtraSoft 57g (48mL), or ExtraExtraSoft 55g (48mL) (1)
b) GC RELINE II PRIMER FOR RESIN 8g (10mL) (1)
c) GC RELINE II FINISHING MATERIAL A and B 10g (9mL) (1 each)
d) GC RELINE II POINT FOR TRIMMING (3)
e) GC RELINE II WHEEL FOR FINISHING (12)
f) GC RELINE II REMOVER (FOR RESIN) 5g (5.5mL) (1)
g) Disposable brush (