

### DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

#### Bactericidal, Tuberculocidal, Fungicidal and Virucidal

This product is a general purpose cleaner, deodorizer and disinfectant to be used where control of cross-contamination hazards is of prime importance. Use this product on all hard, non-porous environmental surfaces/articles including: Chrome, Formica® Brand, glass, painted surfaces, plastics (such as polycarbonate, polyvinylchloride, polypropylene and polystyrene), Plexiglas®, glazed porcelain and stainless steel.

This product is for use in dental offices to clean, deodorize and disinfect such items as dental chairs/ countertops/instruments/implements/OR surfaces, bathroom fixtures, carts, chairs, floors, sinks, tables, telephones, toilets (exterior surfaces, seats and rims) and work stations.

**Cleaning, PreCleaning and Deodorizing:** Prior to disinfection, preclean all surfaces with GC SPRAY-CIDE™, making sure to remove all gross filth and heavy soil. Apply solution to surface and allow surface to remain wet for 30 seconds prior to wiping with a paper towel or clean towel/sponge.

#### To Disinfect Hard, Nonporous Surfaces/Articles

**To Disinfect** *Acinetobacter baumannii*, *Aspergillus niger*, Avian Influenza A Virus, *Candida albicans*, Canine Parvovirus, *Escherichia coli* O157:H7, Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Herpes Simplex Virus Type 2, HIV-1 (AIDS Virus), Influenza Strain A2/HK Virus, *Mycobacterium bovis* BCG, Pandemic 2009 H1N1 Influenza A Virus\*, Polio I Virus, *Pseudomonas aeruginosa*, Rhinovirus, Rotavirus Strain WA, *Salmonella enterica*, *Staphylococcus aureus*, Methicillin-resistant *Staphylococcus aureus* (MRSA), *Trichophyton mentagrophytes* and Vancomycin-resistant *Enterococcus faecalis* (VRE): Using GC SPRAY-CIDE as is (do not dilute), apply solution to thoroughly wet the precleaned surfaces. Allow to remain wet for **6 minutes** prior to wiping or allow surfaces to air dry. For disinfection of Avian Influenza A, use this product only for surfaces that are not conducive to treatment by immersion or excess liquid.

**Note:** Disinfected food contact surfaces must be rinsed with potable water after cleaning and/or disinfection and before reuse. After disinfection, critical and semi-critical devices must undergo appropriate terminal sterilization/high-level disinfection.

# GC SPRAY-CIDE™

## GERMICIDAL SOLUTION

### CLEANER, DISINFECTANT & DEODORIZER

**Net Contents: 1 Gallon (3.8 L)**

**Reorder No. 555128**

#### ACTIVE INGREDIENTS:

n-Alkyl (60%C14, 30%C16, 5%C12, 5%C18)  
dimethyl benzyl ammonium chloride . . . . . 0.154%  
n-Alkyl (68%C12, 32%C14) dimethyl  
ethylbenzyl ammonium chloride. . . . . 0.154%  
Isopropanol . . . . . 21.000%

**OTHER INGREDIENTS** . . . . . 78.692%

**TOTAL** . . . . . 100.000%

Manufactured for:



GC AMERICA INC.  
3737 W. 127th STREET  
ALSIP, IL 60803  
Made in U.S.A.

## KEEP OUT OF REACH OF CHILDREN CAUTION

See side panel for First Aid and additional Precautionary Statements.

This product has demonstrated effectiveness against Influenza Strain A2/HK Virus\* and Avian Influenza A Virus\* and is expected to inactivate all influenza A Viruses including Pandemic 2009 H1N1 (formerly called "swine flu").

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

#### Special Instructions for Cleaning and Decontamination Against HBV, HCV and HIV-1 (AIDS Virus) for Surfaces or Objects Soiled with Blood or Body Fluids:

Kills HBV, HCV and HIV-1 (AIDS Virus) on precleaned environmental surfaces or objects previously soiled with blood or body fluids at room temperature (20-25°C) in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of HBV, HCV and the human immunodeficiency Virus Type 1 (HIV-1) associated with AIDS. **Personal Protection:** When using the germicidal solution, wear disposable latex gloves, protective gown, face mask and eye coverings, as appropriate, when handling HBV, HCV and HIV-1 infected blood or body fluids. **Cleaning Procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the germicidal solution. **Contact Time:** Thoroughly wet surface. Allow surface to remain wet for 6 minutes prior to wiping or allow surface to air dry. **Disposal of Infectious Materials:** Dispose of used solution in accordance with local regulations for infectious waste disposal.

#### STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

**Pesticide Storage:** Store in original container in a cool, well-ventilated (preferably locked) area inaccessible to children and pets. **Pesticide Disposal:** Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility. **Container Disposal:** Non-refillable container; do not reuse or refill this container. To avoid wastes, use all material in this container by application according to label directions. When empty, fill the container about 10 percent full of water. Agitate vigorously and dispose of contents on site or at an approved waste disposal facility. Offer the container for recycling, if available; otherwise dispose of container in trash.

#### FIRST AID

<b>IF INHALED</b>	Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.
<b>IF ON SKIN OR CLOTHING</b>	Take off contaminated clothing. Immediately rinse skin with plenty of water for 15-20 minutes.
<b>IF IN EYES</b>	Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-424-9300 for emergency medical treatment information.

#### PRECAUTIONARY STATEMENTS

##### Hazards to Humans and Domestic Animals

**CAUTION:** Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum or using tobacco. Remove contaminated clothing and wash clothing before reuse.

**Physical or Chemical Hazards:** Do not use or store near heat or open flame.

EPA Reg. No. 1130-15-10214 EPA Est. 1130-IL-1



\* + D6585551280D \*