GRADIA® DIRECT X LIGHT-CURED COMPOSITE RESTORATIVE

For use only by a dental professional in the recommended indications. RECOMMENDED INDICATIONS

- 1. Direct restorative for Class I, II, III, IV, V cavities. 2. Direct restorative for wedge-shaped defects and root surface cavities. 3. Direct restorative for veneers and diastema closure.

- CONTRAINDICATIONS 1. Pulp capping. 2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Shade Selection Clean the tooth with pumice and water. Shade selection should be made prior to isolation. Select the appropriate GRADIA DIRECT X shades by referring to the GRADIA DIRECT shade guide.

Cavity Preparation Prepare cavity using standard techniques. Dry thoroughly blowing with oil free air.

Note: For pulp capping, use calcium hydroxide. Bonding Treatment For bonding GRADIA DIRECT X to enamel and /or dentin, use a light-cured bonding system such as GC Fuji BOND LC, GC UniFil® Bond or GC G-BOND™ (Fig. 1).

Placement of GRADIA DIRECT X 1) Dispensing from a Unitip Insert the GRADIA DIRECT X Unitip into the Unitip APPLIER or equivalent. Remove the cap and extrude material directly into the prepared cavity. Use steady pressure (Fig. 2). Maintain pressure on the applicator handle while removing the Unitip APPLIER and Unitip from the mouth. This prevents the Unitip from coming loose from the applicator.

2) Dispensing from a syringe Remove syringe cap and dispense material onto a mixing pad. Place the material into the cavity using a suitable placement instrument. After dispensing, screw syringe plunger anticlockwise by a half to full turn to release residual pressure inside the syringe. Replace cap immediately after use.

3. After dispensing, avoid too long an exposure to ambient light. Ambient light can shorten the manipulation time.

- Clinical Hints a. In the case of small cavities Restore using a one shade technique. In most cases the use of one standard shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Outside special shades can be selected. b. In the case of large and/or deep cavities In most cases a multi shade layering technique will give the best aesthetic results. To block out shine through from the oral cavity or to mask discolored dentine, place X-AO2 and continue to build up with a standard shade. For optimal aesthetics use a Outside special shade as the final composite layer. In the case of deep posterior cavities, a flowable composite such as GC GRADIA DIRECT Flo /LoFlo or a glass ionomer cement such as GC Fuji LINING LC (Paste Pak) or GC Fuji IX GP can be used on the cavity floor instead of an Opaque shade. See also Examples of Clinical Applications and/or Shade Combination Chart.

Examples of Clinical Applications (Clinical Hint a, b) (a) Standard (b) Outside special (c) X-AO2 or Flowable composite Multi shade layering technique

5. Contouring before Light Curing Contour using standard techniques. 6. Light Curing Light cure GRADIA DIRECT X using a light curing unit (Fig. 3). Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Table with 4 columns: Irradiation time, Plasma arc, GC-Light, Halogen /LED. Rows for X-A1, X-B1, X-WT; X-A2, X-A3, X-B2, X-C2, X-D2; X-A3.5, X-AO2, X-XBW, X-BW.

Table with 4 columns: Shade, Plasma arc, GC-Light, Halogen /LED. Rows for X-A1, X-B1, X-WT; X-A2, X-A3, X-B2, X-C2, X-D2; X-A3.5, X-AO2, X-XBW, X-BW.

Note: 1. Material should be placed and light cured in layers. For maximum layer thickness, please consult the tables. 2. Lower light intensity may cause insufficient curing or discoloration of the material.

Cleaning and Sterilizing of the Unitip APPLIER

1. Clean the applicator with gauze soaked in alcohol. 2. The applicator may be steam autoclaved at 121-135°C (250-275°F) for 10 to 30 minutes or disinfected using alcohol.

Cleaning of the Shade Guide The shade guide can be cleaned with water and mild soap. It should not be sterilized. The use of various chemical disinfectants / sterilants may cause damage to the shade guide.

SHADES Standard shade : X-XBW*, X-BW*, X-A1, X-A2, X-A3, X-A3.5, X-B1*, X-B2*, X-C2*, X-D2* Inside special shade : X-AO2 Outside special shade : X-WT *Not available in Europe.

PACKAGES I. Unitips 1. Pack of 20 tips (each in 5 shades) (0.16mL per tip) X-A1, X-A2, X-A3, X-A3.5, X-WT 2. Pack of 10 tips (each in 7 shades) (0.16mL per tip) X-XBW, X-BW, X-B1, X-B2, X-C2, X-D2, X-AO2

Option a. Unitip APPLIER b. Shade guide c. Mixing pad (No.14B) II. Syringes 1. 1 syringe (in 12 shades) (2.7mL per syringe) Note: Weight per syringe : 5.0g 2. Option a. Shade guide b. Mixing pad (No.14B)

CAUTION 1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.

2. In case of contact with eyes, flush immediately with water and seek medical attention. 3. Take care to avoid ingestion of the material. 4. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity. 5. For infection control reasons, Unitips are for single use only. 6. Wear protective eyeglasses during light curing. 7. When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust. 8. Do not mix with other similar products. 9. Avoid getting material on clothing. 10. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing. 11. Do not use GRADIA DIRECT X in combination with eugenol containing materials as eugenol may hinder GRADIA DIRECT X from setting.

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