Prior to use, carefully read the instructions for use. **EN**

GRADIA® DIRECT X

LIGHT-CURED COMPOSITE RESTORATIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Direct restorative for Class I, II, III, IV, V cavities.
- Direct restorative for wedge-shaped defects and root surface cavities.
- Direct restorative for veneers and diastema closure.

CONTRAINDICATIONS

- Pulp capping.
- Use of this product in patients with known allergies to methacrylate monomer or methacrylic polymer.

DIRECTIONS FOR USE




- Shade Selection
Clean the tooth with pumice and water. Shade selection should be made prior to isolation.
- Cavity Preparation
Prepare cavity using standard techniques.
Note:
For pulp capping, use calcium hydroxide.
- Bonding Treatment
For bonding GRADIA DIRECT X to enamel and / or dentin, use a light-cured bonding system such as G-Premio BOND, G-aenial Bond, G-BOND or GC FUJI BOND LC (Fig. 1). Follow manufacturer's instructions.
- Placement of GRADIA DIRECT X
1) Dispensing from a Unitip Insert the GRADIA DIRECT X Unitip into the Unitip APPLIER or equivalent. Remove the cap and extrude material directly into the prepared cavity. Use steady pressure (Fig. 2).
2) Push the Unitip on the applicator handle while removing the Unitip APPLIER and Unitip from the tooth. This prevents the Unitip from coming loose from the applicator.
3) Dispensing from a syringe
Mix the material and dispense material onto a mixing pad.
Place the material into the cavity using a suitable placement instrument. After releasing, screw syringe plunger anticlockwise by a half to full turn to release residual pressure inside the syringe. Replace cap immediately after use.
Note:
1. Material can be applied in a single layer to achieve aesthetic restorations using Standard shades. For details, refer to the Clinical Hints.
2. Material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at normal room temperature.
3. After dispensing, avoid too long exposure to ambient light. Ambient light can shorten the manipulation time.

Clinical Hints

- In the case of small cavities
Restore using a one shade technique. In most cases the use of one Standard shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Outside special shades can be selected.
In most cases a multi shade layering technique will give the best aesthetic results. To block out shine through from the oral cavity or to mask discolorated dentine, place X-AO2 and continue to build up with a Standard shade. For optimal aesthetics use a Outside special shade as the final composite layer. In the case of deep posterior cavities, a flowable composite such as GC GRADIA DIRECT Fio / LoFlo or a glass ionomer cement such as GC FUJI LINING LC (Paste Pak) or GC FUJI XG GP can be used on the cavity floor instead of an opaque shade. See also Examples of Clinical Applications.

Examples of Clinical Applications (Clinical Hint a, b)

Irradiation time	Plasma arc (2000mW/cm²)	3 sec.		6 sec.	
		High power LED (more than 1200mW/cm²)	10 sec.	20 sec.	40 sec.
Shade	Halogen / LED (700mW/cm²)	20 sec.	20 sec.	40 sec.	40 sec.

	3.0 mm	3.5 mm
X-A1, X-B1, X-WT	3.0 mm	3.5 mm
X-A2, X-A3, X-B2, X-C2, X-D2	2.5 mm	3.0 mm
X-A3.5, X-AO2, X-XBW, X-BW	2.0 mm	3.0 mm

- Finishing and Polishing
Finish and polish using diamond burs, polishing points and discs. To obtain a high gloss, polishing pastes can be used.

- Cleaning and Sterilizing of the Unitip APPLIER
1. Clean the applicator with gauze soaked in alcohol.
2. The applicator may be steam autoclaved at 121-135°C (250-275°F) for 10 to 30 minutes or disinfected using alcohol.
3. Do not attempt to sterilize in a unit at a temperature that exceeds 135°C (275°F) during the drying step of the autoclave procedure.
Otherwise, the applicator may be damaged.

- Cleaning of the Shade Guide
The shade guide can be cleaned with water and mild soap. It should not be sterilized. The use of various chemical disinfectants / sterilants may cause damage to the shade guide.

- SHADES
Standard shade : X-XBW*, X-BW*, X-A1, X-A2, X-A3, X-A3.5, X-B1*, X-B2*, X-C2*, X-D2*
Inside special shade : X-AO2
Outside special shade : X-WT
*Not available in Europe.

- A, B, C, D, AO shades are based on Vita® Shade.
*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

- STORAGE
Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77.0°F) away from high temperatures or direct sunlight. (Shelf life: 3 years from date of manufacture)

- PACKAGES
1. Unitips
a. Pack of 20 tips (each in 5 shades) (0.16mL per tip)
X-A1, X-A2, X-A3, X-A3.5, X-WT
2. Pack of 10 tips (each in 7 shades) (0.16mL per tip)
X-XBW, X-BW, X-B1, X-B2, X-C2, X-D2, X-AO2
Note:
Weight per Unitip : 0.30g
3. Option
a. Unitip APPLIER b. Shade guide c. Mixing pad (No.14B)

- Syringes
1. Syringe (in 12 shades) (2.7mL per syringe)
Note:
Weight per syringe : 5.0g
2. Option
a. Shade guide b. Mixing pad (No.14B)
Note:
All packaging types may not be available in all countries.

- CAUTION
1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Take care to avoid ingestion of the material.
4. Wear plastic or rubber gloves during operation to avoid direct contact with air insoluble resin layers in order to prevent possible sensitivity.
5. For infection control reasons, Unitips are for single use only.
6. Wear protective eyeglasses during light curing.
7. When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
8. Do not mix with other similar products.
9. Avoid getting material on clothing.
10. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
11. Do not use GRADIA DIRECT X in combination with eugenol containing materials as eugenol may hinder GRADIA DIRECT X from setting.
12. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
13. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

- Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<http://www.gceurope.com>
or for The Americas
<http://www.gcamerica.com>
They can also be obtained from your supplier.

- CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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