GC KALORE

LIGHT-CURED RADIOPAQUE UNIVERSAL COMPOSITE RESTORATIVE

For use only by a dental professional in the recommended indications.

- RECOMMENDED INDICATIONS**
 - Direct restorative for Class I, II, III, IV cavities.
 - Direct restorative for wedge-shaped defects and root surface cavities.
 - Direct restorative for veneers and diastema closure.

- CONTRAINDICATIONS**
 - Pulp capping.
 - Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.


DIRECTIONS FOR USE

- Shade Selection
Prepare the tooth with pumice and water. Shade selection should be made prior to isolation. Select the appropriate shades by referring to the GC KALORE Shade Guide or Multi Shade Build-Up Guide.
- Cavity Preparation
Prepare the cavity using standard techniques. Dry by gently blowing with oil free air. Note:
 - For pulp capping, use calcium hydroxide.
- Bonding Treatment
For bonding GC KALORE to enamel and / or dentin, use a light-cured bonding system such as G-BOND, G-aenial Bond or G-Premio BOND (Fig. 1). Follow manufacturer's instructions.
- Placement of GC KALORE
 - Dispensing from a Unipit
Insert the GC KALORE Unipit into a commercially available applicator (Unipit APPLIER II is recommended). Refer to the applicator manufacturer's instructions for use. Remove the cap and extrude material directly into the prepared cavity. Use steady pressure (Fig. 2).
 - Remove syringe
Remove syringe cap and dispense material onto a mixing pad. Place the material into the cavity using a suitable placement instrument. After dispensing, screw syringe cap clockwise until it is tight to fully turn to release residual pressure inside the syringe. Replace cap immediately after use.
- Material can be applied in a single shade layer to achieve aesthetic restorations using Universal shades. For details, refer to the Clinical Hints.
- Material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at normal room temperature.
- After dispensing, minimize exposure to ambient light. Ambient light can shorten the manipulation time.

Clinical Hints

- Finishing small cavities
Restore using a one shade technique. In most cases the use of one Universal shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Translucent shades can be selected.
- In the case of large and/or deep cavities
In most cases a multi shade layering technique will give the best aesthetic results. To block out shine through from the oral cavity or to mask discolorated dentin, use an appropriate Opaque shade and continue to build up with a Universal shade. For optimal aesthetics use a Translucent shade as the final composite layer.
- In the case of deep posterior cavities, a flowable composite such as GRADIA DIRECT Fio / LoFio or G-aenial Fio / G-aenial Fio X / Universal Fio or a glass ionomer cement such as GC Fuji IX, GC Fuji II LC or GC Fuji TRIAGE can be used on the cavity floor instead of an Opaque shade. See also Examples of Clinical Applications and/or Shade Combination Chart.

Examples of Clinical Applications (Clinical Hint)

GC KALORE shade combination chart for multiple layers in deep and/or large cavities

	A1	A2	A3	A3.5	A4	B1	B2	B3	C2	C3	D2	CV	CVD	BW	XBW	
Opaque	OBW	AO2	AO3	AO4	AO4	OBW	AO2	AO3	AO4	AO4	OBW	AO2	AO3	AO4	OBW	OXBW
Universal	A1	A2	A3	A3.5	A4	B1	B2	B3	C2	C3	D2	CV	CVD	BW	XBW	
Translucent (Enamel)	WT	WT	DT	DT	DT	WT	DT	DT	DT	DT	WT	DT	DT	WT	WT	

For details of shades, refer to the following section of SHADES.

- Contouring before Light Curing
Contour using standard techniques.
- Light Curing
Light cure GC KALORE using a light curing unit (Fig. 3). Keep light guide as close as possible to the surface. Refer to the following chart for irradiation Time and Effective Depth of Cure.

GC KALORE : Irradiation Time and Effective Depth of Cure

	Irradiation time	Plasma arc (2000mW/cm ²)	3 sec.	6 sec.
Shade	High power LED (more than 1200mW/cm ²)	Halogen / LED (700mW/cm ²)	10 sec.	20 sec.
			20 sec.	40 sec.
CT, NT, WT, GT, CVT			3.0 mm	3.5 mm
A1, A2, B1, B2, C2, D2, XBW, BW, DT			2.5 mm	3.0 mm
A3, A3.5, A4			2.0 mm	3.0 mm
A4, C3, AO2, AO3, AO4, CV, CVD, OBW, OXBW			1.5 mm	2.5 mm

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STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F) away from high temperatures or direct sunlight. (Shelf life: 3 years from date of manufacture)

PACKAGES

- Unipits
a. Pack of 20 tips (each in 11 shades) (0.16mL per tip)
A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D2
b. Set of 10 tips (each in 15 shades) (0.16mL per tip)
XBW, BW, CV, CVD, AO2, AO3, AO4, OBW, OXBW, WT, DT, CT, NT, GT, CVT
Note:
Weight per Unipit : 0.3g
- Option
a. Shade guide b. Mixing pad (No.14B)
- Syringes
a. Syringe (in 26 shades) (2.0mL per syringe)
Note:
Weight per syringe : 4g
- Option
a. Shade guide b. Mixing pad (No.14B)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingestion of the material.
- Prior to starting the filling process, partially extrude paste outside of the patient's mouth and away from the patient to ensure paste is flowing properly.
- Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers in order to prevent possible sensitivity.
- For infection control reasons, Unipits are for single use only.
- Wear protective eye glasses during light curing.
- When polishing the polymerized material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
- Do not mix with other similar products.
- Avoid getting material on clothing.
- In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge, or cotton pellet before light curing.
- Do not use GC KALORE in combination with eugenol containing materials as eugenol may hinder GC KALORE from setting.
- Personal protective equipment (PPE) such as gloves, face masks and safety glasses should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <http://www.gcgpore.com> or for The Americas <http://www.gcamericas.com>

THEY CAN ALSO BE OBTAINED FROM YOUR SUPPLIER.

CLEANING AND DISINFECTING

USE OF DELIVERED SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSER. Immerse clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional guidelines.

Last revised: 11/2017

© 2017 GC Corporation

GC CORPORATION

GC AMERICA INC.

GC EUROPE N.V.

GC DENTAL PRODUCTS CORP.

RESPONSIBLE MANUFACTURER IN CANADA:

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