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Light-cured resin for temporary crown, bridge, inlay & onlay  
REVOITEK LC is a new light-cure single-component sculptable composite resin for temporary inlays, onlays, crowns and bridges.

For use only by a dental professional in the recommended indications.  
**RECOMMENDED INDICATIONS**  
Fabrication of temporary crowns, bridges, inlays and onlays.

### CONTRAINDICATIONS

In rare cases the product may cause sensitivity in some people. If such any reactions are experienced, discontinue the use of the product and refer to a physician.  
**DIRECTIONS FOR USE**  
I. DIRECT METHOD (on the patient's mouth)  
A. Temporary Crowns & Bridges

1. Transfer of putty stick to the storage case  
Transfer the putty stick from the aluminum foil pack into the storage case provided in the Introductory Set.  
Note:  
1) When putty stick is transferred to the storage case or is dispensed, handle with care in order to prevent contamination.

2) If the storage case is used repeatedly, fine abrasion fragments may be produced from the lid because of the structure of the case. When the lid is replaced, make sure that the inside of the case is clean. If it is dirty, clean with cotton moistened with alcohol.

### Dispensing material

Dispense the required amount of material using the spatula. Adjust the shape of the margins in such a way so that easy placement onto the abutment. To soften the material, lightly knead several times with the fingers.

Note:  
1) Do not warm up before use. The material will become sticky, making it much more difficult to handle.  
2) Do not knead excessively, or the material will become sticky.  
3) After dispensing, immediately close the storage case to protect the material from light.

3. Preparations for pressing material onto the abutment  
Press the putty stick onto the abutment roughly contour with the fingers or spatula.

Note:  
1) Apply GC COCOA BUTTER or Vaseline to the fingers or spatula to help shape the material more easily and to make the surfaces glossier. For additional application of material, temporarily light cure and remove the coated surface using a laboratory carbide bur.  
2) Any material remaining after shaping should not be returned for storage because it has been exposed to ambient light. Dispose of it according to normal practice.

5. Shaping material - Step 1  
Let the patient bite softly to register the occlusal surface and adjust occlusion. Also adjust and contour the margin of the buccal surface.

6. Shaping material - Step 2  
Adjust the margins as necessary, and contour the proximal and lingual surfaces. Using an appropriate instrument, remove any excess material, particularly from interproximal spaces.

Note:  
1) If shaping the material takes a long time, it may become sticky. In this case, apply GC COCOA BUTTER or Vaseline to the fingers or instrument. In cases where additional application of material is needed, temporarily light cure and remove the coated surface using a laboratory carbide bur.

2) Any material remaining after shaping should not be returned for storage because it has been exposed to ambient light. Dispose of it according to normal practice.

7. Temporary light curing  
Temporarily light cure the restoration in the patient's mouth in order to prevent possible deformation during removal from the mouth. Light cure all surfaces of the crown or unit for a combined total of 10 seconds with Halogen/LED or 3-5 seconds with a plasma arc.

Note:  
1) Make sure that no excess material remains in interproximal spaces before light curing. Any material remaining will make it difficult to remove the temporary restoration from the mouth.

2) Remove the trial fit temporarily light cured restoration perpendicular to the abutment.

8. Final light curing  
Perform final polymerization outside the patient's mouth. Light cure each of the buccal, proximal, occlusal and lingual surfaces for 20 seconds with Halogen / LED or 3-5 seconds with a plasma arc. When using a tabletop fluorescent light curing unit (GC LABOLIGHT LV-III or equivalent), light cure for at least 3 minutes.

Note:  
1) When using a hand held type of light curing unit, be sure to light cure the internal surface of each crown or unit for the above designated time.  
2) When using a hand held type of light curing unit, a range of effective polymerization depends on the cross sectional area of the light guide. If the resin restoration is larger than that, light cure in sections for effective polymerization.

3) Too short an irradiation time will result in incomplete hardening and possible discoloration.

Irradiation time and Depth of cure

	1mm	2mm	3mm	4mm
Halogen / LED	-	-	10 sec.	30 sec.
Plasma Arc	-	-	3 sec.	5 sec.
G-Light	-	5 sec.	10 sec.	20 sec.
GC LABOLIGHT LV-III	30 sec.	1 min.	3 min.	15 min.

Note: 20 sec. with Halogen light - 3.5mm

9. Correction of contour, finishing and polishing  
Correct the temporary restoration using a laboratory carbide bur or silicone point. For an external addition, use REVOITEK LC. For an internal addition, use a self-curing resin (GC UNIFAST Trad, ALIKE), a light-cured resin (GC UNIFAST LO) or a flowable composite (GC UNIFAST LC) or a flowable composite with a lot or charms which to provide a glossy and beautiful appearance.

Note:  
1) It is recommended to use a laboratory carbide bur to grind the hardened restoration. When using a fine diamond bur such as GC SMOOTH CUT, or a steel bur, first remove the unpolymerized layer with a laboratory carbide bur or alcohol to prevent the diamond or steel bur from chipping.

2) Before GC UNIFAST Trad, ALIKE or GC UNIFAST LC is additionally poured onto the inside of the restoration, drill a hole into the occlusal or buccal surface for the flow of excess resin so that the temporary restoration can be properly seated.

10. Completion  
Cement the completed temporary restoration using GC FREEGENOL TEMPORARY PACK or equivalent temporary cement.

11. Repair  
If needed, the cemented temporary restoration can be repaired in the patient's mouth using REVOITEK LC or a self-curing resin (GC UNIFAST Trad) or a flowable composite (GC UNIFAST LC) or a flowable composite according to its instructions for use.

a. Remove a layer of the surface to be repaired using a laboratory carbide bur. For optimal results, apply a resin bonding agent (GC RELINE HARD BONDING AGENT) to the area and dry with an air syringe.

b. Place material under the area, shape, cure and polish.

### Temporaries Inlays

1. Transfer of putty stick to the storage case  
Follow instructions as described under A1 above.

2. Dispensing material  
Follow instructions as described under A2 above.

3. Preparations for applying material into the inlay cavity  
Block out the occlusal surface with a V.L.C. or composite filling agent. Apply GC COCOA BUTTER or Vaseline to all surfaces of the preparation, including the V.L.C. glass ionomers base or V.L.C.

4. Pressing the material into the cavity and roughly contour with the fingers or spatula.  
Note:  
1) Apply GC COCOA BUTTER or Vaseline to the fingers or spatula to help shape the material more easily and to make the surfaces glossier. For additional application of material, temporarily light cure and remove the coated surface using a laboratory carbide bur, etc.

2) Let the patient bite softly to register the occlusal surface and adjust occlusion.

5. Shaping material - Step 1  
Adjust the margins as necessary. Using an appropriate instrument, remove any excess material.

6. Shaping material - Step 2  
Adjust the margins as necessary. Using an appropriate instrument, remove any excess material.

Note:  
1) If shaping the material takes a long time, it may become sticky. In this case, apply GC COCOA BUTTER or Vaseline to the fingers or instrument. In cases where additional application of material is needed, temporarily light cure and remove the coated surface using a laboratory carbide bur, etc.

2) Any material remaining after shaping should not be returned for storage because it has been exposed to ambient light. Dispose of it according to normal practice.

3) For ease of removal, you can insert a pin into the inlay.

7. Temporary light curing  
Temporarily light cure the restoration in the patient's mouth in order to prevent possible deformation during removal from the mouth. For curing time recommendations and notes see above under A7.

8. Final light curing  
Perform final polymerization outside the patient's mouth. For curing time recommendations and notes see above under A8.

9. Correction of contour, finishing and polishing  
For instructions and notes see above under A9.

10. Completion  
Cement the completed temporary restoration using GC FREEGENOL TEMPORARY PACK or equivalent temporary cement.

11. Repair  
For instructions see above under A11.

### II. INDIRECT METHOD (on the stone model)

1. When a REVOITEK LC restoration is made on the stone model, apply a small amount of GC COCOA BUTTER or Vaseline to the area which will come in contact with the patient's mouth.

2. Follow instructions described in DIRECT METHOD.

### III. COMBINED USE WITH READY-MADE PLASTIC AESTHETIC TEMPORARY CROWNS

1. When a plastic full temporary crown is used  
a. Apply a resin bonding agent (GC RELINE HARD BONDING AGENT) to the internal surface of the plastic full temporary crown.  
Note:  
More stable adhesion can be obtained by roughening the internal surface of the full temporary crown with a paper coarse shaped point or other instrument before application of the resin bonding agent.

b. Fill the temporary crown with REVOITEK LC and press it onto the abutment. c. Apply light to the outside of the crown for temporary polymerization as described in DIRECT METHOD, A7 Temporary light curing.

d. After removing the crown from the patient's mouth, follow instructions described in DIRECT METHOD, A8 concerning final light curing, correction of contour, finishing, polishing and completion.

### IV. COMBINED USE WITH READY-MADE TRANSPARENT PLASTIC STRIP CROWN FORMS

1. Adjust the cervical area of the plastic strip crown form using crown scissors and bars, and make an escape hole in the occlusal or buccal surface. Trial fit it in the patient's mouth.

2. Lightly knead REVOITEK LC several times. Fill the strip crown form with material and press it to the abutment. Remove excess material.

3. For temporary polymerization, apply light to all surfaces of the crown form as described in DIRECT METHOD, A7 Temporary light curing.

4. After removing the crown from the patient's mouth, for final polymerization apply light to each surface as described in DIRECT METHOD, A8 Final light curing. Remove the crown form using a sharp instrument, trim and finish temporary restoration. Seat in the patient's mouth using the temporary cement.

### STORAGE

Store in a cool and dark place away from high temperature and direct sunlight or intense ambient light.  
(Shelf life: 2 years from date of manufacture)

### SHADE

52 (based on Vita® shades)  
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany

### PACKAGE

1. Introductory Set  
Pack of putty stick 16g, storage case (1), plastic spatula No. 2 (1)

2. Refill  
3 packs of putty stick 16g

### CAUTION

1. In case of contact with eyes, flush immediately with water and refer to a physician.

2. Do not use or store near fire.

3. Take care not to swallow the paste.

4. Do not mix with other materials.

5. When using a light curing unit, wear eye protection glasses. Do not look directly into the light.

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