Prior to use, carefully read the instructions for use. **EN**

# UNIFAST™ LC


LIGHT CURED ACRYLIC RESIN

For use only by a dental professional in the recommended indications.

- RECOMMENDED INDICATIONS**  
Direct technique using a hand held visible light curing device (470 nm wavelength).  
Indirect technique using a table top curing unit (470 nm wavelength).  
By pressure and/or paint on techniques.  
1. Temporary crowns, bridges, inlays and onlays.  
2. Fixing the position of metal cast crowns for soldering.  
3. Repair of dentures.  
4. Individual tooth trays, individual trays.  
5. Orthodontic plates.  
6. Check bite, splint, rebase, night guard, etc.

**CONTRAINDICATIONS**  
Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

- DIRECTIONS FOR USE**  
a) The standard powder/liquid ratio is 1.0g/0.5mL (first graduation on the powder measure/21 drops of liquid).  
b) Add powder to liquid in the mixing cup and mix quickly for 10-15 seconds. The following illustrates the lapsed time stages in polymerisation:

- c) Use the slurry or dough according to application standard techniques.  
d) When all necessary corrections or adjustments have been made, light cure the resin using an appropriate visible light curing device (470 nm wavelength).  
Note:  
1. Quality of cure, and depth of cure, vary depending on the light source used, and on the translucency or opacity of the shade of resin.

	Depth of cure (mm)			
	Translucent	Ivory	Pink # 3/E	Live Pink # 8
Hand held: 20 seconds 60 seconds	6-9 13-16	3-4 5-6	2-3 3-4	3-4 5-7
Table top: 1 minute 10 minutes	10 >25	3 16	2 16	5 24

- The above are typical values and are given for guidance only.  
2. Light cure within 1 hour from starting the mix.  
3. UNIFAST LC can also be cured chemically in 5 minutes in water at 60°C (140°F).  
4. Do not use powder or liquid with components of other brands of acrylic resin.

**STORAGE**  
Recommended for optimal performance, store in a cool dark place away from direct sunlight (4-25°C) (39.2-77.0°F).

- PACKAGES**  
1. Bottle of 50g powder.  
2. Box of 30g liquid (2 x 15g or 14.3mL).  
3. 1-2 package\* : 1 x 50g powder, 2 x 15g (14.3mL) liquid, accessories.  
\* Not available in Europe.

- SHADES**  
1. Crown colour Ivory\*, A2, A3, B2, B3, C2  
2. Gingival colour\* Pink # 3, Live Pink # 8 (veined), Pink E  
3. Translucent Translucent  
\* Not available in Europe.

- CAUTION**  
1. The liquid is volatile and highly flammable. Close the bottle immediately after use and keep away from sources of ignition.  
2. In case of contact with eyes, flush immediately with water and seek medical attention.  
3. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.  
4. In rare cases the product may cause sensitivity in some people.  
If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:  
<http://www.gceurope.com> or for The Americas  
<http://www.gcamerica.com>  
They can also be obtained from your supplier.

**CLEANING AND DISINFECTING**  
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.  
DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised : 04/2018

MANUFACTURED by  
GC DENTAL PRODUCTS CORP.  
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

RESPONSIBLE MANUFACTURER IN CANADA:  
GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900

DISTRIBUTED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8565, Japan

EU : GC EUROPE N.V.  
Researchpark Haasrode-Luven 1240, Interleuvenlaan 33,  
B-3001 Luven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900  
www.gcamerica.com

GC SOUTH AMERICA  
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL  
CEP: 02022-051 - TEL: +55-11-2925-0965  
CNPJ: 08.279.999/0001-61  
RESP. TEC: Mayara de Santis Ribeiro - CRO/SP 105.982

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