Prior to use, carefully read the instructions for use.

EN

UNIFAST LC

LIGHT CURED ACRYLIC RESIN

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

Direct technique using a hand held visible light curing device (470 nm wavelength).

Indirect technique using a table top curing unit (470 nm wavelength).
By pressure and/or paint on techniques.


- Temporary crowns, bridges, inlays and onlays.
- Fixing the position of metal cast crowns for soldering.
- Repair of dentures.
- Individual tooth trays, individual trays.
- Orthodontic plates.
- Check bite, splint, rebase, night guard, etc.

CONTRAINDICATIONS

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

- The standard powder/liquid ratio is 1.0g/0.5mL (first graduation on the powder measure/21 drops of liquid).
- Add powder to liquid in the mixing cup and mix quickly for 10-15 seconds. The following illustrates the lapsed time stages in polymerisation :

- Use the slurry or dough according to application using standard techniques.

- When all necessary corrections or adjustments have been made, light cure the resin using an appropriate visible light curing device (470 nm wavelength).

Note :

- Quality of cure, and depth of cure, vary depending on the light source used, and on the translucency or opacity of the shade of resin.

	Depth of cure (mm)			
	Translucent	Ivory	Pink # 3/E	Live Pink # 8
Hand held: <p>20 seconds</p> <p>60 seconds</p>	6-9	3-4	2-3	3-4
	13-16	5-6	3-4	5-7
Table top: <p>1 minute</p> <p>10 minutes</p>	10	3	2	5
	>25	16	16	24

The above are typical values and are given for guidance only.

- Light cure within 1 hour from starting the mix.
- GC UNIFAST LC can also be cured chemically in 5 minutes in water at 60°C (140°F).
- Do not use powder or liquid with components of other brands of acrylic resin.

STORAGE

Store in a cool dark place away from direct sunlight.
(Shelf life : 3 years from date of manufacture)

PACKAGES

- Bottle of 50g powder.
- Box of 30g liquid (2 x 15g or 14.3mL).
- 1-2 package* : 1 x 50g powder, 2 x 15g (14.3mL) liquid, accessories.

* Not available in Europe.

SHADES

- Crown colour Ivory*, A2, A3, B2, B3, C2
- Gingival colour* Pink # 3, Live Pink # 8 (veined), Pink E
- Translucent Translucent

* Not available in Europe.

CAUTION

- The liquid is volatile and highly flammable. Close the bottle immediately after use and keep away from sources of ignition.
- In case of contact with eyes, flush immediately with water and seek medical attention.

Last revised : 04/2009

CE 0086

'GC'

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900.

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6546 7588

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