Shade	A2, A3, B2, B3, C2	Translucent
Time (sec)	20	40
Depth (mm) (avg)	3.4	5.4

GC AMERICA INC.
ALSIP, IL 60803

MANUFACTURED by:
GC CORPORATION
76-1 HASUNUMA-CHO, ITABASHI-KU, TOKYO JAPAN

DISTRIBUTED by:
GC AMERICA INC.
ALSIP, IL 60803 U.S.A.

3910900

UNIFAST™ LC

LIGHT QUICK-CURING ACRYLIC RESIN

GC UNIFAST LC gradually changes in consistency from slurry to dough-like to rubber-like. Outstanding elasticity in its rubber-like phase enables easy removal of undercuts. Also, prior to light irradiation undercuts or excess material can be trimmed easily using a scissors or cutting knife.

APPLICATIONS

- Temporary Crowns, bridges, inlays and onlays.
- Fixing the position of metal cast crowns for soldering.
- Individual tooth trays, individual trays.
- Orthodontic plates.
- Checkbite, bite plate (splint), rebase, night guard.

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
CURING TIME AND DEPTH

Curing time and depth of cure will vary depending on the brand and type of light being used. Below are average curing times and depths.

Hand Held Curing Lights.

Shade	A2, A3, B2, B3, C2			Translucent		
	20	40	60	20	40	60
Time (sec)	20	40	60	20	40	60
Depth (mm) (avg)	3.4	4.5	5.4	7.6	11.7	15.2

NOTES:

1. Photocure the mixed resin within one hour from the start of mix to avoid deterioration of its original efficiency. The mixture will change the balance of its composition as time passes even without irradiation of light.
2. GC UNIFAST LC can also be fully cured in 5 minutes in water (60°C).
3. Non-eugenol type temporary cement such as GC FREEGENOL TEMPORARY PACK or NOGENOL TEMPORARY CEMENT is recommended for temporary cementation.
4. Since the liquid is flammable, keep it away from fire. Use in a well ventilated place. If it should spill on table or floor, wipe off with a dry cloth.
5. The liquid bottle should be tightly closed after use and stored in a cool and dark place, away from direct sunlight.
6. Neither the liquid nor powder component of UNIFAST LC should be mixed with any other brand of self-curing acrylic resin.
7. Handle the liquid with care and avoid direct contact with the skin. In case of contact with the skin, immediately wash off with soap and water. The liquid may cause skin irritation.
8. Do not return powder and liquid left in the rubber cup to their containers in order to keep good quality of the contents.

SHADES

A2, A3, B2, B3, C2, Translucent

PACKAGES

1. 1 Bottle (30g) each - A2-A3-B2-B3-C2- Translucent Powder
2. 2 Bottles (15g each) Liquid
3. 2 Rubber Mixing Bowls, 1 Brush, 1 Spatula, 2 Measuring Vials (1 Liquid - 1 Powder)

UNIFAST LC

• TECHNIQUES

1. Direct - using a gun type light curing unit.
2. Indirect - using a table top light curing unit.
3. Pressure and/or brush on technique.

• APPLICATIONS


1. Temporary crowns, bridges, inlays and onlays.
2. Fixing the position of metal cast crowns for soldering.
3. Individual tooth trays, individual trays.
4. Orthodontic plates.
5. Checkbite, bite plate (splint), rebase, night guard.

• DIRECTIONS FOR USE

Take an impression using alginate or silicone impression material prior to cutting the preparation. Follow the directions below.

1. For smaller restorations such as inlays, onlays or single crowns use the standard powder/liquid ratio (1g powder to .5ml of liquid). This will be the first gradation on both the powder and liquid measuring vials.
2. Mix quickly for 10-15 seconds with the plastic spatula provided.
3. Pour Unifast LC mixture into the impression (be sure to cut relief slots or grooves to allow excess material to escape).
4. Seat the impression back into its original position. Any excess mixture will escape through the relief slots.
5. Hold the impression in place for two minutes.
6. Remove the impression from the mouth. Carefully remove the rubber-like Unifast from the impression, trim the flash and replace the flexible Unifast back in the mouth.
7. Partially photocure the resin while having the patient bite gently into occlusion. This assures accurate occlusion and eliminates any excessive grinding while finishing. Cure until the material begins to get firm.
8. Remove the partially cured material from the mouth and cure any area that may have been unexposed to the curing light. This can be done with either a hand held curing light or a table top light box.
9. Adjust occlusion if needed, shape and finish in the usual manner.

• ELAPSED TIME AND STAGES IN POLYMERIZATION

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