## General Information, Precautions & Cleaning Instructions

### The Ultimate F&P Kit:

The carbide burs and polishing instruments of the **Ultimate F&P Kit** are shipped in non-sterile packaging and **MUST** be sterilized prior to use. **PLEASE NOTE:** These rotary instruments are capable of multiple use and **MUST** be cleaned and sterilized between each clinical patient application. Please read this instruction sheet in its entirety prior to proceeding.

### General Product Description:

The **carbide burs** of the Ultimate F&P Kit are manufactured as "one piece," tungsten carbide instruments - which ensure better concentricity when compared to "two piece" instruments (made up of shank and cutting end). These rotary instruments are designed to be used with industry standard, high-speed handpieces with "push button", auto-chucking mechanism. These carbide burs can be used to shape, contour and create composite restorations in both the anterior and posterior region of a patient's dentition. Recommended speed for carbide burs is **150,000-200,000 RPM**.

The **polishing instruments** of the Ultimate F&P Kit are made with high-quality, stainless steel shanks and diamond-impregnated rubber material available in disc, cup and wheel shapes for all types of polishing applications. These polishing instruments are designed to be used with industry standard, slow-speed handpieces with RA (right angle) latch-type chucking mechanism. Recommended speed for the polishing instruments is **5,000-7,500 RMP**.

### Clinical Application:

- The rotary instruments of the Ultimate F&P Kit are only to be used by a licensed dentist where patients are involved. These instruments may be used by a dental laboratory technicians, dental students and/or dental auxiliary "only" where no patient is involved. Any other use of these products outside of the aforementioned settings is strictly prohibited.
- These instruments may be used by a dental student (on a live patient) only under the direct supervision of a license dental instructor in a clinical teaching setting.

### Intended Use:

- These rotary instruments should only be used with properly maintained/functioning high-speed and low-speed handpieces following the manufacturer's guidelines and indication for use. Never use any rotary instrument beyond the suggested RPM's (listed above) or outside of handpiece manufacturer's indications for use.
- Using a carbide bur and/or polishing piece at higher RPM may generate unnecessary heat, fatigue/breakage of the instrument leading to patient discomfort and/or injury.
- Proper irrigation should always be used when operating these, or any rotary instruments. Lack of proper irrigation during operation of these rotary instruments may result in the generation of unnecessary heat that may generate unnecessary heat, fatigue/breakage of the instrument leading to patient discomfort and/or injury.

### Precautions:

- The Tungsten Carbide Burs within this kit will last for approximately **15 uses**; however, general wear is dependent upon hand pressure, duration of clinical procedure, air pressure, RPM's, etc. Additionally, the Diamond Impregnated Polishing Instruments will last for approximately **10 uses** and general wear may vary with the aforementioned conditions: hand pressure, duration of procedure, air pressure, RPM's, etc.

- **Do not** apply excessive hand pressure on the handpiece or bur during operation as this may generate unnecessary heat, fatigue/breakage of the instrument leading to patient discomfort and/or injury.
- Properly clean and sterilize all carbide burs and polishing instruments prior to any use.
- **Do not** use chemical or dry heat to sterilize these rotary instruments. (Please see sterilization procedure below).
- Proper cleaning and sterilization is required after every use of every rotary instrument to prevent cross-contamination.
- Failure to properly remove any composite debris from a carbide bur may cause the instrument to fatigue and break during use, which may result in the generation of unnecessary heat that could cause necrosis of both hard & soft tissue and severe patient discomfort.
- It is recommended to use a rubber dental dam when using any rotary instrument to prevent accidental aspiration or swallowing of the bur or polishing instrument.

### General Operating Instructions:

Inspect every carbide bur and polishing instrument prior to any use:

- The use of worn-out or dull carbide burs may generate unnecessary heat, fatigue/breakage of the instrument leading to patient discomfort and/or injury.
- The use of bent carbide burs will cause unwanted rotational vibration that may generate unnecessary heat, fatigue/breakage of the instrument leading to patient discomfort and/or injury.
- During clinical application, be sure to move the bur in a continuous motion to avoid localized heating and/or damage caused from unnecessary heat, fatigue/breakage of the instrument.
- Use caution when removing the handpiece and carbide bur from the patient's oral environment. The sharp edges of the carbide burs can cut both patient and operator leading to irritation and injury.
- Be sure all handpieces are in "good working condition" to ensure maximum effectiveness of the bur and dental procedure.
- Failure to properly maintain handpieces may lead to clinical failures: delay in procedure, injury to the patient or user, aspiration or swallowing of rotary instrument or damage to the preparation site.
- Always ensure the carbide bur or polishing instrument is fully seated into the handpiece's chucking mechanism (for both high-speed and low-speed handpieces). Failure to properly seat a rotary instrument may lead to ejection (also called "walk out") of the carbide bur or polishing instrument causing injury to the patient or user.
- Never force a rotary instrument into a handpiece's chucking mechanism as this could cause damage to the handpiece's collet, which could result instrument "walk out" and other handpiece damage.

### General Warnings & Cautions:

1. **Do not** use cleaning agents with **chlorine** or **chloride** as the active ingredient when cleaning any GC America rotary instrument as these agents are corrosive to stainless steel and must never be used. Always use cleaning agents with **neutral pH** balances.
2. **Do not** use Cold Sterilizing Methods (including Glutaraldehyde, Activated Dialdehyde, Ortho-Phthaladehyde, etc.) for the treatment of carbide burs or polishing instruments. These agents often contain strong oxidizing chemicals that may adversely affect the properties of these rotary instruments.
3. Dental care providers should always **wear a protective mask** when working with rotary instruments to avoid unnecessary inhalation of aerosol and/or hard/soft tissue pieces and particles.
4. **Always** read and review package labels and instructions to ensure use of the appropriate device for the desired clinical application. Failure to follow instructions-for-use may cause the following: preparation site damage, injury to the patient or user, or possible aspiration or swallowing of the rotary instrument.

**PLEASE NOTE:** Tungsten Carbide Burs contain nickel and should not be used on individuals who have metallurgy and/or hypersensitivity to this substrate. Additionally, the State of California requires that any product that contains nickel must be identified as an agent that may cause cancer, birth defects or other reproductive harm.

### General Cleaning & Sterilizing Instructions:

Always wear protective gloves when handling contaminated instruments to avoid possible infection/cross-contamination. These instructions are intended to be followed regarding GC America carbide dental burs and polishing instruments. They are applicable before initial use and after each subsequent use.

**PLEASE NOTE:** All GC America carbide burs and polishing instruments are provided in clean condition, but are not sterilized. **All GC America rotary instruments should be sterilized before any/all use.** All GC America carbide burs and polishing instruments have been manufactured to specific tolerances. When placing them back into the Bur Block, be careful not to force or misdirected the "shank" portion into the positioning hole. If an instrument "does not" fit easily into a positioning hole, be sure to inspect to insure no debris or material is lodged or attached to the shank (preventing the instrument from seating back into the block).

**Always wear proper infection control barriers (gloves, mask, eye protection, etc.) during any instrument cleaning process.**

- Always inspect all carbide burs and polishing instruments prior to any use and during the cleaning process.
- Carbide burs can be handled "wet" or "dry" and should always be protected from damage by comingling with other instruments. Do not allow these rotary instruments to remain "wet" for long periods of time as this may increase the chance of staining or corrosion. Additionally, avoid prolonged exposure to disinfectant solutions which may result in damage to the instrument metal and/or rubber.
- Always clean any rotary instrument in a designated infection control sink reserved for cleaning instruments with a designated instrument brush.
- After cleaning off loose debris, replace the rotary instruments back into the Bur Block and rinse it under cold water for at least one (1) minute.
- Prepare a fresh batch of neutral-pH cleaning solution in a clean and appropriate tub. Following the manufacturer's instructions, immerse the Bur Block (containing both the carbide burs and polishing instruments) into the bath for at least ten (10) minutes.
- After this ten (10) minute immersion, brush off any remaining debris into a sink or appropriate receptacle using a neutral cleaning agent - being careful not to get any material/liquid on you. Brush and clean for approximate one (1) minute. Care should be taken to avoid spreading contaminant debris onto you, others or general surfaces. Always use caution during this process as any brushing motion can cause particles to spray or project in a variety of directions.
- After this cleaning process, visually inspect each instrument to ensure the removal of general debris. Repeat this series of steps if necessary.
- Thoroughly rinse the Bur Block under warm, running water for at least one (1) minute and until all cleaning solution is visibly washed away.
- Dry the Bur Block using a dry wipe and/or clean, compressed air.

### Ultrasonic Cleaning Instructions:

- Prepare a fresh batch of pH-neutral cleaning solution and submerge the Bur Block into Ultrasonic

bath following the manufacturers' instructions.

- Expose the submerged Bur Block for at least fifteen (15) minutes. Remove and thoroughly rinse the Bur Block under warm, running for at least (1) minute.
- Visually inspect to confirm the removal of any/all debris and repeat the cycle if needed. Dry the Bur Block using a dry wipe and/or clean, compressed air.
- Visually inspect each rotary instrument for damage and/or wear that would prevent proper operation.
  - Do not use if any part of the rotary instrument is broken.
  - Do not use if any part of the rotary instrument shows signs of corrosion or unusual wear.

**Sterilization Guidelines:**

Use the following cycle for sterilization:

<b>Cycle Type:</b>	<b>Minimum Exposure Time</b>	<b>Minimum Exposure Temperature</b>	<b>Minimum Dry Time</b>
Gravity	10 Minutes	135°C (275°F)	30 Minutes
Pre-Vacuum (4 pulses)	3 Minutes	134°C (273°F)	30 Minutes

Please be sure that the maximum load for the sterilizer has not been exceeded. Please be sure to allow for ample “dry time” and that the instruments are never left “wet”. Leaving rotary instruments in wet condition may cause corrosion and eventual failure.

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