SALIVA-CHECK MUTANS

Streptococcus mutans RAPID DETECTION KIT
For use only by a dental professional in the recommended indications.

The bacterial species *Streptococcus mutans* have a significant role in the initiation of dental caries. Salivary level of this bacterium is useful for caries risk assessment in individual patients, when used in conjunction with other clinical information. The SALIVA-CHECK MUTANS test kit provides a semi-quantitative evaluation of the level of *Streptococcus mutans* in the saliva in 15 minutes by using monoclonal antibodies.

RECOMMENDED INDICATIONS
For detection of *Streptococcus mutans* level in saliva.

- FEATURES**
- Patients with high levels (>5x10⁵ CFU/mL saliva) of *Streptococcus mutans* can be quickly and easily identified within 15 minutes without any specialized equipment and without cultivation of the bacteria.
 - SALIVA-CHECK MUTANS detects *Streptococcus mutans* in saliva using a highly specific immunochromatography process and does not rely upon growth of bacteria. Therefore problems such as contamination by other bacteria or pathogens are avoided.
 - The test result is shown visually with a line, so outcome can be easily determined.
 - By using a swab (available separately), a saliva sample can be collected from a patient (such as infants or medically handicapped patients) who has difficulty expectorating in the saliva mixing container.

- PRECAUTIONS**
- The test device, paraffin gum, pipette and mixing container are for single use only.
 - It is also recommended to use a plaque disclosing agent and/or a disclosing tablet to gather a caries risk profile and to educate the patient about their caries risk status.
 - Do not use a plaque disclosing agent before collecting saliva sample for the SALIVA-CHECK MUTANS test kit because of the risk of false results.
 - If the saliva is contaminated with the plaque disclosing agent, the saliva will not change to a light green colour after adding the reagent #1 and #2.
 - The patient rinses the mouth to remove the disclosing agent. *Streptococcus mutans* may also be removed, and the test result may not be accurate.
 - Prior to any visit where a saliva diagnostic procedure is planned, instruct the patient to avoid tooth brushing, flossing, or brushing the teeth or use a mouth wash for at least one hour prior to the scheduled appointment time.

- DIRECTIONS FOR USE**
- Open the foil package and take out the test device, paraffin gum, mixing container and pipette.
 - Instruct the patient to chew the gum provided for 1 minute to stimulate their secretion of saliva.
- NOTE:**
Gum should not be swallowed.
Collect the patient's saliva sample in the mixing container. The volume obtained must reach line A. Remove any excess above this line.
- NOTE:**
For all children or patients who have difficulty expectorating in the mixing container, it is recommended to use a collection swab which is available separately. Refer to the respective manufacturer's instructions for use.
- Fill the bottle of reagent #1 vertically, add 1 drop of reagent #1 to the saliva. Hold the opening of the mixing container tightly to avoid spilling of saliva. Tap the mixing container 15 times over a period of 10 seconds with a finger to mix the saliva and reagent #1 thoroughly.

- NOTE:**
Reagent #1 is an alkaline solution which will break down and dissolve viscous components in the saliva. This is an important step in the test process. If the saliva sample is able to easily flow through the test device.
- Ensure saliva bubbles do not collect around the neck of the mixing container as they can obstruct the addition of Reagent #1. If small bubbles in the pack, remove these bubbles with a tissue prior to addition of Reagent #1.
 - Add 4 drops of reagent #2 (yellow) to the mixing container. Shake the device for several seconds to mix. Check that the saliva sample has changed to a light green colour (alkaline to neutral pH).

- Using the graduated pipette, take sufficient saliva from the mixing container to fill to line 3 on the pipette and dispense into the test window at the end of the test device.
- Leave sitting on the bench for 15 minutes (room temperature). A red thin line should be observed in the control (C) window of the test device, indicating that the test is working properly. At the same time, the control (T) window. The result is positive if a thin red line appears in the T window, indicating that salivary levels of *Streptococcus mutans* are high (>5 x10⁵ CFU/mL saliva) and thus the patient has a potential high risk of future caries activity. If no line can be observed after 15 minutes, this indicates a low salivary level of *Streptococcus mutans* and a low potential risk of caries at this time.

- NOTE:**
- Make sure you check the result at 15 minutes after dispensing saliva into the sample window. The results read before or after 15 minutes may be inaccurate.
 - The result is regarded as positive even if the line at the test (T) window is very pale. Refer to picture #14 on the instructions for use.
 - If the red line does not appear in the control (C) window, the test is invalid. Please start the test again using a new test device.

- The test will be positive when the number of *Streptococcus mutans* bacterial cells in 1 ml of stimulated saliva is > 5 x 10⁵ CFU/mL (greater than 500,000 CFU per mL of saliva).
- A negative result indicates that the number of *Streptococcus mutans* is < 5 x 10⁵ CFU/mL (less than 500,000 CFU per mL of saliva).

- MECHANISM**
SALIVA-CHECK MUTANS detects *Streptococcus mutans* in saliva using a highly specific immunochromatography process. When high levels of *Streptococcus mutans* are present in saliva, they react with a colloidal gold-labeled anti-S. mutans monoclonal antibody which is contained in the test device. This gold colloid particles attach to the surface of *Streptococcus mutans* bacteria. The reaction between the antibody and bacteria to form the red line on the T window. The colloidal gold-labeled anti-S mutans monoclonal antibody which did not react with *Streptococcus mutans* reacts with an anti-mouse immunoglobulin in the C window and forms the red line.

- The SALIVA-CHECK MUTANS test kit includes reagents #1 and #2 to overcome difficulties in using saliva as a test sample. Whole human saliva does not easily flow into the test device due to its high viscosity. *Streptococcus mutans* bacteria are often covered with glucans which may inhibit the reaction with the antibody. SALIVA-CHECK MUTANS reagents #1 and #2 are designed to remove these impediments to obtain accurate results. Reagents #1 is NaOH solution, #2 is organic acid solution.

- PERFORMANCE**
- Detectability: To check the intra-lot accuracy, the same positive samples and a detection limit has been evaluated by diluting the pure *Streptococcus mutans* solution (1x 10⁸ CFU/mL). Results show that the number of *Streptococcus mutans* detected is of 9.4x10⁴ CFU in the test device (300 µl of 5x10⁵ CFU/mL solution).
 - Sensitivity-Specificity (Correlation): An evaluation has been conducted on 89 saliva samples. Results of MUTANS test device were compared with those of real-time PCR method.

MUTANS Test Device	Positive	Real-Time PCR	
		Positive	Negative
	10	2	12
	1	76	77
Total	11	78	89

Sensitivity	Specificity	Echtzeit-PCR	
		Positiv	Negativ
90.9%	97.4%	10	2
83.3%	98.7%	76	77
Negative Predictive Value	98.7%	76	77

- C. Reproducibility:**
To check the intra-lot accuracy, the same positive samples and a dilution buffer solution have been processed 15 times on test devices of the same production lot in the same experimental conditions. All observed results were correct as expected. To check the inter-lot accuracy, same samples (positive and dilution buffer) were processed on three different production lots. All results were correct as expected.

- D. Interference:**
Crossed the ability to samples positive for the following oral bacteria was tested and found to be negative: *Streptococcus sobrinus*, *Streptococcus sanguis*, *Streptococcus salivarius*, *Streptococcus mitis*, *Lactobacillus casei*, *Streptococcus gordonii*, *Streptococcus ratti*, *Streptococcus faecalis*, *Streptococcus acidophilus*, *Lactobacillus plantarum*, *Lactobacillus salivarius*, *Porphyromonas gingivalis*, *Prevotella intermedia*, *Tannerella forsythensis*, *Eikenella corrodens*, *Fusobacterium nucleatum*, *Treponema denticola*, *Campylobacter rectus*, *Myxococcus xanthus*, and *Escherichia coli*.

- PACKAGE CONTENTS**
SALIVA-CHECK MUTANS
10 test kits - each kit contains MUTANS Test Device (1), Paraffin Gum (1), Pipette (1), Mixing container (1)
1 bottle of Reagent #1 (2 mL)
1 bottle of Reagent #2 (4 mL)
Optional (sold separately):
Collection Swab (10 pieces)

- STORAGE**
Store the kit in a cool dark place away from high temperatures and direct sunlight. Do not store at or below, or at high temperatures for a long time, as this may cause inaccurate results. Use the test kit immediately after opening the foil bag. (Shelf life: two and a half years from date of manufacture.)

- CAUTION**
- Avoid contact of reagent #1 or #2 with the eyes. In case of contact, flush with copious amounts of water and seek medical attention.
 - Avoid contact of reagent #1 or #2 with the skin. In case of contact, immediately flush with water.
 - A negative test result indicates that the patient has low levels of *Streptococcus mutans* in the mouth and is considered at low risk of developing dental caries at this time. This does not mean that they cannot develop dental caries in the future, since risk factors may change in later years.
 - For accurate diagnosis, the test results should be compared with other available clinical and laboratory information.
 - All components used for each test are intended for single use only. Do not reuse.
 - Dispose of all the components according to local regulations.
 - Reagent #1 is alkaline, reagent #2 is acidic. When disposing of these liquids, mix two together and flush with water.

- Instructions for use: Last updated on October 15, 2007, Rev. # 1.

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