
KIT CONTAINS: SET-INHALT: CONTENUTO DEL KIT: CONTENIDO DEL KIT:

Test device Tester Dispositif d'essai (test) Test Dispositivo per il test Test

- Control Window C Kontrollfenster C Fenêtre de contrôle C Finestra di controllo C Ventana de control C
- Test Window T Testfenster T Fenêtre T Finestra controllo T Ventana de control T
- Sample Window Probenfenster Fenêtre Echantillon Vaschetta campione Ventana de muestra

- Paraffin gum Paraffin Paraffine Gomma di paraffina Gomma de parafina

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
- Fig. 1 Abb. 1

MUTANS Test Device	Real-Time PCR		Total
	Positive	Negative	
	10	2	12
	76	77	77
	11	78	89

Sensitivity 90.9% (10/11)
Specificity 97.4% (76/78)
Positive Predictive Value 83.3% (10/12)
Negative Predictive Value 98.7% (76/77)

Reproducibility : To check the intra lot accuracy, the same positive samples and a dilution buffer solution have been processed 15 times on test devices of the same production lot in the same experimental conditions. All observed results were correct as expected. To check the inter lot accuracy, same samples (positive and dilution buffer) were processed on three different production lots. All results were correct as expected.


D. Interference: Cross-reactivity to samples positive for the following oral pathogens was tested and found to be negative: *Streptococcus sobrinus*, *Streptococcus sanguis*, *Streptococcus salivarius*, *Streptococcus mitis*, *Lactobacillus casei*, *Streptococcus gordonii*, *Streptococcus mutans*, *Actinomyces naeslundii*, *Lactobacillus acidophilus*, *Lactobacillus plantarum*, *Lactobacillus salivarius*, *Porphyromonas gingivalis*, *Prevotella intermedia*, *Tannerella Forsythensis*, *Eikenella corrodens*, *Fusobacterium nucleatum*, *Treponema denticola*, *Campylobacter rectus*, *Myxococcus xanthus* and *Escherichia coli*.

PACKAGE CONTENTS
SALIVA-CHECK MUTANS
 10 test kits - each kit contains MUTANS Test Device (1), Paraffin Gum (1), Pipette (1), Mixing container (1)
 1 bottle of Reagent #1 (2 mL)
 1 bottle of Reagent #2 (4 mL)
 Optional (sold separately):
 Collection Swab (10 pieces)

STORAGE
 Store the kit in a cool dark place away from high temperatures and direct sunlight. Do not store at 0°C or below, or at high temperatures for a long time, as this may cause inaccurate results. Use the test kit immediately after opening the foil bag. (Shelf life: two and a half years from date of manufacture.)

CAUTION
 1. Avoid contact of reagent #1 or #2 with the eyes. In case of contact, flush with copious amounts of water and seek medical attention.
 2. Avoid contact of reagent #1 or #2 with the skin. In case of contact, immediately flush with water.
 3. A negative test result indicates that the patient has low levels of *Streptococcus mutans* in the mouth and is considered at low risk of developing dental caries at this time. This does not mean that they cannot develop dental caries in the future, since risk factors may change in later years.
 4. For accurate diagnosis, the test results should be compared with other available clinical and laboratory information.
 5. All components used for each test are intended for single use only. Do not reuse.
 6. Dispose of all the components according to local regulations.
 7. Reagent #1 is alkaline, reagent #2 is acidic. When disposing of these liquids, mix two together and flush with water.

Instructions for use: Last updated on October 15, 2007, Rev. # 1.

DISTRIBUTED BY
GC CORPORATION
 76-1 Hasunuma-cho, Tabahi-ku, Tokyo 174-8685, Japan

EU: GC EUROPE N.V.
 Research Park Hasselt-Luven 1240, Interleuvenlaan 33,
 B-3001 Leuven, Belgium Tel: +32 16 39 80 50

GC AMERICA INC.
 3737 West 127th Street, Alhambra, IL 60802 U.S.A.
 Tel: +1 708-597-0900

GC ASIA DENTAL PTE. LTD.
 110, Linyang Way, #09-27 Singapore 508724
 Tel: +65 6546 7588

MANUFACTURED BY
**SHENZHEN KANG SHENG
 BAO BIO-TECHNOLOGY CO., LTD.**
 Block 2, Bao An Industrial Park Phase 1, Da Yuan Development Zone, Sub Zone II,
 Fu Yong, Bao An District, Shenzhen, People's Republic of China, 518103

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MECHANISME
 SALIVA-CHECK MUTANS detects *Streptococcus mutans* in saliva using a highly specific immunochromatography process. When high levels of *Streptococcus mutans* are present in saliva, they react with a colloidal gold-labeled anti-S. mutans monoclonal antibody which is contained in the test device. Thus, gold colloid particles attach to the surface of *Streptococcus mutans* bacteria. This reaction causes the red line to form in the control window (C) of the test device. The colloidal gold-labeled anti-S. mutans monoclonal antibody which did not react with *Streptococcus mutans* reacts with an anti-mouse immunoglobulin in the C window and forms the control red line.

Die SALIVA-CHECK MUTANS test kit includes reagents #1 and #2 to overcome difficulties in using saliva as a test sample. Whole human saliva does not easily flow into the test device due to its high viscosity. *Streptococcus mutans* bacteria are often covered with glycans which may inhibit the reaction with the antibody. SALIVA-CHECK MUTANS reagents #1 and #2 are designed to remove these impediments to obtain accurate results. Reagents #1 is NaOH solution, #2 is organic acid solution.

PERFORMANCE
A. Detectability:
 The detection limit has been evaluated by diluting the pure *Streptococcus mutans* sswell (1×10^4 CFU/mL). Results show that the number of *Streptococcus mutans* detected is of 9.4×10^4 CFU in the test device (300 µl of 5×10^5 CFU/mL solution).
B. Sensitivity-Specificity (Correlation):
 An evaluation has been conducted on 89 saliva samples. Results of MUTANS test device were compared with those of real-time PCR method.

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Sensitivity 90.9% (10/11)
Specificity 97.4% (76/78)
Value predictive positive 83.3% (10/12)
Value predictive negative 98.7% (76/77)

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