# **G-Cera ORBIT VEST**

Revised October 2018

S3305  
950166

**EN Refractory die material**

**DE Feuerfestes Stumpfmateriale**

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**SV Eldfast stansmateriale**

**DA Ildfast modelmateriale**

The logo for G-Cera, featuring the letters 'G', 'C', and 'E' in a bold, italicized, sans-serif font. The 'G' and 'E' are larger and more prominent, with the 'C' nestled between them. The letters are black with a slight shadow effect.

EN

Prior to use, carefully read  
the instructions.

## **INSTRUCTIONS FOR USE**

G-CERA ORBIT VEST is a phosphate-bonded refractory die material developed for production of all-ceramic crowns, inlays, onlays and veneers. G-Cera Orbit Vest allows various porcelains to be built directly onto the refractory model; has outstanding dimensional stability and provides sufficiently high strength to withstand repeated firings.

### **FEATURES**

#### **1. TWO TYPES OF POWDER OFFERED**

G-Cera Orbit Vest has two types of powder with different expansions for crowns/veneers, and for inlays. This provides for a wide range of use in producing both external ceramic restorations such as crowns or veneers, and internal ceramic restorations such as inlays.

#### **2. SUPERIOR STRENGTH**

G-Cera Orbit Vest has sufficiently high strength to endure repeated firings. It has been formulated to curb the expansion rate to a minimum in a given range of temperature rise, causing no crack or damage to its model.

#### **3. SMOOTH SURFACE**

Thanks to its unique granular formula, G-Cera Orbit Vest provides a smooth model surface.

#### **4. EXCELLENT COMPATIBILITY WITH VARIOUS PORCELAINS**

G-Cera Orbit Vest is designed for use with most metal bonding porcelains which have average coefficients of thermal expansion of  $11.5-13.5 \times 10^{-6}$  at 100-500°C (212-932°F). It will not cause the porcelain to crack or separate from the model during firing.

### **PHYSICAL DATA**

<b>Item</b>	<b>Value</b>
Powder / Liquid ratio	33 g / 5 ml (6.5g)
Pouring time (min. sec.) (after completion of mix)	2'00"
Setting time (min. sec.)	7'00"
Compressive strength (MPa) (after degassing)	30.6

Test conditions: Temperature (23±2°C) (73.4±3.6°F) - Relative humidity (50±5%).

## **DIRECTIONS FOR USE:**

### **1. THE POWDER / LIQUID RATIO**

The powder/liquid ratio is 33g (one pack of powder) / 5ml (6.5g).

*Note:*

*Powder and liquid should be correctly dispensed. Otherwise proper physical properties such as thermal and setting expansions can not be obtained.*

### **2. MIXING**

After hand-mixing thoroughly, vacuum mix uniformly for 60 seconds.

*Note:*

*The smallest possible mixing container (5-6cm internal diameter) should be used. The use of a larger container will not allow for a uniform consistency, resulting in improper physical properties such as thermal and setting expansion.*

### **3. POURING THE IMPRESSION**

Prior to pouring the refractory, a wetting agent can be applied to the impression in order to obtain a bubble-free refractory model. After applying the wetting agent, lightly rinse under running water and dry well. After mixing the refractory, immediately pour the mixture into the impression while using a gypsum vibrator (maximum pouring time after mix: about 2 minutes).

*Note:*

*Some wetting agents may roughen the surface of the refractory model.*

### **4. REMOVING FROM THE IMPRESSION**

Remove the model from the impression 1-1.5 hour after pouring.

*Note:*

*Avoid early removal as it will adversely effect the strength, dimensional precision and surface conditions of the model. Also do not allow the model to stay in the impression for over 2 hours, or you will have difficulty with separation.*

### **5. DEGASSING THE MODEL**

Gradually heat the model in the ring furnace from room temperature up to 700°C (1,292°F) over 1 hour. Maintain (700°C / 1,292°F) temperature for 10 minutes.

*Note:*

*Too rapid a rise in temperature may fracture the model due to too rapid evaporation of moisture inside the refractory model.*

Next, place the model in the preheated (700°C / 1,292°F) porcelain furnace and raise to 1,050°C (1,922°F) at a rate of 55°C (131°F) per minute under atmospheric pressure. Maintain this temperature for 5 minutes. Remove the model and bench-cool to room temperature. It is not necessary to bench-cool in the porcelain furnace.

*Notes:*

- 1) Do not degass the model at a temperature higher than 1,050°C (1,922°F). Otherwise it may affect its dimensional precision.
- 2) Degassing the model can also be performed in a porcelain furnace, but phosphate-bonded refractory materials give off small amounts of ammonia gas, which may affect the porcelain furnace. Therefore a ring furnace is recommended to heat the model up to 700°C (1,292°F).

#### **SAFETY RECOMMENDATIONS AND HAZARD WARNINGS:**

1. Before use read the powder and liquid products Safety Data Sheets communicated by your dealer.  
These are also available at <http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>
2. The solid product contains quartz and cristobalite. Causes damage to the lung through prolonged or repeated exposure.  
Route of exposure: Inhalative. Do not breathe dust. In case of inadequate ventilation wear respiratory protection.
3. Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.  
Rinse the empty investment material bag with water before disposal.
4. Remove dust from your working place only when it is wet.
5. To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
6. When sandblasting the cast object, always use a fine-dust filter extraction system.
7. During heating up the material, use of a fume hood is mandatory. Working in closed places is to be avoided.  
Do not inhale fumes when heating the material.
8. Never touch hot materials by hand, always use appropriate tools like muffle grippers in order to avoid burns.
9. Take care of the heat caused by the chemical setting of the material, in order to avoid harm.
10. This product is to be used only by a dental professional.
11. This product should be used only for the applications described in the instructions.
12. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

#### **NOTES**

1. If the liquid should be contaminated with extraneous matter such as plaster, acid, table salt, etc., or should be placed below 5°C (41°F), it will turn into a gel and can not be used again.
2. Working temperature: Keep both powder and liquid away from direct sunlight. Recommended for optimal performance, store at room temperatures (20-25°C / 68-77°F).
3. After opening the powder pack, use material up as soon as possible.
4. Some components contained in the powder will delay the setting time of gypsums or gypsum bonded investments. Designate a mixing bowl to be used exclusively with G-Cera Orbit Vest.

#### **PACKAGES:**

1. Box of 30 x 33g powder for crowns/veneers.
2. Box of 30 x 33g powder for inlays.
3. Bottle of 200ml liquid.

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
TAISEI DENTAL MFG. CO., LTD.