COE-FLEX® LEAD-FREE

(P)

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- Misture e preencha a moldeira com COE-FLEX LEAD-FREE Regular ou COE-FLEX LEAD-FREE Heavy.

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(GB)

Polysulfide Impression Material

- Injection - Type 3 Light-bodied consistency
- Regular - Type 2 Medium-bodied consistency
- Regular - Type 2 Fast Set, Medium-bodied consistency
- Heavy - Type 1 Heavy-bodied consistency

Prior to use, carefully read the instructions for use.

Do not intermix the Catalyst or Base from the original COE-FLEX Impression Material with the COE-FLEX LEAD-FREE Catalyst or Base.

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS:

For taking impressions for full dentures, partial dentures, inlays, onlays, and crowns and bridges.

CONTRAINDICATIONS:

In rare cases this material may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

STORAGE:

Store COE-FLEX LEAD-FREE Impression Material in a dry location at room temperature (70° to 77°F; 21° to 25°C). Material stored in such conditions will have a shelf life of 36 months.

CAUTION:

- If unmixed material comes in contact with oral tissue or skin, remove and flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingestion of the material.

TO AVOID STAINED CLOTHING:

Let the mixed COE-FLEX LEAD-FREE Impression Material cure completely before attempting to remove it from skin or clothing; it may be peeled off in some instances. NEVER attempt to wipe or rub a spot freshly stained with mixed COE-FLEX LEAD-FREE Impression Material.

If some of the catalyst or base paste (while yet unmixed) stains a garment, it can be removed in some instances with standard dry cleaning solvent. Apply solvent liberally with a soft cloth while successive pieces of clean, white absorbent paper or blotter are held under the spotted area. After most of the stain is removed, wash out the remaining spot with soap and water.

DIRECTIONS FOR USE:

1. PREPARATION:

Use a custom tray constructed of COE TRAY PLASTIC™ or stock metal or plastic tray. Trays must be absolutely rigid. Paint COE-FLEX Adhesive evenly over entire surface of tray. When adhesive is dry, mix and apply the COE-FLEX LEAD-FREE Impression Material.

TIMING:

At normal room temperature (74°F / 23°C ±4°F / 2°C and 50±10% relative humidity). Elevated temperature and/or humidity will accelerate the setting reaction.

Physical and Mechanical Properties:

Consistency / Type	Light-Bodied Type 3	Medium-Bodied Type 2	Medium-Bodied Fast Set – Type 2	Heavy-Bodied Type 1
Mixing Time	45 sec – 1 min	45 sec – 1 min	30 sec	45 sec – 1 min
Total Working Time	4 min	3 min 30 sec	2 min 30 sec	3 min
Mouth Removal Time Post Insertion (Minimum Time in Mouth)	8 – 10 Minutes	8 – 10 Minutes	4 minutes	8 – 10 Minutes
Max Strain in Compression	20%	20%	20%	20%
Recovery from Deformation	min 96.5	min 96.5	min 96.5	min 96.5
Linear Dimensional Change	max 0.5%	max 0.5%	max 0.5%	max 0.5%

Test Conditions: Temperature (23±2)°C, Relative humidity (50±10)%

To Speed up Setting Time - Use more catalyst (blue paste) in proportion to the amount of white paste.

To Retard Setting Time - Use less catalyst (blue paste) in proportion to the amount of white paste.

PROPORTIONS:

When using a new tube: Remove the cap and then completely remove the foil tamper evident seal. It is important to completely remove the seal so that the flow of material is not restricted. Failure to do so may result in incorrect mixing ratios and poor or non performance of the product.

By Weight: 5 parts Catalyst with 9 parts Base.

By Volume: 1 part Catalyst with 2 parts Base.

(Catalyst and Base tube orifice diameters has been designed accordingly.)

Extrude equal lengths of blue and white pastes on a mixing pad. Recap each tube after use. The diameter of the ropes of paste on the mixing pad must be the same diameter as the orifice of the tubes. Thinner and thicker ropes give improper proportions.

NOTICE: The conversion to plastic laminate tubes resulted in a change in the orifice size of the base and catalyst vs. the orifice sizes of the original formula COE-FLEX base and catalyst that was packaged in metal tubes. You will now need to extrude longer lengths (approximately 1.5 times longer) of base & catalyst to achieve the same volume as the old packing provided.

2. MIXING:

Coat the spatula with the blue paste before incorporating the blue and white pastes together. Spatulate the two lengths of paste together with broad sweeps. A proper mixture is free from white or blue streaks. High energy mixing accelerates the setting reaction.

- Fill an impression syringe (COE Syringe or GC Plastic Syringe) with a sufficient amount of COE-FLEX LEAD-FREE Injection material to inject onto and around prepared teeth.
- Mix and fill the tray with COE-FLEX LEAD-FREE Regular or COE-FLEX LEAD-FREE Heavy.
- When using the syringe, inject prepared teeth with a steady, even pressure from the posterior to anterior region. Inject over the occlusal surface of a non-involved tooth to complete injection procedure. Keep tip of syringe as close to preparations as possible.

3. IMPRESSION TAKING:

Firmly place the filled tray over the impression area. Do not apply pressure against the tray while COE-FLEX LEAD-FREE is setting. This may cause movement and distortion of the impression. Do not remove impression prematurely. Touch excess material extruded from the tray with a blunt instrument. COE-FLEX LEAD-FREE is cured when no indentation remains.

4. REMOVAL:

Break seal with a slight lateral movement and remove impression with a straight pull.

Rinse off adhering saliva under running lukewarm water.

5. POURING OF MODELS AND DIES:

Pour casts as soon as possible after the impression has been removed from the mouth; no later than 8 hours after the impression has been taken. Recommended gypsum products: COECAL Dental Stone or FujiRock EP Die Stone.

DISINFECTION AND STERILIZATION OF THE IMPRESSION:

COE-FLEX LEAD-FREE impression can be disinfected by immersion in a glutaraldehyde solution (2.5 to 3.4% concentration). Refer to the manufacturer's instructions for recommended immersion times and temperatures.

PACKAGE:

Standard Package*: Contains 1 Base (90 mL), 1 Catalyst (45 mL), 1 Adhesive

Refill Package*: Contains 1 Base (90 mL), 1 Catalyst (45 mL)

COE-FLEX LEAD-FREE Catalyst* Only (45 mL)

COE-FLEX LEAD-FREE Base* Only (90 mL)

*Consistencies: Injection - Type 3 Light-Bodied
Regular - Type 2 Medium-Bodied
Regular Fast Set - Type 2 Medium-Bodied
Heavy - Type 1 Heavy-Bodied

GC®

MANUFACTURED by
GC AMERICA INC.
3737 West 127th Street
Alsip, Illinois 60803 U.S.A.
Tel: +1-708-597-0900

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