Prior to use, carefully read the instructions for use.

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EXAJET™ 370

HYDROPHILIC VINYL POLYSILOXANE IMPRESSION MATERIAL

For use only by a dentist or physician.

INDICATIONS FOR USE

Tray material for both double mix and putty wash impression techniques. To be used in a dynamic mixing device.

CONTRAINDICATIONS

In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

PROPERTIES

	EXAJET 370 Heavy Body Type 1: High consistency Normal set	EXAJET 370 Heavy Body Type 1: High consistency Fast set	EXAMIX Regular Type 2: Medium consistency	EXAMIX Injection Type 3: Low consistency
Total working time (min.sec.)	2'00"	1'30"	2'00"	2'15"
Setting time (min. sec.)	4'00"	2'15"	4'00"	4'00"
Minimum time in mouth (min. sec.)	3'00"	2'30"	4'00"	4'00"
Recovery from deformation (%)	99.4	99.4	99.7	99.7
Maximum strain in compression (%)	2.5	2.5	6.0	7.0
Linear dimensional change after 24 hrs. (%)	0.2	0.2	0.2	0.2

DIRECTIONS FOR USE CARTRIDGE LOADING AND DISPENSING

NOTE: The cartridges are heavy and dropping may result in damage to the cartridges rendering it unusable. The cartridge and mixing tip provided are compatible with a dynamic mixer. Also refer to the manufacturers instructions.

- 1) Hold and support the cartridge with the extrusion nozzles facing upwards.
- 2) Grasp the tab on the cap and bend in an upward 90° direction. Do not place excessive stress on tab.
- 3) While holding the cartridge securely, grasp the tab with thumb and forefinger and remove the cap from the cartridge.
- 4) Insert the cartridge into the mixer according to the manufacturers instructions.
- 5) Start the mixer according to the mixer manufacturer's instructions. Start mixing the material, extrude and discard the initial unevenly mixed paste. Wipe off the extruded material in a vertical direction to avoid cross contamination of the nozzles.
- 6) If the mixing tip does not seat easily, be sure the internal central hex of the mixing tip is properly aligned with the drive shaft of the mixer.
NOTE: Do not use other manufacturer's mixing tips. They will not fit properly and cause the material to leak from the nozzles.
- 7) Once the mixing tip is properly placed, place the nozzle locking ring over the mixing tip and turn it 1/4 turn clockwise until ring is locked in place.
- 8) Continue to extrude the material until the mix is homogenous (extrude 3-6 cm). Do not use if streaks (improper mixing) appear.
- 9) Repeat steps 5, 6, 7 and 8 each time the machine is in use for impression taking.

IMPRESSION MAKING PUTTY WASH TECHNIQUE

- 1) Coat the impression tray (e.g. Coe trays) with a thin layer of adhesive and allow to dry.
- 2) Load EXAJET 370 material into the tray.
- 3) Take the pre-impression.
- 4) Remove after the material has set, clean thoroughly and dry.
- 5) Prepare the pre-impression for final impression.
- 6) Inject EXAMIX Injection or Regular onto the prepared teeth.
- 7) Seat tray with pre-impression into the mouth.
- 8) Once the setting time has elapsed, remove the impression from the mouth.
- 9) The resulting impression should be cleaned and then disinfected according to the manufacturer's recommendations.

DOUBLE MIX TECHNIQUE

- 1) Coat the impression tray (e.g. Coe trays) with a thin layer of adhesive and allow to dry.
- 2) Inject EXAMIX Injection or Regular onto the prepared teeth.
- 3) Load EXAJET 370 material into the tray.
- 4) Seat tray into the mouth.
- 5) Once the setting time has elapsed, remove the impression from the mouth.
- 6) The resulting impression should be cleaned and then disinfected according to the manufacturer's recommendations.

Note: Make sure to complete step 2 when the tray is fully loaded.

NOTES

Wait one hour before pouring the model.
EXAJET 370 impressions can be silver- or copper-plated.

CAUTION

- 1) When extruding EXAJET 370, care should be taken to avoid mixing or contact with the following materials. They may delay the setting time.
 - Condensation cured silicone impression materials
 - Polysulfide impression materials
 - Polyether impression materials
 - Eugenol containing materials
 - Certain types of latex gloves or retraction cords
 - Also avoid moisture and glycerol when extruding
- 2) Care should be taken to avoid getting material on clothing. It is hard to remove when it sets up on clothes.
- 3) Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

STORAGE

Recommended for optimal performance, store at normal temperature 15-25°C (59-77°F).
Shelf-life:
2 years from date of manufacture.

PACKAGES

- 1) Intro package Heavy Body – Normal or Fast set
 - 1 cartridge
 - 10 dynamic mixing tips
- 2) Refill package Heavy Body – Normal or Fast set
 - 2 cartridge
- 3) Clinic package Heavy Body – Normal or Fast set
 - 6 cartridge
- 4) 50 dynamic mixing tips

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available for the Americas at:
<http://www.gcamerica.com/downloads>

Last Revised: 05/2017

Rx Only

'GC'

MANUFACTURED for
GC AMERICA INC.
3757 West 127th Street
Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

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