EN
Prior to use, carefully read the instructions for use.

G-aenial™ Bond

ONE COMPONENT SELF-ETCHING LIGHT-CURED ADHESIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Bonding of light cured composites and acid modified composites (compomers) to tooth structure.
2. Bonding of dual cured luting and core build up composites to tooth structure as long as these materials are light cured.

CONTRAINDICATIONS

1. Pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED

1. In combination with chemically cured composite resin.
2. In combination with eugenol containing materials as eugenol may hinder G-aenial Bond from setting or bonding properly.
3. In combination with desensitizers as desensitizers may hinder G-aenial Bond from setting or bonding properly.
4. In combination with dual cured luting and core build up composites unless both G-aenial Bond and these materials are individually light cured.

DIRECTIONS FOR USE

1. **CAVITY PREPARATION**
Prepare cavity using standard techniques. Make sure that the temporary material is thoroughly removed from the surface. Use rubber dam for protection. Dry the prepared tooth surfaces by gently blowing with an air syringe.
Note:
For pulp capping, use calcium hydroxide.

2. TECHNIQUE SELECTION

1. Select out of the following 2 technique options
 - a) Self-etch technique – apply G-aenial Bond to bur cut enamel and dentin, without a separate etching step.
 - b) Selective etching of enamel – prior to applying G-aenial Bond to enamel and dentin, etch (bur cut) enamel with 35-40% phosphoric acid gel for 10 seconds, rinse for 5 seconds and gently dry.
Note:
Non bur cut enamel should always be treated with a 35-40% phosphoric acid gel for 10 seconds, rinsed with water for 5 seconds and gently dried.

3. APPLICATION

- a) Prior to dispensing, shake the bottle of G-aenial Bond thoroughly (Fig. 1). Dispense a few drops into a clean dispensing dish (Fig. 2). Replace bottle cap immediately after use (Fig. 3).
- b) IMMEDIATELY apply to the prepared enamel and dentin surfaces using the dispensable applicator (Fig. 4).
- c) Rinse undisturbed for 10 seconds after the end of application (Fig. 5).
- d) Then, dry thoroughly for 5 seconds with oil free air under MAXIMUM air pressure. Use vacuum suction to prevent splatter of the adhesive (Fig. 6). The final result should be a thin adhesive film with the appearance of frosted glass and which does not visibly move under further air pressure.
Note:
When removing from refrigerator after long-term storage, leave G-aenial Bond to stand at room temperature for several minutes prior to use.

- 2) Apply G-aenial Bond immediately as the material contains a volatile solvent.
- 3) Remove excess material remaining on the tooth (other than the surfaces to be bonded) using a sponge or cotton pellet as residue material is hard to remove after light curing.
- 4) Should the applied material be contaminated with water, blood or saliva prior to light curing, wash and dry the tooth and repeat the procedure by re-applying material.

4. LIGHT CURING

- Light cure using a visible light curing unit (Fig. 7).
Irradiation time
Halogen / LED (700 mW/cm²) : 10 seconds
High power LED (more than 1200 mW/cm²) : 5 seconds

In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following time:
Halogen / LED (700 mW/cm²) : 20 seconds
High power LED (more than 1200 mW/cm²) : 10 seconds

- Note:
1) Light cure completely for an effective bonding strength. Lower light intensity may cause insufficient adhesion.
2) Use a protective light shield or similar protective eye wear during light curing.

- 5A. **PLACEMENT OF LIGHT CURED COMPOSITES AND COMPOMERS**
After having light cured the adhesive, follow manufacturer's instructions for use in order to place, contour and light cure the composite resin material of choice.
- 5B. **PLACEMENT OF DUAL CURED COMPOSITES**
After having light cured the adhesive, make sure to light cure the dual cured composite material separately. Self curing only will result in insufficient bonding strength.
6. **FINISHING**
Adjust restoration, finish and polish using standard techniques.
7. **STORAGE**
Store at room temperature (1-28°C) (33.8-82.4°F). If not in use for a prolonged period of time, store in refrigerator.

1. **PACKAGE**
G-aenial Bond Bottle Kit: 5 mL liquid (1), DISPOSABLE DISPENSING DISH (20), DISPOSABLE APPLICATOR II (50)
2. G-aenial Bond Bottle Refill: 5 mL liquid (1)
3. G-aenial Bond 3-bottle Pack: 5 mL liquid (3)
1. **CAUTION**
G-aenial Bond is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. G-aenial Bond is volatile. Use in a well ventilated place. Advise patients to breath through nose.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
5. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until a mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply cocoa butter to the area where rubber dam cannot cover.
6. Avoid inhalation or digestion of material.
7. When spilt on the table or floor, immediately wipe off with a dry cloth.
8. Do not mix with other products.
9. Dispose of all waste according to local regulations.
10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
11. In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

- 5A. **PLACEMENT OF LIGHT CURED COMPOSITES AND COMPOMERS**
After having light cured the adhesive, follow manufacturer's instructions for use in order to place, contour and light cure the composite resin material of choice.
- 5B. **PLACEMENT OF DUAL CURED COMPOSITES**
After having light cured the adhesive, make sure to light cure the dual cured composite material separately. Self curing only will result in insufficient bonding strength.
6. **FINISHING**
Adjust restoration, finish and polish using standard techniques.
7. **STORAGE**
Store at room temperature (1-28°C) (33.8-82.4°F). If not in use for a prolonged period of time, store in refrigerator.

1. **PACKAGE**
G-aenial Bond Bottle Kit: 5 mL liquid (1), DISPOSABLE DISPENSING DISH (20), DISPOSABLE APPLICATOR II (50)
2. G-aenial Bond Bottle Refill: 5 mL liquid (1)
3. G-aenial Bond 3-bottle Pack: 5 mL liquid (3)
1. **CAUTION**
G-aenial Bond is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. G-aenial Bond is volatile. Use in a well ventilated place. Advise patients to breath through nose.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
5. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until a mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply cocoa butter to the area where rubber dam cannot cover.
6. Avoid inhalation or digestion of material.
7. When spilt on the table or floor, immediately wipe off with a dry cloth.
8. Do not mix with other products.
9. Dispose of all waste according to local regulations.
10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
11. In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

- Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com
or for The Americas
http://www.gcamerica.com
They can also be obtained from your supplier.

- CLEANING AND DISINFECTING:**
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERGE. Thoroughly clean device to prevent dry coating and label with the following link:
http://www.gcamerica.com
healthcare-grade infection control product according to regional/national guidelines.

Last revised: 10/2018

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Tomimatsu-cho, Kasugai, Aichi 466-0844, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
EU : GC EUROPE N.V.
Researchpark Haarrodde-Leuven 1240, Internevenlaan 33,
8-3001 Leuven, Belgium TEL : +32 16 74 10 90
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL : +1-708-597-0900
www.gcamerica.com
GC SOUTH AMERICA
Rua Heliópolis 399, Sítio Heliópolis, São Paulo, SP, BRASIL
CEP: 02022-051 - TEL.: +55-11-2925-0985
CNPJ: 08.391999/0001-61
RESP. TEC. Mayara de Santis Ribeiro - CROSP/15.009.2
GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05 Singapore 529729.
TEL.: +65 6464 7698
GC AUSTRALASIA DENTAL PTY. LTD.
1723 Botany Rd, Bankstown, NSW 2019, Australia
TEL.: +61 2 9301 8200
PRINTED IN JAPAN

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