Prior to use, carefully read the instructions for use.

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ES

# GC Fuji LINING LC

RADIOPAQUE LIGHT CURED GLASS IONOMER LINING CEMENT

For use only by a dental professional in the recommended indications.

## RECOMMENDED INDICATIONS

As a base or liner in prepared cavities.

## CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of the product and refer to a physician.

## DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	1.4 / 1.0
Mixing Time (sec.)	20"-25"
Working Time (min., sec.)	3'00"
Light Curing Time (sec.)	30"
Depth of Cure (mm)	1.6

Test conditions : Temperature (23 +/-1°C)  
Relative humidity (50 +/-10%)

## 1. POWDER AND LIQUID DISPENSING

- a) The standard powder to liquid ratio is 1.4/1.0g. 1 level scoop of powder to 1 drop of liquid.
- b) For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- c) Hold the liquid bottle vertically and squeeze gently.
- d) Close bottles immediately after use.

## 2. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into two equal parts. Mix the first portion with all of liquid for 5-10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 10-15 seconds.

## 3. RESTORATIVE TECHNIQUE

- a) Prepare tooth using standard techniques. Extensive mechanical retention is not necessary.
- b) Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer. For pulp capping, use calcium hydroxide.
- c) Mix the required amount of cement. Working time is 3 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- d) Remove surface moisture, but DO NOT DESICCATE.
- e) Transfer cement to the preparation using a syringe or suitable placement instrument, covering dentine up to the dentino-enamel junction.
- f) Light cure with a suitable visible light curing device (470nm wavelength) for 30 seconds. Place light source as closely as possible to the cement surface.  
Note :  
If application exceeds 2 mm in depth, use a layering technique.
- g) Proceed with standard composite or amalgam restorative techniques.

## STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life : 3 years from date of manufacture.)

## PACKAGES

1-1 pkg. : 5g powder, 4g (3.4mL) liquid, powder scoop.

## CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. DO NOT mix powder or liquid with any other glass ionomer components.
4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at:

<http://www.gcamerica.com>

They can also be obtained from your supplier.

## CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised : 02/2016

Rx Only

MANUFACTURED by  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan  
RESPONSIBLE MANUFACTURER IN CANADA:  
GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED by  
GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900

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