Prior to use, carefully read the instructions for use.

(EN)

Miracle Mix®

RADIOPAQUE SILVER ALLOY-GLASS IONOMER CEMENT RESTORATIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Core build-up
- Restoration of primary teeth
- Cases where radiopacity is required
- As a base or liner
- Class I, limited Class II, temporary fillings, etc.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced,discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

- POWDER AND LIQUID DISPENSING

Pour the Miracle Mix alloy into the bottle of Miracle Mix powder and shake thoroughly.

 - The recommended powder to liquid ratio is 2 or 3 level scoops of powder to 2 drops of liquid.
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.

Shake the powder bottle occasionally.
- MIXING

Dispense powder and liquid onto pad or glass slab. Divide the powder into 2 equal parts. Mix the first portion with all of liquid for 15-20 seconds. Incorporate the remaining portion and mix the whole thoroughly for 20 seconds.

Note :

A heavy consistency is recommended.
- RESTORATIVE TECHNIQUE
 - Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
 - Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer.
 - Mix the required amount of cement. Working time is 1 min. 30 seconds at 23°C (73.4°F). Higher temperatures will shorten working time.
 - Remove surface moisture but DO NOT DESICCATE.
 - Transfer cement to the preparation using a syringe or suitable placement instrument. Avoid air bubbles.
 - Form the contour.
 - Setting time is 5'00" to 5'30" from start of mixing.
- FINISHING

Use standard techniques.

STORAGE

Store at temperature of 4-25°C (39,2-77,0°F). (Shelf life : powder 3 years, liquid 2 years from date of manufacture)

PACKAGES

- Package : 15g powder with measure, 10g (8,0mL) liquid, 17g alloy.
- Bottle of 15g powder with measure.
- Bottle of 10g (8,0mL) liquid.
- Bottle of 17g alloy.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.

Last revised : 12/2008

CE 0086

GC

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6546 7588

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