EXABITE™ II NDS

Prior to use, carefully read the instructions for use.

VINYL POLYSILOXANE BITE REGISTRATION CREME

EXABITE II NDS is an improved version of EXABITE, an addition reaction polysiloxane impression material, providing sharper setting and higher strength. It is available in a new cartridge which must be used with the new GC CARTRIDGE DISPENSER 2. For use only by a dental professional in the recommended indications.

INDICATIONS

Checking occlusal registrations.


CONTRAINDICATIONS

Patients with known allergy to silicone materials: In case of allergy refer to a physician.

PHYSICAL PROPERTIES

COMPARISSON DES PROPRIÉTÉS	EXABITE I	EXABITE II
Working Time (sec)*	25**	20**
Minimum retention from mouth (sec)	45**	60**
Maximum strain in compression (%)	0.9	1.0
Recovery from deformation (%)	99.6	99.9
Linear dimensional change (%)	0.3-1.0	0.3-1.0
Hardness (Shore A)	67	69

EXABITE II NDS - Working Time (at 74°/23°C) and Minimum Time to Remove From Mouth:

A = Working time B = Holding time in the mouth

DIRECTIONS FOR USE - New Cartridge Dispensing

- Lift the release lever of the CARTRIDGE DISPENSER 2 (referred to as the dispenser hereafter) and pull the piston plunger all the way back into the dispenser. Lift the cartridge holder of the dispenser and load the cartridge, ensuring that the V-shaped notch on the flange of the cartridge is facing down. Push the cartridge holder down to hold the cartridge firmly in place.
- Remove the cartridge cap by rotating 1/4 turn anti-clockwise. Tilt the cap downward and peel it away from the cartridge. Gently squeeze the dispenser handle to extrude a small amount of material from the two openings at the end of cartridge. Make sure that base and catalyst come out evenly.
- Align the V-shaped notch on the rim of MIXING TIP to the V-shaped notch between the cartridge barrels. Push firmly to attach the mixing tip. Then, rotate the colored collar of the mixing tip 1/4 turn clockwise to the end of the cartridge.
- Squeeze the handle several times to extrude the material. After use, do not remove the MIXING TIP as this will become the storage cap until next use. When replacing the MIXING TIP, rotate the collar on the mixing tip 1/4 turn anti-clockwise to align the V-shaped notch on the cartridge. Tilt the mixing tip downward and peel it away from the cartridge.
- Remove and replace the old mixing tip, immediately prior to next use. Before attaching a new tip, gently extrude a small amount of material to ensure that base and catalyst are flowing evenly from both openings. If the materials should fail to extrude, remove any hardened materials from the end of the cartridge.
- To replace the cartridge, lift the release lever and retract the piston plunger fully. Remove the empty cartridge by lifting the cartridge holder, and load a new one into the dispenser.

CLINICAL CASE

Occlusal Registration for Protheses.

- Instruct the patient to occlude into the centric position. Inject material into a space between the opposing tooth and the abutment tooth from the buccal side.
- When taking an occlusal registration, inject material onto the occlusal surface, then let the patient occlude.

NOTES

- To reduce the potential for improper model articulation, inject material onto the occlusal surfaces on both sides of the arch.
- The working time may be as little as 20 seconds when the material is injected directly into the mouth (at 95°/35°C).
- Remove the registration from the mouth. Trim undercuts and excess material using a sharp knife. After trimming, make sure that the registration fits the prepared working models and articulate the case as necessary.

NOTE: Use only a sharp knife to trim the material. The use of scissors might fracture the bite registration.

Registration of Checkbite


Using EXABITE II NDS, it is possible to record bite registration in eccentric positions with the jaw moved forward, right, left, and to the center respectively.

Caution

- In case of contact with eyes, flush immediately with water and seek medical attention.
- When extruding EXABITE II NDS, care should be taken to avoid mixing or contact with the following materials. They may delay setting time.
 - Condensation type silicone impression material
 - Polyoxide impression material
 - Eugenol type materials (Ex. GC Impression Paste)
 - Others (Water, Glycerine)
- Sulfur -Oil -Latex -Acrylates
- The obtained bite registration should be cleaned then disinfected, utilizing a 2.5% or 3.4% glutaraldehyde, or other appropriate disinfectant, according to the manufacturer's label recommendations.
- Avoid getting material on clothing. It is difficult to remove once it has set.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be used.

IDENTIFICATION OF PARTS

V-shaped Notch between the cartridge barrels

Colored Collar for Mixing Tip 2L

STORAGE
Recommended for optimal performance store at room temperature (74°/23°C, 50± 10%, Relative humidity).

PACKAGES

Cartridge 83g (48mL) x 2
Mixing Tip x 12
Cartridge 83g (48mL) x 8
Mixing Tip x 48
(Minimum Mixed Volume 96mL)

Option:
1. CARTRIDGE DISPENSER II 1 piece
2. MIXING TIP 60 pieces

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at: <http://www.gccorporation.com>

or for the Americas: <http://www.gcamerica.com>

They can also be obtained from your supplier.

UNDESIRED EFFECTS - REPORTING
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not included in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.
Last revised: 03/2020

MANUFACTURED BY GC CORPORATION

76-1 Hasegawa-cho, Saitama-shi, Tokyo 174-8585, Japan
Tel: +81-3-3965-1221

RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC.
3737 West 127th Street, Atsip LK, 60803 U.S.A.

DISTRIBUTED BY GC AMERICA INC.
3737 West 127th Street, Atsip LK, 60803 U.S.A.
Tel: +1-708-991-0900
www.gcamerica.com

EU: GC EUROPE N.V.
Researchpark Hasselt-Luven 1240, Interleuvenlaan 33, B-3001 Luven, Belgium
Tel: +32-30-212121

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01, Tampines Plaza, Singapore 529541
Tel: +65-6546-7888

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia TEL: +61 2 9301 8200

GC SOUTH AMERICA
Rua Itaipava, 991, Santana - São Paulo, SP - Brasil
RESP: (0020) 9051 - TEL: +55-11-2925-0965 - CNPJ: 08.279.999/0001-61
RESP: TEL: Mayara De Santos Ribeiro - CNP/SP: 105.982

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