


G-ænial™ A'CHORD

Universal Light-Cured Radiopaque Composite Restorative

Simplicity, esthetics and performance in your hands



Since 1921

100 years of Quality in Dental




Achieve the Esthetics of 16 Classic VITA®+ Shades with Just 5 Core Shades

GC's proprietary Full-Coverage Silane Coating (FSC) and High-Performance Pulverized CERASMART® (HPC) technology provides long-lasting, natural looking restorations with simplified handling and delivery.

Superb Handling


- Sculpts easily with excellent prep adaptation
- Bis-MEPP monomer allows for non-sticky placement of restorative
- May be used in either anterior or posterior, with no slumping



Courtesy of Dr A. Mak, Australia

Unishade Simplicity

- Mimics natural light reflection for invisible restorations
- Polishes with ease and retains high gloss finish
- Natural fluorescence in any light
- Reduced inventory need saves time and cost




High Performance

- High wear and stain resistance
- Long lasting gloss achieved with minimal polishing
- High radiopacity (318%) for easy identification on radiograph


Impress Your Patients With This Simplified Solution

G-aenial™ A'CHORD is a universal composite that delivers long lasting, invisible restorations. Designed to save time and therefore costs, it brings together superior technology in a straightforward product. Now you can feel truly confident using a material with effortless handling even in challenging cases.

Time-saving results: The remarkable high gloss achieved with minimal polishing.



Source: GC R&D, Japan, 2020. Data on file.



Source: GC R&D, Japan, 2020. Data on file.



G-aenial™ A'CHORD Syringe Refills (2.1ml/4g)

952230 A1	952233 A3.5	952236 A6	952240 AO3	952237 BW
952231 A2	952234 A4	952238 AO1	952242 JE	952241 BOW
952232 A3	952235 A5	952239 AO2	952243 AE	

G-aenial™ A'CHORD Unitip Refills (15 x 0.16ml)

952244 A1	952247 A3.5	952250 A6	952254 AO3	952251 BW
952245 A2	952248 A4	952252 AO1	952256 JE	952255 BOW
952246 A3	952249 A5	952253 AO2	952257 AE	

Scan for a Free Lunch & Learn!

