



PRESS RELEASE:  
Available 2/26/10

## NEW! GC Initial IQ – Lustré Pastes NF

GC Initial IQ Lustré Pastes NF, are 3-dimensional ceramic pastes developed to create color depth and lifelike translucency with a single paint on application. Thanks to the “New Formula” you can now really bring your crowns and bridges to life with this dedicated Lustré Paste Set, with an optimal match to both **low** and **high** CTE ceramics. Available in a ready to use consistency, the Lustré Pastes NF are based on fine ceramic particles and can be applied in a thicker layer where they exhibit unsurpassed vitality and a very natural glaze to your restorations.

- With an optimal match to both “low” & “high” CTE range, GC Initial IQ Lustré Pastes NF are compatible to almost all types of dental ceramic.
- GC Initial IQ Lustré Pastes NF are based on fine ceramic particles allowing a thicker application compared to conventional stains & glazes.
- Due to the unique paste medium, it has a fine thixotropic property for easy and exact applications.
- GC Initial IQ Lustré Pastes NF are low fusing, ceramic pastes that do not wear like some stains, glazes, or surface protectants.
- Just one dedicated kit for all different indications for the complete GC Initial ceramic range.
- Available in a ready to use consistency (1 x Neutral, 1 x Value, 4 x Body Shades, 8 x Enamel Effect Shades).

For more information call Dan Boskocevic at 800-323-7063 ext. 4043