



**Contact:**

Dan Boskocevic  
Senior Dental Laboratory Product Manager  
GC America Inc.  
800-323-7063 ext. 3012

**FOR IMMEDIATE RELEASE:**

**[GC America Inc. Introduces Additional GC Initial® LiSi Press LT Shades]**

Alsip, IL – February 2019 - GC America Inc. introduces GC Initial LiSi Press LT additional shades. A revolutionary new pressable ceramic, it combines unparalleled strength and amazing aesthetics. Best of all, it is fast with processing, optimized to be used with GC Initial® LiSi Veneering Ceramic and leaves virtually no reaction layer making your laboratory more productive.

Finally! A lithium disilicate ceramic with the aesthetics and strength technicians demand without being low in value.

GC Initial LiSi Press is the first lithium disilicate ceramic ingot with High Density Micronization (HDM), a technology unique to GC that provides exceptional physical properties and the most natural, life-like aesthetics. HDM uses equally dispersed lithium disilicate micro-crystals to fill the entire glass matrix rather than using traditional larger size crystals that do not take full advantage of the matrix structure. The result is the ultimate combination of strength and aesthetics, making GC Initial LiSi Press perfectly suitable for all types of restorations through all levels of transparency. Most importantly, HDM technology helps ensure the product remains extremely stable, without distortion or drop in value, even after multiple firings.

**Product will be available for Canada after February 25, 2019**

- Excellent flexural strength (>500 MPa) with HDM Technology
- Save up to 20 minutes per restoration
- Amazing aesthetics
- Rich, warm, and bright colors with excellent fluorescence
- Virtually no reaction layer when divested

Item Number/SKU	Product Description	
012835	GC Initial LiSi Press	LT-B00 (3gx5)
012837	GC Initial LiSi Press	LT-B0 (3gx5)
012838	GC Initial LiSi Press	LT-A1 (3gx5)
012839	GC Initial LiSi Press	LT-A2 (3gx5)
012840	GC Initial LiSi Press	LT-A3 (3gx5)
012843	GC Initial LiSi Press	LT-B1 (3gx5)
012844	GC Initial LiSi Press	LT-B2 (3gx5)
012847	GC Initial LiSi Press	LT-C1 (3gx5)
012848	GC Initial LiSi Press	LT-C2 (3gx5)
012851	GC Initial LiSi Press	LT-D2 (3gx5)

###