



GC AMERICA INC.
3737 W 127th STREET
ALSIP, ILLINOIS 60803
TEL (800) 323-3386
FAX (708) 897-4062
www.gcamerica.com

Contact:

Dan Boskocevic
Dental Laboratory Product Manager
GC America Inc.
800-323-7063. ext. 3012

FOR IMMEDIATE RELEASE:

Alsip, IL – December 2014 – GC America introduces GC Initial Zirconia Disk is now available for purchase in Canada.

An extremely high quality zirconium disk is now available to the GC Initial portfolio. The Zirconium Disk is Cold Isostatic Pressed (CIP) for optimized physical properties. There are two types to choose from: Standard Translucency (ST) is ideal for veneering and High Translucency (HT) is ideal for full contour, monolithic restorations.

Characteristics and benefits:

- Extremely high quality zirconium
- Cold Isostatic Pressed (CIP) for optimized physical properties
- Ideal for veneering
- Ideal for full contour monolithic restorations

The GC Initial Zirconia Disk will be available February 20, 2014.

Catalog Code No.	Product Description
878114	GC Initial Zirconia Disk ST 98.5mm x 12mm
878115	GC Initial Zirconia Disk ST 98.5mm x 14mm
878116	GC Initial Zirconia Disk ST 98.5mm x 18mm
878117	GC Initial Zirconia Disk ST 98.5mm x 25mm
878118	GC Initial Zirconia Disk HT 98.5mm x 12mm
878119	GC Initial Zirconia Disk HT 98.5mm x 14mm
878120	GC Initial Zirconia Disk HT 98.5mm x 18mm
878121	GC Initial Zirconia Disk HT 98.5mm x 25mm

###

