

**Contact:**

Toni Cusimano  
Associate Dental Lab Product Manager  
GC America Inc.  
800-323-7063 ext. 3006

**FOR IMMEDIATE RELEASE:**

**[Alsip, IL – February 2019 – GC America Inc. introduces GRADIA® Plus]**

**GRADIA® Plus** is a Modular Composite System for Indirect Restorations. Based on the latest ceramic polymer technology this advanced, high-strength, nano-hybrid, light-curing composite offers brightness, translucency, chroma and a natural opalescence in the oral environment that is similar to porcelain.

Its unique modular concept has fewer standard shades but uses a more individual mixing and layering approach making it more compact and cost-effective. Nevertheless, it will meet all demands for indications or technique used from classic or multi-chromatic build-up to monolithic approach.

Innovation meets indication: GRADIA Plus functions as a wide range of clinical applications from metal free inlays, veneers and crowns to frame supported crown and bridges along with implant superstructures. And GRADIA Plus has improved aesthetics and superior mechanical properties.

- State-of-the art ceramic polymer technology
- Red & white aesthetics in perfect harmony
- Indication related consistency of pastes
- Rational and modular
- Curing at the speed of light

Product is available to ship March 29, 2019.

SKU#:	Product Description:
901048	GRADIA Plus Layer Set
901409	GRADIA Plus Layer Pro Set
901050	GRADIA Plus Paint Set
901051	GRADIA Plus Gum Shade Set
901052	GRADIA Plus One Body Set
901053	GRADIA Plus Accessory Set

###





**GRADIA Plus kit contents and descriptions**

Product SKU#	Product Description:	Kit Contents:	Individual Replacement/SKU#	
901048	GRADIA Plus Layer Set	GRADIA Plus Opaque O-Base, 2.0 mL each	901054	
		GRADIA Plus Opaque OA, 2.0 mL each	901055	
		GRADIA Plus Opaque OB, 2.0 mL each	901056	
		GRADIA Plus Opaque OC, 2.0 mL each	901057	
		GRADIA Plus Opaque OD, 2.0 mL each	901058	
		GRADIA Plus Paste HB-DA1, 3.3 mL each	901069	
		GRADIA Plus Paste HB-DA2, 3.3 mL each	901070	
		GRADIA Plus Paste HB-DA3, 3.3 mL each	901071	
		GRADIA Plus Paste HB-DA3.5, 3.3 mL each	901072	
		GRADIA Plus Paste HB-DB1, 3.3 mL each	901073	
		GRADIA Plus Paste HB-DB3, 3.3 mL each	901074	
		GRADIA Plus Paste HB-DC3, 3.3 mL each	901075	
		GRADIA Plus Paste HB-DD2, 3.3 mL each	901076	
		GRADIA Plus Paste HB-EL, 3.3 mL each	901077	
		GRADIA Plus Paste HB-ED, 3.3 mL each	901078	
		GRADIA Plus Paste HB-CLF, 3.3 mL each	901079	
		GRADIA Plus Paste HB-PE, 3.3 mL each	901080	
		GRADIA Plus Paste HB-ODA, 3.3 mL each	901081	
		GRADIA Plus Paste HB-ODB, 3.3 mL each	901082	
			GRADIA Plus Dispensing Tip Needle Type Small & Light Protective Cover, 5 pc.	901134
	GRADIA Plus Dispensing Tip Needle Type & Light Prot.	901135		
	Mixing Pads #22 (1 pad in a kit)	001414		
	Plastic Spatula Blue	001769		
901409	GRADIA Plus Layer Pro Set	GRADIA Plus Paste HB-ODC, 3.3 mL each	901083	
		GRADIA Plus Paste HB-ODD, 3.3 mL each	901084	
		GRADIA Plus Paste HB-ODW, 3.3 mL each	901085	
		GRADIA Plus Paste HB-DW, 3.3 mL each	901086	
		GRADIA Plus Paste LB-Base E, 2 mL each	901091	
		GRADIA Plus Paste LB-Base CLF, 2 mL each	901092	
		GRADIA Plus Paste LB-Base D, 2 mL each	901093	
		GRADIA Plus Paste LB-DW, 2 mL each	901094	
		GRADIA Plus Paste LB-Base OD, 2 mL each	901095	
		GRADIA Plus Paste LB-ODW, 2 mL each	901096	
		GRADIA Plus Paste LB-Base Opal, 2mL each	901097	
		GRADIA Plus Paste LB-Orange, 2 mL each	901098	
		GRADIA Plus Paste LB-Yellow, 2 mL each	901099	
		GRADIA Plus Paste LB-Red, 2 mL each	901100	
		GRADIA Plus Paste LB-Grey, 2 mL each	901101	
		GRADIA Plus Paste LB-Blue, 2 mL each	901102	
		GRADIA Plus Paste LB-Milky, 2 mL each	901103	
		GRADIA Plus Paste LB-Inlay E, 2 mL each	901104	
		GRADIA Plus Paste LB-Inlay TD, 2 mL each	901105	
			GRADIA Plus Mixotip, (10 pcs.)	901136
			GRADIA Plus Dispensing Tip Wide Plastic Type & Light Protective Cover (5 pcs.)	901133
			GRADIA Plus Dispensing Tip Needle Type & Light Prot.	901135
	Mixing Pads #22 (1 pad in a kit)	001414		
	Plastic Spatula Blue	001769		



Product SKU#	Product Description:	Kit Contents:	Individual Replacement SKU#
901050	GRADIA Plus Paint Set	GRADIA Plus LP-A, 0.8 mL each	901115
		GRADIA Plus LP-B, 0.8 mL each	901116
		GRADIA Plus LP-C, 0.8 mL each	901117
		GRADIA Plus LP-D, 0.8 mL each	901118
		GRADIA Plus LP-CLF, 0.8 mL each	901119
		GRADIA Plus LP-Blue, 0.8 mL each	901120
		GRADIA Plus LP-Grey, 0.8 mL each	901121
		GRADIA Plus LP-Cream, 0.8 mL each	901122
		GRADIA Plus LP-Lavender, 0.8 mL each	901123
		GRADIA Plus LP-CL, 0.8 mL each	901124
		GRADIA Plus LP Diluting Liquid 3 mL each	901127
		GRADIA Plus Dispensing Tip Needle Type Small & Light Protective Cover (5 pcs.)	901134
		GRADIA Plus Dispensing Tip Needle Type & Light Prot.	901135
		Brush Round No.1 (10 pcs.)	001533
		Brush Flat No.1 (10 pcs.)	001534
		Brush Holder (2 pcs.)	001535
		Disposable Palette with Cover (5 pcs.)	N/A
		Mixing Pads No.14B (1 pad in kit)	423003
901051	GRADIA Plus Gum Shades Set	GRADIA Plus Opaque O-Base, 2.0 mL each	901054
		GRADIA Plus Opaque GO-1, 2.0 mL each	901067
		GRADIA Plus Opaque GO-2, 2.0 mL each	901068
		GRADIA Plus LP-CL, 0.8 mL each	901124
		GRADIA Plus GLP-Bright Red, 0.8 mL each	901125
		GRADIA Plus GLP-Violet, 0.8 mL each	901126
		GRADIA Plus Light Body GLB-1, 2 mL each	901111
		GRADIA Plus Light Body GLB-2, 2 mL each	901112
		GRADIA Plus Light Body GLB-3, 2 mL each	901113
		GRADIA Plus Light Body GLB-CL, 2 mL each	901114
		GRADIA Plus Heavy Body GHB-2, 3.3 mL each	901088
		GRADIA Plus Mixotip, (10 pcs.)	901136
		GRADIA Plus Dispensing Tip Wide Plastic Type & Light Protective Cover (5 pcs.)	901133
		GRADIA Plus Dispensing Tip Needle Type Small & Light Protective Cover, (5 pcs.)	901134
		GRADIA Plus Dispensing Tip Needle Type & Light Prot.	901135
		Brush Round No.1 (10 pcs.)	001533
		Brush Flat No.1 (10 pcs.)	001534
		Brush Holder (2 pcs.)	001535
		Disposable Palette with Cover (5 pcs.)	N/A
		Mixing Pads (1 pad in a kit)	001414
		Plastic Spatula Blue	001769



Product SKU#	Product Description:	Kit Contents:	Individual Replacement SKU#
901052	GRADIA Plus One Body Set	GRADIA Plus Opaque O-Base, 2.0mL each	901054
		GRADIA Plus Opaque OA, 2.0mL each	901055
		GRADIA Plus Opaque OB, 2.0mL each	901056
		GRADIA Plus Opaque OC, 2.0mL each	901057
		GRADIA Plus Opaque OD, 2.0mL each	901058
		GRADIA Plus Paste LB-Base OD, 2mL each	901095
		GRADIA Plus One Body LB-A, 2.0mL each	901106
		GRADIA Plus One Body LB-B, 2.0mL each	901107
		GRADIA Plus One Body LB-C, 2.0mL each	901108
		GRADIA Plus One Body LB-D, 2.0mL each	901109
		GRADIA Plus One Body LB-W, 2.0mL each	901110
		GRADIA Plus Dispensing Tip Wide Plastic Type & Light Protective Cover (5 pcs.)	901133
		GRADIA Plus Disposable Tip N-type & Light Protective Cover (5 pcs.)	901135
		Mixing Pads #22 (1 pad in a kit)	001414
Plastic Spatula Blue	001769		
901053	GRADIA Plus Accessory Set	GRADIA Plus Air Barrier, 10mL each	901128
		GRADIA Plus Modeling Liquid, 3mL each	901129
		GRADIA Plus Separator, 5mL each	901130
		GRADIA Plus Die-Hardner, 5mL each	901131
		Ceramic Primer II, 3 mL	008551
		Metal Primer Z, 5 mL	009295
		GC Acrylic Primer, 6 mL each	901138
		Diapolisher Paste Syringe	004415
		Brush Round No.1 (10 pcs.)	001533
		Brush Flat No.1 (10 pcs.)	001534
		Brush Holder (2 pcs.)	001535
		Brush No.7 Bulk	001536
		Mixing Pads #22 (1 pad in a kit)	001414

###