



PRESS RELEASE:
Available 10/10/11

GC America Introduces: Three New State of the Art Gypsums.

NEW FUJIROCK
NEW FUJIROCK IMP
NEW FUJIROCK FAST SET

NEW FUJIROCK- One of the world's most popular dental die stone manufacturers now brings you NEW FUJIROCK. NEW FUJIROCK is an improved type IV die stone in its physical properties and workability. This new die stone is ideal for all kinds of prosthetics, with high precision, outstanding edge hardness and high pressure stability. Thus a precise model as a basis for precise prosthetics is achieved.

- Easy to mix. Flows well yet runs less from the impression when pouring
- Low setting expansion (**0.08**) provides accurate models with smooth surfaces
- Increased compressive strength results in no fractures during removal from impression
- Available in Golden Brown & Alabaster Ivory - (Ideal for ceramic restorations)

NEW FUJIROCK IMP- One of the world's most popular dental die stone manufacturers now brings you NEW FUJIROCK IMP. NEW FUJIROCK IMP is an extremely precise and hard type IV die stone with very low expansion (0.06). Because of these very special properties NEW FUJIROCK IMP is an ideal stone for Implant cases, CAD CAM, and long-span superstructures which require a long span of accuracy and detail. A very precise die stone for very precise work.

- Easy to mix. Flows well yet runs less from the impression when pouring. Exhibits superior handling characteristics
- Low setting expansion (**0.06**) provides accurate models with smooth surfaces
- Increased compressive strength results in no fractures during removal from impression
- Available in Golden Brown & Alabaster Ivory - (Ideal for ceramic restorations)
- Specifically designed for, Implant, CAD CAM, and very precise long span restorations

NEW FUJIROCK FAST SET- One of the world's leading dental die stone manufacturers now brings you NEW FUJIROCK FAST SET. NEW FUJIROCK FAST SET is everything you would expect from the FUJIROCK line of high quality gypsum materials. NEW FUJIROCK FAST SET has all of the strong and precise qualities as NEW FUJIROCK but sets in only **5 minutes**.

- Sets in 5 minutes, yet provides proper compressive strength in a shorter time. Enables quality work to be done very quickly
- Easy to mix. Flows well yet runs less from the impression when pouring
- Low setting expansion (**0.08**) provides accurate models with smooth surfaces
- Increased compressive strength results in no fractures during removal from impression
- Available in Golden Brown
- Provides ample working time with a sharp set

For more information call Dan Boskocevic at 800-323-7063 ext. 3012