

Contact:

Toni Cusimano
Associate Dental Lab Product Manager
GC America Inc.
800-323-7063 ext. 3006



FOR IMMEDIATE RELEASE:

[GC America Inc. introduces PATTERN RESIN™ LS Mini]

Alsip, IL – February 2019, PATTERN RESIN LS Mini is a self-curing, general purpose acrylic resin used for making patterns for fabricating copings, post and core build-ups, lingual and palatal bars, implant attachments, adhesion bridges, clasps, telescoping crowns and certain pre-soldering techniques. The stability and handling properties of the resin allow the dental technician to create restorations with high detail and precision.

- Ash-free; leaves no ash or residue after burn out
- Results in smoother cast surfaces
- Low polymerization shrinkage
- Great accuracy and dimensional stability
- High strength; ultra smooth surface
- Virtually no distortion
- Excellent fluidity
- Quick setting time

PATTERN RESIN LS Mini will be available for purchase at Lab Day Chicago, March 1, 2019. Canadian launch to follow soon.

SKU#:	Product Description:
335207	PATTERN RESIN LS 1:1 Package Mini
335208	PATTERN RESIN LS Mini Liquid Refill 30mL
335209	PATTERN RESIN LS Mini Powder Refill 30g

###